



Introduction to
Extraction Chromatography Resins
and their
Use in Rapid Methods

Lawrence Jassin and Terence O'Brien

At the FERN National Training Conference
June 30th, 2010

Variety is the **spice** of life



However, it is a **nightmare** in the Lab

eichrom We Need Sample Preparation



eichrom



eichrom

Your separation resin drawer!



Hows

Versatility

Results

eichrom


Stops in the Lab/Kitchen along the way

- Sample Preparation Concepts
- Extraction Chromatography (EXC) Fundamentals
- Separation Tools – EXC Uptake Curves
- Some example separations and RESULTS
- Questions ?





Your goals for Sample Preparation

- Minimize sample prep time and waste generation
- Trade offs between
 - sample size, detection level, sample preparation rigor
- Data quality objectives
- Detection instrument selection
- Instrument's capabilities





eichrom


Measurement of Gamma Energy



Measurement of Beta Energy (LSC)




Measurement of Alpha Energy




Americium Spectrum after TRU Resin Separation

presence of rare earths degrades spectrum- self absorption issues



Am Spectrum after TEVA Resin Separation

Lanthanide elements removed - cleaner spectrum



Radiologic Screening Counters

- Gas-Flow Proportional Counters
 - Alpha/Beta Weight Attenuation Curves
 - Alpha/Beta Cross Talk Calibration Curves
 - Alpha/Beta Isotope Calibration
- Liquid Scintillation Counters
 - Alpha/Beta Efficiency Determination
 - Alpha/Beta Cross Talk
 - Alpha/Beta Quench/Color Correction





Bulk Sample Preparation Steps

- Volume Reduction
 - Precipitation
 - Evaporation
 - Ion Exchange
- Sample Modification
 - Digestion
 - Leach
 - Fusion





Actinides/Sr in Fish Method

- 200 g fish
- Wet ash
 - aqua regia/ $\text{HNO}_3/\text{H}_2\text{O}_2$
- Furnace
 - 550°C
- Dissolve in 12 ml 6M HNO_3 + 12 ml 2M $\text{Al}(\text{NO}_3)_3$ + 3M HNO_3 as needed (~40-45 ml load solution)



10. 10. 2006 13:53

FlowChart of with time frames Radiostrontium Emergency Method At SRS





Maxwell III, SL and Culligan, B.K., "Rapid Method for Determination of Radiostrontium in Emergency Milk Samples", Journal of Radioanalytical and Nuclear Chemistry, Vol. 279, 3 (2009) 757-760

Separation of target analyte(s) from interferences

- Precipitation or Co-Precipitation
 - Adjustment of sample conditions that causes target analytes to fall out of solution
- Ion Exchange
 - Retention of target analytes on a activated support
- Solvent Extraction
 - Liquid to Liquid Extraction Separation
- Extraction Chromatography
 - Retention of target analytes in a supported liquid to liquid system.

Metal Anion Complex Formation



Absorption of elements from HCl solutions by Anion Resin




Fig. 6. Anion exchange distribution coefficients in HCl solutions. (Dowex 1-X10) (Ref. 3)

Absorption of elements from HNO₃ solutions by Anion Resin




Fig. 7. Anion Exchange distribution coefficients in HNO₃ solutions. (Dowex 1-X10) (Ref. 4)

eichrom

Extraction Chromatography



Solvent Extraction




Column Chromatography

eichrom

Extraction Chromatographic Resin

Surface of Porous Bead



Inert support =

Macroporous Acrylic Resin

Example Stationary Phases

- **Crown Ether (Sr)**
- **CMPO (TRU)**
- **DAAP (UTEVA)**



eichrom
TECHNOLOGIES INC.
8005 South Coast Pkwy, Suite 111
Darien, IL 60561 USA
Tel: 630.963.0300
Fax: 630.963.1928

TEVA Resin
Part #: TE-8000
Lot: TEA0002
Amount: 200 grams

eichrom
TECHNOLOGIES INC.
8005 South Coast Pkwy, Suite 111
Darien, IL 60561 USA
Tel: 630.963.0300
Fax: 630.963.1928

eichrom
TECHNOLOGIES INC.
8005 South Coast Pkwy, Suite 111
Darien, IL 60561 USA
Tel: 630.963.0300
Fax: 630.963.1928

Sr Resin
Part #: SR-B100-A
Lot: SRA03032
Amount: 100 grams

eichrom
TECHNOLOGIES INC.
8005 South Coast Pkwy, Suite 111
Darien, IL 60561 USA
Tel: 630.963.0300
Fax: 630.963.1928

eichrom
TECHNOLOGIES INC.
8005 South Coast Pkwy, Suite 111
Darien, IL 60561 USA
Tel: 630.963.0300
Fax: 630.963.1928

TRU Resin SPS
Part #: TR-8000
Lot: TRS0000
Amount: 50 grams

eichrom
TECHNOLOGIES INC.
8005 South Coast Pkwy, Suite 111
Darien, IL 60561 USA
Tel: 630.963.0300
Fax: 630.963.1928

eichrom

Batch Reaction / MnO₂ Resin for Ra Analysis



eichrom

Column Format / Gravity Flow



Gravity Flow vs. Vacuum Assisted



TEVA Discs/ Gravity or Vacuum Assisted Flow



eichrom

Selecting
the
right
utensil
for
the
job



Knife to peel



Vegetables knife



Puntilla knife



Cutlet knife



Cleaner knife



Knife to fillet



Knife to carve



Cook knife



Salmon knife



Jam knife



Spatula



Bread knife



Fork to carve



Santoku knife



Kitchen axe

Our Products

Eichrom commercialized its breakthrough chemical separation techniques in the 1990s. Today, accredited laboratories and regulatory agencies worldwide trust Eichrom's line of analytical-chemistry products, all developed through the expertise of our team of dedicated chemists.

Product Name	Color Code	Applications*		
Accessories		Plastic accessories for columns and cartridges	product info	part numbers
Resolve™ Filters		Alpha spectroscopy source preparation	technical info	part numbers
Resolve™ PTFE Filters		- NEW Radiological Air Monitoring	technical info	part numbers
Actinide Resin	● Yellow	Group actinide separations/gross alpha measurements	technical info	part numbers
Beryllium Resin		Be	technical info	part numbers
DGA Resin		Actinides, Lanthanides, Y, Ra	technical info	part numbers
Diphonix® Resin		Actinides and transition metals	technical info	part numbers
Ion Exchange Resins		Analytical grade cation and anion exchange resins	technical info	part numbers
Ln Resin	● Purple	Lanthanides, Ra-228	technical info	part numbers
MnO ₂ Resin		Ra	technical info	part numbers
Nickel Resin	● Pink	Ni	technical info	part numbers
Pb Resin	● Black	Pb	technical info	part numbers
Pre-filter Material		Organics removal	technical info	part numbers
RE Resin		Th, U, Np, Pu, Am, Cm, rare earth elements	technical info	part numbers
Sr Resin	● Red	Sr, Pb	technical info	part numbers
TEVA® Resin	● Green	Tc, Th, Np, Pu, Am/Lanthanides	technical info	part numbers
Tritium Column		³ H	technical info	part numbers
TRU Resin	● Blue	Fe, Th, Pa, U, Np, Pu, Am, Cm	technical info	part numbers
UTEVA® Resin	● Orange	Th, U, Np, Pu	technical info	part numbers


*Primary applications shown in blue.

EXC Resin Extractant Choice

- Ionic Recognition Extractant based resin
 - Analyte retention related to charge and ionic radii
- Neutral and Anionic Extractant based resin
 - Analytes are directly hydrogen ion dependant, although anionic extractant based resins tend to show reduced analyte uptake as the acid competes
- Acidic Extractant based resin
 - Analyte retention is inverse hydrogen ion dependant

Sr Resin

di-t-butylcyclohexano 18-crown-6




Diluent: 1-octanol

Ionic
Recognition,
size and charge




Acid dependency of k' for various ions at 23-25°C.

Sr Resin



Anionic Extractant TEVA Resin


Trialkylmethylammonium Nitrate or Chloride



R = C₈H₁₇ and C₁₀H₂₁

Acid dependency of k' for various ions at 23°C.


TEVA Resin




eichrom

TRU Resin: Neutral Extractant



DGA Resin: Neutral extractant/ ionic recognition




TRU (CMPO)




DGA




Actinides on TRU vs DGA out of HNO_3






Actinide Resin uptake of various actinides with DIPEX® extractant (Liquid Chelating Exchanger)



Data developed at Argonne National Laboratory, USA

Pb and Po in milk, crop and animal samples (Dell)



Actinides/Sr in Fish Method



3 + 3M

- **Load to TEVA+TRU+DGA**
 - after valence adjustment
 - using sulfamic acid, iron (if Np-237 needed), ascorbic acid, followed by sodium nitrite
- **Collect load/rinse (evaporate and redissolve later in 8M HNO₃ for Sr Resin)**

Actinides and Sr-90 in Fish Data

Tracer/carrier Recovery	Avg. Recovery	MS
Pu-236	99.8%	100% (Pu-238) 90.0% (Np-237)
Am-243	109%	94.1% (Am-241) 94.3% (Cm-244)
U-232	97.1%	91.1% (U-235)
Sr carrier	84.9%	97.7% (Sr-90)

Typical Actinide Tracers and Sr Carrier Recoveries for Animal Tissue Matrices

<u>Matrix</u>	<u>Pu-236</u>	<u>Am-243</u>	<u>U-232</u>	<u>Sr Carrier</u>
Beef (N=6)	98.7% ±5.7%	97.1% ± 8.4%	93.4% ±4.7%	96.3% ±0.5%
Deer (N=59)	99.3% ±12%	93.4% ±10%	90.4% ±8.0%	83.4%±3.5%
Fish-Bass (N=72)	96.2% ±14%	101.8%±13%	95.1% ±8.1%	89.0% ±16%
Fish-Bream (N=57)	96.6% ±12%	98.4% ±7.7%	91.1% ±6.3%	91.7% ±10%
Fish-Catfish (N=69)	98.3% ±12%	103.7% ±7.6%	89.4% ±12%	89.4% ±17%
Fish-Mullet (N=6)	96.2% ±6.8%	100.4% ±8.9%	91.0% ±8.1%	85.6% ±17%
Fish-Red Fish (N=6)	99.5% ±11%	105.2% ±8.6%	95.7% ±3.2%	77.7% ±21%
Fish-Sea Trout (N=6)	100.5% ±5.0%	102.2% ±7.6%	83.5% ±20%	74.4% ±25%
Hog (N=17)	93.0% ±20%	96.4% ±9.7%	86.4% ±15%	86.0%±7.1%
Shellfish (N=5)	101.3 ±2.2%	97.4% ±7.1%	81.7 ±3.2%	97.5% ±0.89%

Composite tissue samples

100 gram-deer, hog, bream, shellfish

200 gram-catfish, bass, red drum, mullet, sea trout

25 gram- nonedible fish samples including bones

Sr-89/90 in Milk Column Extraction

- Redissolve in 10 ml 8M HNO_3 -1M $\text{Al}(\text{NO}_3)_3$
- Perform typical Sr Resin Separation using 3 ml Sr resin
 - (2 ml +1 ml cartridges)
- Rinses:
 - 15 mL of 8M HNO_3
 - 5 ml 3M HNO_3 -0.05M oxalic acid
 - 7 ml 8M HNO_3
- Sr Elution: 13 ml 0.05M HNO_3



Performance of New Radiostrontium in Milk - 60 minute Count

⁹⁰ Sr Added (Bq/L)	⁹⁰ Sr Measured (Bq/L)	Uncertainty (%, K=2)	Difference (%)
0	0.11	130	N/A
0	0.27	59	N/A
2.86*	3.09	13.2	+8.0
2.86*	3.11	16.7	+8.7
2.86*	2.67	13.6	-6.6
2.86*	2.67	11.3	-6.6
5.70	5.85	10.4	+2.6
5.70	5.75	8.3	+0.9
5.70	6.04	8.2	+5.9
14.3	13.6	6.1	-4.9
14.3	14.0	6.1	-2.1
14.3	14.2	6.1	-0.7

* Added using NRIP water standard

Average +0.52

eichrom

Your separation resin drawer!



Hows

Versatility

Results

eichrom

¿ Questions ?

