




# Extraction Chromatography Resins and their Use in Food Methods

Lawrence Jassin and Terence O'Brien

Winchester Engineering and Analytical Center  
(WEAC)

July 15<sup>th</sup>, 2010

Variety is the **spice** of life



However, it is a **nightmare** in the Lab

# eichrom We Need Sample Preparation





eichrom



eichrom

# Your separation resin drawer!



Hows

Versatility

Results

# eichrom


## Stops in the Lab/Kitchen along the way

- Why Sample Preparation / Instrumentation Choices
- Matrix Elimination Approaches
- Extraction Chromatography (EXC) Separations
- Some example separations and RESULTS
- Questions ?



# Your goals for Sample Preparation

- Minimize sample prep time and waste generation
- Trade offs between
  - sample size, detection level, sample preparation rigor
- Data quality objectives
- Detection instrument selection
- Instrument's capabilities





# Radiologic Screening Counters

- Gas-Flow Proportional Counters
  - Alpha/Beta Weight Attenuation Curves
  - Alpha/Beta Cross Talk Calibration Curves
  - Alpha/Beta Isotope Calibration
- Liquid Scintillation Counters
  - Alpha/Beta Efficiency Determination
  - Alpha/Beta Cross Talk
  - Alpha/Beta Quench/Color Correction






# Measurement of Beta Energy (LSC)




# Measurement of Alpha Energy




# Americium Spectrum after TRU Resin Separation

presence of rare earths degrades spectrum- self absorption issues



# Am Spectrum after TEVA Resin Separation



Lanthanide elements removed - cleaner spectrum





eichrom

# Measurement of Gamma Energy



# Actinides by ICP-MS: The issues

## Polyatomic ions

$^{230}\text{ThH}$  and  $^{231}\text{Pa}$

$^{236}\text{UH}$  and  $^{237}\text{Np}$

$^{237}\text{NpH}$  and  $^{238}\text{U}$

$^{238}\text{UH}$  and  $^{239}\text{Pu}$

$^{240}\text{PuH}$  and  $^{241}\text{Am}$

$^{242}\text{PuH}$  and  $^{243}\text{Am}$


## Isobars

$^{236}\text{U}$  and  $^{236}\text{Np}$

$^{238}\text{U}$  and  $^{238}\text{Pu}$

$^{241}\text{Pu}$  and  $^{241}\text{Am}$

$^{242}\text{Pu}$  and  $^{242}\text{Am}$





## Matrix Elimination Approaches

- Sample Modification
  - Digestion
  - Leach
  - Fusion
- Volume Reduction
  - Precipitation
  - Evaporation
  - Ion Exchange



# Actinides/Sr in Fish Method

- 200 g fish
- Wet ash
  - aqua regia/ $\text{HNO}_3/\text{H}_2\text{O}_2$
- Furnace
  - 550°C
- Dissolve in 12 ml 6M  $\text{HNO}_3$  + 12 ml 2M  $\text{Al}(\text{NO}_3)_3$  + 3M  $\text{HNO}_3$  as needed (~40-45 ml load solution)






# Sr in Milk

## Calcium Phosphate Precipitation

- Milk Sample
  - 500 ml aliquot routine
  - 100 mL emergency + 50 mL Water to facilitate separation
- Add 2 mL 1.25M  $\text{Ca}(\text{NO}_3)_2$  and 5 mL  $(\text{NH}_4)_2\text{HPO}_4$ 
  - Ca added so water blanks will precipitate (not really needed for milk)

500 ml milk



# Calcium Phosphate Precipitation

- Add phenolphthalein indicator
- Add  $\text{NH}_4\text{OH}$  to dark pink (pH 10)
- Centrifuge 10 minutes at 3500 rpm
- Discard the supernatant

# Acidification of Precipitant

Separate the fat/protein from the Sr

- Add 20 mL of 3M HNO<sub>3</sub> into each tube
- Ca<sub>3</sub>(PO<sub>4</sub>)<sub>2</sub> dissolves
- Fat/protein coagulate
- Centrifuge
- Transfer supernatant to beaker (contains the Sr)
- Rinse solids with 10-15 ml 3M HNO<sub>3</sub>
- Transfer supernatant to beaker (remaining Sr)
- Evaporate beaker to dryness


# Destruction of the remaining matrix

- Wet ash
  - 15 ml concentrated  $\text{HNO}_3$  and 5 ml 30 wt%  $\text{H}_2\text{O}_2$
- Heat beakers in a furnace
  - at 550C for 30-60 minutes to turn the solids white
- Wet ash
  - 10-15 ml concentrated  $\text{HNO}_3$  and 5 ml 30 wt%  $\text{H}_2\text{O}_2$
- Redissolve in 10 ml 8 M  $\text{HNO}_3$  -1M  $\text{Al}(\text{NO}_3)_3$



# Sample Preparation

500 ml sample ppt. shown



$\text{Ca}_3(\text{PO}_4)_2$  and  
fat/protein

Add 3M  $\text{HNO}_3$




# Sample Preparation




Centrifuge


Most fat/protein is removed



# Sample Preparation



Heat on hot plate



# Sample Preparation



Heat to dryness






# Sample Preparation





Heated at 550C 30 min.







# Flowchart of with time frames Radiostrontium Emergency Method At SRS



Maxwell III, SL and Culligan, B.K., "Rapid Method for Determination of Radiostrontium in Emergency Milk Samples", Journal of Radioanalytical and Nuclear Chemistry, Vol. 279, 3 (2009) 757-760

# Metal Anion Complex Formation



eichrom

# Extraction Chromatography



Solvent Extraction




Column Chromatography

eichrom

# Extraction Chromatographic Resin

Surface of Porous Bead



Inert support =

**Macroporous Acrylic Resin**

## Example Stationary Phases

- **Crown Ether (Sr)**
- **CMPO (TRU)**
- **DAAP (UTEVA)**



**eichrom**  
TECHNOLOGIES INC.  
8005 South Coast Pkwy, Suite 111  
Darien, IL 60561 USA  
Tel: 630.963.0300  
Fax: 630.963.1928

**TEVA Resin**  
Part #: TE-8000  
Lot: TEA0002  
Amount: 200 grams

**eichrom**  
TECHNOLOGIES INC.  
8005 South Coast Pkwy, Suite 111  
Darien, IL 60561 USA  
Tel: 630.963.0300  
Fax: 630.963.1928

**eichrom**  
TECHNOLOGIES INC.  
8005 South Coast Pkwy, Suite 111  
Darien, IL 60561 USA  
Tel: 630.963.0300  
Fax: 630.963.1928

**Sr Resin**  
Part #: SR-B100-A  
Lot: SRA03032  
Amount: 100 grams

**eichrom**  
TECHNOLOGIES INC.  
8005 South Coast Pkwy, Suite 111  
Darien, IL 60561 USA  
Tel: 630.963.0300  
Fax: 630.963.1928

**eichrom**  
TECHNOLOGIES INC.  
8005 South Coast Pkwy, Suite 111  
Darien, IL 60561 USA  
Tel: 630.963.0300  
Fax: 630.963.1928

**TRU Resin SPS**  
Part #: TR-8000  
Lot: TRS0000  
Amount: 50 grams

**eichrom**  
TECHNOLOGIES INC.  
8005 South Coast Pkwy, Suite 111  
Darien, IL 60561 USA  
Tel: 630.963.0300  
Fax: 630.963.1928



eichrom

## Batch Reaction / MnO<sub>2</sub> Resin for Ra Analysis



eichrom

## Column Format / Gravity Flow




eichrom

## Cartridge Format / Vacuum Assisted Flow





# Gravity Flow vs. Vacuum Assisted



# TEVA Discs/ Gravity or Vacuum Assisted Flow





eichrom

Selecting  
the  
right  
utensil  
for  
the  
job



Knife to peel



Vegetables knife



Puntilla knife



Cutlet knife



Cleaner knife



Knife to fillet



Knife to carve



Cook knife



Salmon knife



Jam knife



Spatula



Bread knife



Fork to carve



Santoku knife



Kitchen axe



Radiochemistry Products



Dioxin Test Kits



Ion Exchange Resins



Nuclear Power Outfitters



European Radioactivity Testing Laboratories

**LATEST NEWS:** Nuclear Power Outfitters Featured Product: "T-Flex" Tungsten Shielding System

**Technical Info**

[View Application Notes](#)

[View Methods](#)

[Search Bibliography](#)



**Order Online**

Free ground shipping now available for all orders placed online. [Click Here](#)



**Resources**

- > [MSDS](#)
- > [Certificates of Analysis](#)
- > [ISO 9000](#)
- > [About Us](#)
- > [Contact Us](#)



Product Name	Color Code	Applications*		
Accessories		Plastic accessories for columns and cartridges	<a href="#">product info</a>	<a href="#">part numbers</a>
Resolve™ Filters		Alpha spectroscopy source preparation	<a href="#">technical info</a>	<a href="#">part numbers</a>
Resolve™ PTFE Filters		• <b>NEW</b> Radiological Air Monitoring	<a href="#">technical info</a>	<a href="#">part numbers</a>
Actinide Resin	● Yellow	<b>Group actinide separations/gross alpha measurements</b>	<a href="#">technical info</a>	<a href="#">part numbers</a>
Beryllium Resin		<b>Be</b>	<a href="#">technical info</a>	<a href="#">part numbers</a>
DGA Resin		Actinids, Lanthanides, Y, Ra	<a href="#">technical info</a>	<a href="#">part numbers</a>
Diphonix® Resin		Actinides and transition metals	<a href="#">technical info</a>	<a href="#">part numbers</a>
Ion Exchange Resins		Analytical grade cation and anion exchange resins	<a href="#">technical info</a>	<a href="#">part numbers</a>
Ln Resin	● Purple	<b>Lanthanides, Ra-228</b>	<a href="#">technical info</a>	<a href="#">part numbers</a>
MnO <sub>2</sub> Resin		<b>Ra</b>	<a href="#">technical info</a>	<a href="#">part numbers</a>
Nickel Resin	● Pink	<b>Ni</b>	<a href="#">technical info</a>	<a href="#">part numbers</a>
Pb Resin	● Black	<b>Pb</b>	<a href="#">technical info</a>	<a href="#">part numbers</a>
Pre-filter Material		<b>Organics removal</b>	<a href="#">technical info</a>	<a href="#">part numbers</a>
RE Resin		Th, U, Np, Pu, Am, Cm, <b>rare earth elements</b>	<a href="#">technical info</a>	<a href="#">part numbers</a>
Sr Resin	● Red	<b>Sr, Pb</b>	<a href="#">technical info</a>	<a href="#">part numbers</a>
TEVA® Resin	● Green	<b>Tc, Th, Np, Pu, Am/lanthanides</b>	<a href="#">technical info</a>	<a href="#">part numbers</a>
Tritium Column		<sup>3</sup> H	<a href="#">technical info</a>	<a href="#">part numbers</a>
TRU Resin	● Blue	<b>Fe, Th, Pa, U, Np, Pu, Am, Cm</b>	<a href="#">technical info</a>	<a href="#">part numbers</a>
UTEVA® Resin	● Orange	<b>Th, U, Np, Pu</b>	<a href="#">technical info</a>	<a href="#">part numbers</a>


\*Primary applications shown in **blue**.

## EXC Resin Extractant Choice

- Ionic Recognition Extractant based resin
  - Analyte retention related to charge and ionic radii
- Neutral and Anionic Extractant based resin
  - Analytes are directly hydrogen ion dependant, although anionic extractant based resins tend to show reduced analyte uptake as the acid competes
- Acidic Extractant based resin
  - Analyte retention is inverse hydrogen ion dependant

# Sr Resin

di-t-butylcyclohexano 18-crown-6



**Diluent:** 1-octanol


Ionic  
Recognition,  
size and charge






Acid dependency of  $k'$  for various ions at 23-25°C.

Sr Resin





# Sr-89/90 in Milk Column Extraction

- Redissolve in 10 ml 8M HNO<sub>3</sub>-1M Al(NO<sub>3</sub>)<sub>3</sub>
- Perform typical Sr Resin Separation using 3 ml Sr resin – (2 ml +1 ml cartridges)
- Rinses:
  - 15 mL of 8M HNO<sub>3</sub>
  - 5 ml 3M HNO<sub>3</sub>-0.05M oxalic acid
  - 7 ml 8M HNO<sub>3</sub>
- Sr Elution: 13 ml 0.05M HNO<sub>3</sub>




## Performance of New Radiostrontium in Milk - 60 minute Count

<sup>90</sup> Sr Added (Bq/L)	<sup>90</sup> Sr Measured (Bq/L)	Uncertainty (%, K=2)	Difference (%)
0	0.11	130	N/A
0	0.27	59	N/A
2.86*	3.09	13.2	+8.0
2.86*	3.11	16.7	+8.7
2.86*	2.67	13.6	-6.6
2.86*	2.67	11.3	-6.6
5.70	5.85	10.4	+2.6
5.70	5.75	8.3	+0.9
5.70	6.04	8.2	+5.9
14.3	13.6	6.1	-4.9
14.3	14.0	6.1	-2.1
14.3	14.2	6.1	-0.7


\* Added using NRIP water standard

Average +0.52


# Acid dependency of $k'$ for various ions at 23-25°C. Sr Resin








## Pb and Po in milk, crop and animal samples (Dell)



# Pb, Po, U, Th, Ra


2.0 mL cartridge Sr Resin, 22(1)°C, 50-100 μm, 2 mL/min




## UTEVA Resin / Neutral Extractant

Diamyl Amylphosphonate (DAAP)

a.k.a. Dipentyl Pentylphosphonate (DPPP)



Acid dependency of  $k'$  for various ions at 23-25°C.  
UTEVA Resin






# UTEVA

*CJ. M. Oliveira and F. P. Carvalho*

Sequential extraction procedure for determination of U, Th, Ra, Pb and Po radionuclides by alpha spectrometry in Environmental Samples

Isotopic Tracer Recovery Yields (mean  $\pm 1\sigma$ ) for 19 Biota Samples

Analyte	Yield	Uncertainty
U	0.88	$\pm 0.12$
Th	0.47	$\pm 0.18$
Ra	0.22	$\pm 0.12$
Po	0.80	$\pm 0.08$
Pb	0.28	$\pm 0.24$



# Absorption of elements from HCl solutions by Anion Resin




Fig. 6. Anion exchange distribution coefficients in HCl solutions. (Dowex 1-X10) (Ref. 3)

# Absorption of elements from HNO<sub>3</sub> solutions by Anion Resin





Fig. 7. Anion Exchange distribution coefficients in HNO<sub>3</sub> solutions. (Dowex 1-X10) (Ref. 4)


Anionic Extractant  
TEVA Resin

Trialkylmethylammonium Nitrate or Chloride




R = C<sub>8</sub>H<sub>17</sub> and C<sub>10</sub>H<sub>21</sub>

Acid dependency of k' for various ions at 23°C.  
TEVA Resin



eichrom


Neutral Extractant - TRU Resin,  
Octyl(phenyl)-N,N-diisobutylcarbamoylmethylphosphine oxide  
(CMPO)




Diluent: tri-*n*-butyl phosphate



Acid dependency of  $k'$  for various ions at 23-25°C.  
TRU Resin



Acid dependency of  $k'$   
for various ions at 23°C.  
TRU Resin



# Actinides/Sr in Fish Method




- **Load to TEVA+TRU+DGA**
  - after valence adjustment
    - using sulfamic acid, iron (if Np-237 needed), ascorbic acid, followed by sodium nitrite
- **Collect load/rinse (evaporate and redissolve later in 8M HNO<sub>3</sub> for Sr Resin)**

## Actinides and Sr-90 in Fish Data

Tracer/carrier Recovery	Avg. Recovery	MS
Pu-236	99.8%	100% (Pu-238) 90.0% (Np-237)
Am-243	109%	94.1% (Am-241) 94.3% (Cm-244)
U-232	97.1%	91.1% (U-235)
Sr carrier	84.9%	97.7% (Sr-90)

# Actinide Column Separation

- 1) Redissolve in 7 mL warm 6M  $\text{HNO}_3$  and 7mL 2M  $\text{Al}(\text{NO}_3)_3$
- 2) Add 0.5 mL 1.5M Sulfamic Acid + 1.25 mL 1.5M Ascorbic Acid/ 1 mg Fe (if Np-237 analyzed)
- 3) Add 1 mL 3.5 M Sodium Nitrite







eichrom

TRU Resin: Neutral Extractant



DGA Resin: Neutral extractant/ ionic recognition




TRU (CMPO)



DGA



# Actinides on TRU vs. DGA out of $\text{HNO}_3$



- DGA alone:
  - Rinse DGA with 8 ml 0.1M HNO<sub>3</sub>
  - Add rinse to initial load/rinse solution containing Sr and evaporate
    - Contains some Sr
    - Place TRU cartridge above DGA and elute Am/Cm from TRU onto DGA with 15 ml 4M HCL
  - Discard rinse
- Rinse DGA (alone) with 3 ml 1M HNO<sub>3</sub>/10 ml 0.1M HNO<sub>3</sub> to remove interferences
  - Discard rinse
- Elute Am/Cm from DGA with 10 ml 0.25M HCl
  - Cerium fluoride microprecipitation
  - Alpha spectrometry

# eichrom Actinides/Sr in Fish Method

- Rinse TRU with 12 ml 4M HCL-0.2M HF
  - Th removal
- Elute U from TRU with 15 ml 0.1M ammonium bioxalate
  - Cerium fluoride precipitation/alpha spectrometry



## Typical Actinide Tracers and Sr Carrier Recoveries for Animal Tissue Matrices

<u>Matrix</u>	<u>Pu-236</u>	<u>Am-243</u>	<u>U-232</u>	<u>Sr Carrier</u>
Beef (N=6)	98.7% ±5.7%	97.1% ± 8.4%	93.4% ±4.7%	96.3% ±0.5%
Deer (N=59)	99.3% ±12%	93.4% ±10%	90.4% ±8.0%	83.4%±3.5%
Fish-Bass (N=72)	96.2% ±14%	101.8%±13%	95.1% ±8.1%	89.0% ±16%
Fish-Bream (N=57)	96.6% ±12%	98.4% ±7.7%	91.1% ±6.3%	91.7% ±10%
Fish-Catfish (N=69)	98.3% ±12%	103.7% ±7.6%	89.4% ±12%	89.4% ±17%
Fish-Mullet (N=6)	96.2% ±6.8%	100.4% ±8.9%	91.0% ±8.1%	85.6% ±17%
Fish-Red Fish (N=6)	99.5% ±11%	105.2% ±8.6%	95.7% ±3.2%	77.7% ±21%
Fish-Sea Trout (N=6)	100.5% ±5.0%	102.2% ±7.6%	83.5% ±20%	74.4% ±25%
Hog (N=17)	93.0% ±20%	96.4% ±9.7%	86.4% ±15%	86.0%±7.1%
Shellfish (N=5)	101.3 ±2.2%	97.4% ±7.1%	81.7 ±3.2%	97.5% ±0.89%


**Composite tissue samples**

**100 gram-deer, hog, bream, shellfish**


**200 gram-catfish, bass, red drum, mullet, sea trout**

**25 gram- nonedible fish samples including bones**





Actinide Resin uptake of various actinides with DIPEX® extractant (Liquid Chelating Exchanger)



Data developed at Argonne National Laboratory, USA

eichrom

## Your separation resin drawer!



Hows

Versatility

Results

eichrom

¿ Questions ?

