

# GROSS ALPHA DETERMINATION IN DRINKING WATER USING A HIGHLY SPECIFIC RESIN AND LSC

**S. HAPPEL\*<sup>a</sup>, P. LETESSIER<sup>a</sup>, J.H. EIKENBERG<sup>c</sup>, A.H. THAKKAR<sup>b</sup>, E.P. HORWITZ<sup>b</sup>**

<sup>a</sup>EICHROM EUROPE - Parc de Lormandière - Bât. C - Rue Maryse Bastié - Campus de Ker Lann - 35170 Bruz – France

<sup>b</sup>EICHROM TECHNOLOGIES - 8205 S. Cass Ave / Suite 106 - Darien, IL 60561 - USA

<sup>c</sup>PAUL SCHERRER INSTITUTE - Villigen PSI - 5232 - Switzerland

# Introduction

- Drinking Water Directive 98/83/EC
- Planchet Evaporation
- Resin / LSC Approach
  - Extraction
  - Interferences
  - LS counting
- Gross alpha protocol
- Real samples
- Conclusions / Next steps

# Drinking Water Directive 98/83/EC

- Covers water intended for human consumption
  - Water for drinking, food preparation and other domestic use
  - Implementation by latest 25 December 2003
- Radiological aspects implemented
  - Recommended by WHO guidelines
  - H-3 < 100 Bq L<sup>-1</sup>
  - Total indicative dose < 0.1 mSv a<sup>-1</sup>  
(H-3, K-40, Rn-222 & daughters excluded)
- Screening approach (guideline values, WHO):
  - 0.1 Bq L<sup>-1</sup> gross alpha activity
  - 1 Bq L<sup>-1</sup> gross beta activity
  - No further action if determined activities are below

# Planchet Evaporation

➤ Commonly used method

- Evaporation of 100 – 1000 ml water on 50 – 200 mm Ø planchets
- Gas proportional counting

➤ Pro:

- Number of samples counted simultaneously
- $\beta$ - $\alpha$ -spill-over

➤ Drawbacks:

- Evaporation time
- Dissolved solids (self absorption correction)
- Inhomogeneous distribution of solids => precision
- Detection efficiency / counting time

# Resin / LSC Approach


## ➤ Resin approach

- Large number of samples prepared simultaneously
- Batch extraction
- Matrix removal
- Resin and filter are dried and mixed with scintillation cocktail
  - Samples with very similar composition / SQP(E) value
  - No extensive quench correction

## Resin / LSC Approach


- $\alpha$ -/ $\beta$ -discrimination LS-counting
  - High detection efficiency
  - Low background count rates
  - Shorter counting times
  
- Used in routine in The Netherlands (Actinide Resin<sup>®</sup>)
  - Drinking and waste water
  - Drawback: Ra uptake interfered by Ca (>100 ppm)

# Extraction conditions



- pH 1: no Am uptake
- pH 2: overall good extraction
- No further improvement by shaking for 4 h


# Extraction at pH 2



- TracerMix: equal activities of Ra, U, Pu, Am
- pH 2: overall good extraction
- 30 min shaking time sufficient





# Ca-Interference





- U, Pu, Am, Np and Ra uptake high for varying Ca contents (> 90 – 100 %)
- Ra-uptake is good (> 90 %) even at high Ca contents (200 ppm)
  - Ca content of European Drinking water 30 – 150 ppm

# K and Fe interferences



- Only slight interference even for 500 ppm K
  - Drinking water content usually < 30 ppm K
- High Fe(III) amounts (100 ppm) interfere with Pu uptake
- Fe(II) shows no interference
  - 0.2 ppm Fe allowed (Drinking Water Directive 98/83/EC)

# Carbonate and Sulphate interferences




- Carbonate and Sulphate interfere when present in high concentrations (e.g. carbonate > 2000 ppm)
- interference is increased by large amounts of Ca (350 ppm)
- Carbonate content of drinking water usually < 1000 ppm
- Sulphate content restricted to 240 ppm (Drinking Water Directive 98/83/EC)

# Resin approach

- Extraction near quantitative
  - Actinides and Ra
  - pH 2 (commonly used)
  - 30 min shaking time
- 60 samples of different chemical composition
  - Chemical Yield (mean): 0,9
  - Reproducibility  $s_R = 5 \%$  (k=1)
  - Repeatability  $s_r = 5 \%$  (k=1)
- Interferences
  - Ca, Fe(III), carbonate and sulphate do not interfere  
(except when present in contents usually not found or not allowed in water)
  - No K interference
- Resin approach suited

# LS-counting



- $\alpha$ -/ $\beta$ -discrimination LS-counting (1220 Quantulus)
- 18 mL UltimaGold AB, dried Resin and filter (Am-241, Sr-90)
- Good reproducibility of samples SQP(E) value (800 +/- 10)



- High PSA values (>200) best suited
  - Low  $\beta$ -spill-over
  - $\alpha$  detection efficiency still high (60 - 80 % )

# LS-counting

➤ High PSA values (>200) best suited



- FoM ( $\epsilon_{\alpha}^2/\epsilon_{\beta f}$ ) values > 100
- Good detection limits
  - 40 mBq L<sup>-1</sup> (< 4 h, 100 mL sample)
  - 90 mBq L<sup>-1</sup> (< 90 min, 100 mL sample)

# Calculations

$$c_{A_\alpha} = \frac{R_\alpha - \left( \frac{\varepsilon_{\beta f} R_\beta}{\varepsilon_\beta} \right)}{\left( \varepsilon_\alpha - \left( \frac{\varepsilon_{\beta f} \varepsilon_{\alpha f}}{\varepsilon_\beta} \right) \right) \cdot V \cdot 60 \cdot R}$$

**Gleichung 1:** Gesamt-Alpha-Aktivitätskonzentration,  $\alpha$ - $\beta$ -LSC in Bq.L<sup>-1</sup>

$$\frac{U_{c_{A_\alpha}}}{c_{A_\alpha}} = k \cdot \sqrt{\left[ \frac{\sqrt{sR_\alpha^2 + \left( \frac{s\varepsilon_{\beta f} + \frac{s\varepsilon_\beta}{\varepsilon_\beta} + \frac{sR_\beta}{R_\beta} \right) \cdot \frac{\varepsilon_{\beta f} R_\beta}{\varepsilon_\beta}}}{\left( R_\alpha - \left( \frac{\varepsilon_{\beta f} R_\beta}{\varepsilon_\beta} \right) \right)} \right]^2 + \left[ \frac{\sqrt{s\varepsilon_\alpha^2 + \left( \frac{s\varepsilon_{\beta f}}{\varepsilon_{\beta f}} + \frac{s\varepsilon_\beta}{\varepsilon_\beta} + \frac{s\varepsilon_{\alpha f}}{\varepsilon_{\alpha f}} \right) \cdot \frac{\varepsilon_{\beta f} \varepsilon_{\alpha f}}{\varepsilon_\beta}}}{\left( \varepsilon_\alpha - \left( \frac{\varepsilon_{\beta f} \varepsilon_{\alpha f}}{\varepsilon_\beta} \right) \right)} \right]^2 + \left( \frac{u_V}{V} \right)^2 + \left( \frac{s_R}{R} \right)^2 \right]}$$

**Gleichung 2:** Unsicherheit Gesamt-Alpha-Aktivitätskonzentration,  $\alpha$ - $\beta$ -LSC in Bq.L<sup>-1</sup>

$$c_{\min}^* = \left( k_{1-\alpha} \cdot \sqrt{\left( \varepsilon_{\beta f} R_\beta + R_{\alpha,0} \right) \cdot \left( \frac{1}{t_0} + \frac{1}{t_m} \right)} + \frac{1}{4} (k_{1-\alpha} + k_{1-\beta})^2 \left( \frac{1}{t_0} + \frac{1}{t_m} \right) \right)$$


**Gleichung 3:** Erkennungsgrenze  $c_{\min}^*$  Gesamt-Alpha-Aktivitätskonzentration,  $\alpha$ - $\beta$ -LSC in Bq.L<sup>-1</sup>

$$c_{\min} = \frac{1}{\varepsilon_\alpha \cdot V \cdot 60 \cdot R} \left( (k_{1-\alpha} + k_{1-\beta}) \sqrt{\left( \varepsilon_{\beta f} R_\beta + R_{\alpha,0} \right) \cdot \left( \frac{1}{t_0} + \frac{1}{t_m} \right)} + \frac{1}{4} (k_{1-\alpha} + k_{1-\beta})^2 \left( \frac{1}{t_0} + \frac{1}{t_m} \right) \right)$$

**Gleichung 4:** Nachweisgrenze  $c_{\min}$  Gesamt-Alpha-Aktivitätskonzentration,  $\alpha$ - $\beta$ -LSC in Bq.L<sup>-1</sup>

$\varepsilon_\alpha$ : Korrekt klassifizierte  $\alpha$ -particles;  $\varepsilon_{\alpha f}$ : Inkorrekt klassifizierte  $\alpha$ -particles;  $\varepsilon_\beta$ : Korrekt klassifizierte  $\beta$ -particles;  $\varepsilon_{\beta f}$ : Inkorrekt klassifizierte;  $Z_\alpha / Z_\beta$ : Alpha/Beta-Fenster netto Zählrate;  $k_{1-\alpha}$  and  $k_{1-\beta}$ : Quantile der Normalverteilung (3,0 und 1,645),  $Z_0$ : Untergrundzählrate,  $t_0$  and  $t_m$ : Messzeit Untergrund/Probe.

# Gross alpha protocol





# Real samples

**table 1:** determined and reference values, gross-alpha determination, real samples

Probe	reference value / Bq L <sup>-1</sup>	determined value / Bq L <sup>-1</sup>	d/r / %
Marburg	404 +/- 11 <sup>+</sup>	408 +/- 18	100,9 +/- 4,5 %
Dueseldorf	414 +/- 11 <sup>+</sup>	408 +/- 19	98,6 +/- 4,5 %
Kaiserquelle Aachen	401 +/- 11 <sup>+</sup>	391 +/- 18	97,5 +/- 4,4 %
BfS 3/03 MW	4,32 +/- 0,11 <sup>*</sup>	3,94 +/- 0,14	91,2 +/- 4,0 %
BfS 3/03 RW	7,72 +/- 0,27 <sup>*</sup>	7,16 +/- 0,39	92,7 +/- 4,4 %

+ TracerMix-activity added

\* obtained using established evaporation protocol

- Good recoveries for spiked drinking water samples
- Good agreement between LSC and evaporation results

# Conclusions

- Uptake of Actinides is good; high Ca contents do not interfere
- Resin shows high Ra-uptake (> 95%) for low Ca-content; for high Ca-contents (200 ppm) uptake is > 90 %
- Fe(III), carbonate and sulphate do not interfere  
(except when present in contents not found or not allowed in water)
- Direct  $\alpha$ -/ $\beta$ -discrimination LS-counting of dried resin
  - good reproducibility of SQP(E)
  - LoD of 40 mBq L<sup>-1</sup> in less than 4 h counting time
  - LoD of 90 mBq L<sup>-1</sup> in less than 90 min counting time
- Good precision
- Quick results
- First test with real samples showed good results

## Next steps

- Improvement of LS-counting conditions
- Further evaluation with intercomparison/reference samples
- CETAMA Intercomparison
  - Comparison between Resin vs. Planchet Evaporation
  - Mineralwater
- Seawater