

# What Does Eichrom's Quality Control System Do For You?


Sarah Teichma, Joel Williamson, Jill Bryant, Terence O'Brien, Lawrence Jassin,  
E. Philip Horwitz, Daniel McAlister



Eichrom Technologies, LLC

56<sup>th</sup> Annual Radiobioassay & Radiochemical  
Measurement Conference

October 25-28, 2010



# Outline

- QC Goals
- Standard Testing
  - Finished Product Testing
  - Individual Constituent Testing
- QC Progression of a New Resin
- System Improvements
- Upcoming QC Upgrades

# Eichrom's Key QC Goals

1. Ship at least 98% of line items by the customer's required delivery date.

- 97% (To Date)

2. Implement one or more improvements in product quality, measurement systems of product quality, or efficiency quarterly.


(2010 highlights)

- Shelf Life Study Completed for DGA, Sr, TRU, TEVA, & UTEVA
- Reformulation of Tritium Columns

# Finished Product Testing

- Dry Weight Distribution Ratio ( $D_w$ )
  - Since the formation of Eichrom, an essential feature of our quality control system was the measurement of  $D_w$  for each batch of new resin manufactured.
  - Radioactive or Stable Elements

## Dry Weight Distribution Ratio




$$D_w = \frac{A_0 - A_s}{w(g)} \bigg/ \frac{A_s}{v(\text{mL})}$$

# Finished Product Testing (cont.)


## Column Elutions

- DGA  
Am & U
- Sr  
Ba, Ca, Sr, & Y
- TEVA  
Pu & Th
- TRU  
Am & Pu
- UTEVA  
Th & U



# Individual Constituent Testing


- Extractants & Diluents
  - Solvent Extraction
    - CMPO/TBP (TRU Resin)
    - Aliquat-336 (TEVA Resin)
    - DAAP (UTEVA Resin)



- Organic Phase
- Aqueous Phase
- Am 241

# Individual Constituent Testing (cont.)

- Support Resin
  - Flow Rate
    - A-Grade Specs (0.6-0.8mL/min.)
  - Particle Characterization
    - New Lot





# QC Progression of a New Resin

- Development of the Resin
- QC Testing Based on Application of Resin
- Collaborating w/ Input from Client
  - Ensuring the resin will meet their needs.
- Development of a Standard QC
- Setting of Performance Specifications

# QC Improvements

- Constantly Striving to Deliver the Best Product
- New Material
- Change of Material Supplier
- Inconsistent Data

# New Material

- Resolve Filters
  - Supplier Changed Manufacturing Process
  - Did Not Meet Needs
  - Designed a rugged enough QC to differentiate between materials.
  - Finally found a material that was comparable.

# New Material

- Tritium Columns
  - Supply of 100-200 mesh Diphonix was low, almost non-existent.
  - Batch Uptakes
  - Performance Directly Related to Flow Rate
  - 50-100 mesh Diphonix Did Not Perform Well
  - Mixed Particle Diphonix + S-Grade Pre-Filter
  - Reformulation/Collaborations with Clients => Better Product (Better Results/Faster)

# Supplier Change

- CMPO (TRU Resin)
  - New Supplier After 15+ Years
  - Ensure Quality/Purity
  - Th-227 Column Elution
    - Extremely High Specific Activity

# Inconsistent Data

- Sample From Larry's Trunk
- Matt O'Hara Collaboration
  - Inconsistent Results Using DGA
- Original QC: Column Elution w/Eu-152
- New QC: Column Elution w/ Am-241 & U-233
- Inconsistency Observed w/ U-233=> Purifying Extractant
  - Consistent Resin

**Load:** 3M HNO<sub>3</sub> **U Strip:** 0.5M HNO<sub>3</sub> **Am Strip:** 0.5M HCl

# Upcoming QC Upgrades



- LN QC Procedure
- Anion Procedure
  - Increase Ruggedness
  - New Needs
    - Client
    - Field Changes
  - Interferences

# Conclusions


- Standard testing is done to ensure the consistency of a product before the end user receives it.
- Eichrom always has open ears. If there is a comment or concern, please share.
- Specifications associated with the QC performed on the major resins are in the handout.
- CoAs Online



## Standard QC Specifications

<b>DGA-Normal Resin</b>	<sup>241</sup> Am Strip	<sup>241</sup> Am Breakthrough	<sup>233</sup> U Strip	<sup>233</sup> U Breakthrough
Mean (%)	97	1	97	1
Minimum (%)	95		95	
Maximum (%)		2		2
<b>MnO<sub>2</sub> Resin</b>	<sup>133</sup> Ba Strip	<sup>133</sup> Ba Breakthrough		
Mean (%)	95	1		
Minimum (%)	90			
Maximum (%)		2		
<b>Sr Resin</b>	Ba Breakthrough	Ca Breakthrough	Sr Strip	Y Breakthrough
Mean (A-Grade%/S-Grade%)	0.5/0.25	0.25/0.15	85/85	0.25/0.15
Minimum (A-Grade%/S-Grade%)			80/80	
Maximum (A-Grade%/S-Grade%)	1.5/0.5	0.75/0.25		0.75/0.25
<b>TEVA Resin</b>	<sup>239</sup> Pu Strip	<sup>239</sup> Pu Breakthrough	<sup>230</sup> Th Strip	<sup>230</sup> Th Breakthrough
Mean (%)	90	1	93	2
Minimum (%)	85		88	
Maximum (%)		2		5
<b>TRU Resin</b>	<sup>241</sup> Am Strip	<sup>241</sup> Am Breakthrough	<sup>239</sup> Pu Strip	<sup>239</sup> Pu Breakthrough
Mean (%)	95	1	95	1
Minimum (%)	90		90	
Maximum (%)		2		2
<b>UTEVA Resin</b>	<sup>230</sup> Th Strip	<sup>230</sup> Th Breakthrough	<sup>233</sup> U Strip	<sup>233</sup> U Breakthrough
Mean (%)	90	0.5	90	1
Minimum (%)	85		85	
Maximum (%)		0.7		2

### Sr Resin Recoveries



# Online: Certificate of Analysis



Radiochemistry Products



Dioxin Test Kits



Ion Exchange Resins



Nuclear Power Outfitters



European Radioactivity Testing Laboratories

**LATEST NEWS:** Nuclear Power Outfitters Featured Product: "T-Flex" Tungsten Shielding System

## Technical Info

[View Application Notes](#)

[View Methods](#)

[Search Bibliography](#)



## Order Online

Free ground shipping now available for all orders placed online. [Click Here](#)



## Resources

- > [MSDS](#)
- > [Certificates of Analysis](#)
- > [ISO 9000](#)
- > [About Us](#)
- > [Contact Us](#)





WARNING:

BE CAREFUL WHAT YOU GET OUT OF LARRY'S  
TRUNK!

