

# URANIUM IN SOIL

(2 GRAM SAMPLE)

## 1. SCOPE

- 1.1. This is a procedure for the separation of uranium from 2 gram soil samples. After separation of uranium with this method, source preparation for uranium measurement by alpha spectrometry is performed by electrolytic deposition onto stainless steel planchets (Eichrom Method SPA02) or by rare earth fluoride micro precipitation onto polypropylene filters (Eichrom Method SPA01).
- 1.2. This method does not address all aspects of safety, quality control, calibration or instrument set-up. However, enough detail is given for a trained radiochemist to achieve accurate and precise results for the analysis of the analyte(s) from the appropriate matrix, when incorporating the appropriate agency or laboratory safety, quality and laboratory control standards.

## 2. SUMMARY OF METHOD

- 2.1. Uranium is separated on Eichrom UTEVA resin prior to measurement by alpha spectrometry. A  $^{232}\text{U}$  tracer is used to monitor chemical recovery and correct results to improve precision and accuracy. A  $^{232}\text{U}$  self-cleaning tracer is recommended for this method. Preparation of the  $^{232}\text{U}$  self-cleaning tracer is described in Eichrom Method TP01.

## 3. SIGNIFICANCE OF USE

- 3.1. This is a rapid, reliable method for measurement of uranium in soil samples that is more cost-effective and efficient than traditional ion exchange, solvent extraction and precipitation techniques.

## 4. INTERFERENCES

- 4.1. Actinides with un-resolvable alpha energies such as  $^{241}\text{Am}$  and  $^{238}\text{Pu}$  or  $^{237}\text{Np}$  and  $^{234}\text{U}$  must be chemically separated to enable measurement. This method separates these isotopes effectively.
- 4.2.  $^{210}\text{Po}$  has an alpha energy similar to  $^{232}\text{U}$ . For samples containing very high levels of  $^{210}\text{Po}$  (e.g., environmental samples with the natural uranium decay chain intact or dissolved air filter samples) an additional column rinse is recommended.

4.3. This method may not adequately dissolve uranium in soils containing refractory particles. If refractory particles may be present in samples, a soil fusion method, as described in references (4) and (5), is recommended.

## 5. APPARATUS

- Analytical Balance, -0.0001 gram Sensitivity
- Beakers, glass
- Beakers, Teflon
- Centrifuge tubes, 50mL
- Centrifuge, rotor and carriers for 50mL tubes
- Column rack, Eichrom Part: AC-103
- Extension funnels, 25 mL, Eichrom Part: AC-120
- Fume Hood
- Hotplate
- Muffle Oven

## 6. REAGENTS

**Note:** Analytical grade or ACS grade reagents are recommended. Evaluation of key reagents, such as aluminum nitrate and ammonium hydrogen phosphate, for contribution to method background levels from naturally occurring radioactive materials is recommended.

Aluminum nitrate nonahydrate, $\text{Al}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$
Ammonium hydroxide (57% $\text{NH}_4\text{OH}$ or 28% $\text{NH}_3$ ), Concentrated $\text{NH}_4\text{OH}$
Appropriate tracers or standards ( $^{232}\text{U}$ self-cleaning tracer)
Deionized water, All reagents are prepared with deionized water
Hydrochloric acid (37%), concentrated HCl
Hydrofluoric acid (49%), concentrated HF
Nitric acid (70%), concentrated $\text{HNO}_3$
Oxalic acid dihydrate, $\text{H}_2\text{C}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$
UTEVA® Resin - 2mL Prepacked Column, 100-150 $\mu\text{m}$ , Eichrom Part UT-C50-A

6.1. *Hydrochloric Acid (1M)* - Add 83mL of concentrated HCl to 900mL of water. Dilute to 1L with water.

- 6.2. *Hydrochloric Acid (5M) - Oxalic Acid (0.05M)* - Dissolve 6.3g oxalic acid dihydrate in 400mL of water. Add 417mL concentrated HCl. Dilute to 1L with water.
- 6.3. *Hydrochloric Acid (9M)* - Add 750mL of concentrated HCl to 100mL of water. Dilute to 1L with water.
- 6.4. *Nitric Acid (3M)* - Add 188mL of concentrated HNO<sub>3</sub> to 700mL water. Dilute to 1L with water.
- 6.5. *Nitric Acid (8M)* - Add 500mL of concentrated HNO<sub>3</sub> to 450mL of water. Dilute to 1L with water.
- 6.6. *Nitric Acid (3M) - Aluminum Nitrate (1M)* - Dissolve 375g of Al(NO<sub>3</sub>)<sub>3</sub>·9H<sub>2</sub>O in 500mL of water, add 188mL of concentrated HNO<sub>3</sub>. Dilute to 1L with water.

## 7. PROCEDURE

### 7.1. Sample Preparation

- 7.1.1. Using an analytical balance, weigh up to 2g of sample into a 200mL glass beaker.
- 7.1.2. Heat the sample at 110°C until dry.
- 7.1.3. Weigh the sample again to achieve dry weight.
- 7.1.4. Spike samples with the appropriate yield tracer.

**Note: If using a self-cleaning <sup>232</sup>U tracer (Eichrom Method TP01), prior to removing aliquot for spiking soil samples, vortex the tracer for 1-2 minutes to suspend BaSO<sub>4</sub> precipitant and then centrifuge for 5 minutes. This will ensure that <sup>228</sup>Th and its daughters are effectively removed.**

- 7.1.5. Place the dry sample in a muffle furnace and ash overnight at 510°C.
- 7.1.6. Transfer the ashed soil sample to a 125mL Teflon beaker using 10mL of concentrated HNO<sub>3</sub>.
- 7.1.7. Add an additional 10mL of concentrated HNO<sub>3</sub> and 5mL of concentrated HCl to each beaker.
- 7.1.8. Place a watch glass on each beaker, and heat each sample to near boiling on a hot plate for 3 hours.

- 7.1.9. Cool the sample and transfer, including solids, to a 50mL centrifuge tube. Complete transfer with deionized water. Dilute to 50mL with deionized water.
- 7.1.10. Centrifuge and decant the supernate to a clean labeled Teflon beaker, leaving the residue in the bottom of the centrifuge tube. Set the beaker aside, while the residue is treated in the following steps.
- 7.1.11. Transfer the residue to the original 125mL Teflon beaker using 10 mL of concentrated HNO<sub>3</sub>. Add an additional 10mL of concentrated HNO<sub>3</sub> and 15mL of concentrated HF to the beaker.
- 7.1.12. Place a Teflon cover on each beaker and heat on a hot plate until the residue is dissolved.
- 7.1.13. Remove the Teflon cover, and add the solution containing the dissolved residue to the supernate in the beaker from step 7.1.9. (If any residue remains, repeat the HNO<sub>3</sub>/HF digestion from steps 7.1.11. and 7.1.12.)
- 7.1.14. Evaporate the combined solutions in the beaker from step 7.1.13 to dryness.
- 7.1.15. Remove beaker from hot plate and add 5mL of conc. HNO<sub>3</sub>.
- 7.1.16. Evaporate to dryness and dissolve in 10mL of 3M HNO<sub>3</sub>/1M Al(NO<sub>3</sub>)<sub>3</sub>. Transfer to a 50mL centrifuge tube.
- 7.1.17. Rinse the beaker with an additional 5mL of 3M HNO<sub>3</sub>/1M Al(NO<sub>3</sub>)<sub>3</sub> solution. Add rinse to the centrifuge tube, cap the centrifuge tube, and mix by swirling. Cool to room temperature and centrifuge the solution before loading on the columns.

## 7.2.U Separation Using UTEVA Resin

- 7.2.1. For each sample solution, place a UTEVA Resin column (with extension funnel) in the column rack.
- 7.2.2. Place a waste tray below the columns, remove the bottom plugs from each column, push the top frit down to the top of the resin bed, and allow to drain.

**Note: Make sure all reagent solutions have cooled to room temperature before proceeding with the method.**

- 7.2.3. Add 5mL of 3M HNO<sub>3</sub> into each column to precondition resin and allow to drain.
- 7.2.4. Transfer each solution from step 7.1.17. into the appropriate UTEVA Resin column . Allow the solution to drain by gravity flow. Uranium will be retained by the resin.
- 7.2.5. Add 5mL of 3M HNO<sub>3</sub> to rinse to each sample tube. Transfer each rinse solution into the appropriate UTEVA Resin column. Allow solution to drain.
- 7.2.6. Add 5mL of 3M HNO<sub>3</sub> into each column. Allow solution to drain.
- 7.2.7. Add 15mL of 8M HNO<sub>3</sub> to each UTEVA column. Allow solution to drain.

**Note: This rinse removes Po isotopes, including <sup>210</sup>Po which can affect measurement of <sup>232</sup>U by alpha spectrometry.**

- 7.2.8. Add 5mL of 9M HCl into each column. Allow solution to drain.

**Note: This rinse converts the resin to the chloride system. Some Th and Np may be removed here.**

- 7.2.9. Add 20mL of 5M HCl-0.05M oxalic acid into each column. Allow solution to drain. Discard the combined eluate to this point as waste.

**Note: This rinse removes plutonium, neptunium and thorium from the column.**

- 7.2.10. If using electrodeposition to prepare sources for alpha spectrometry, place a clean, labeled beaker below each column. If using rare earth fluoride micro precipitation to prepare sources for alpha spectrometry, place a clean, labeled 50mL polypropylene centrifuge below each column.
- 7.2.11. Add 15mL of 1M HCl into each column to strip the uranium. Allow solution to drain.
- 7.2.12. Prepare sources for uranium measurement by alpha spectrometry using electrodeposition (Eichrom Method SPA02) or rare earth fluoride micro precipitation (Eichrom Method SPA01).

## 8. CALCULATIONS

Calculate the actinide activity as follows:

Calculate tracer yield:

$$\text{Yield} = \frac{(C_s - B_s)}{E_s \times A_s}$$

where:

$C_s$  = measured actinide tracer, cpm

$B_s$  = background, cpm

$E_s$  = counting efficiency for tracer

$A_s$  = tracer activity, dpm

**Note: If any tracer may be present in the sample, a spiked and unspiked sample must be analyzed to determine chemical yield, where:**

$$Y = \frac{(\text{spiked sample tracer cpm} - \text{unspiked sample tracer cpm})}{E \times \text{actinide spike activity, dpm}}$$

$$\text{Percent yield} = \text{Yield} \times 100$$

Calculate actinide isotope activity:

$$\text{Sample dpm/L} = \frac{S - B}{E \times V \times Y}$$

where:

$S$  = sample activity, cpm

$B$  = background, cpm

$E$  = counting efficiency

$V$  = sample weight, g or volume, L

$Y$  = yield

Conversion of dpm/g to pCi/gram:

$$\text{pCi/g} = (\text{dpm/g}) / 2.22$$

## 9. REFERENCES

- 1) ASTM Method C1000-11, "Standard Method for the Radiochemical Determination of Uranium Isotopes in Soil by Alpha Spectrometry."
- 2) Horwitz, E.P., Chiarizia, R., Dietz, M.L., Diamond, H., Essling, A.M., and Gracyk, D.W. "Separation and Preconcentration of Uranium from Acidic Media by Extraction Chromatography," *Analytica Chimica Acta*, 266, 25-37 (1992).
- 3) ASTM Method D3648-14, "Standard Practices for the Measurement of Radioactivity."
- 4) Maxwell, S.L., Hutchinson J.B., McAlister D.R. "Rapid Fusion Method for the determination of refractory Th and U in soil samples," *Journal of Radioanalytical and Nuclear Chemistry*, in press, (2015).
- 5) Eichrom Application Note AN-1430, "Rapid Determination of Actinides in Soil Samples".
- 6) Eichrom Method TP01, "Self-cleaning U-232 Tracer."

