

THORIUM IN WATER

1. SCOPE

- 1.1. This is a method for the separation of thorium from water samples. After completing this method, source preparation for measurement of thorium by alpha spectrometry is performed by electrolytic deposition onto stainless steel planchets (Eichrom Method SPA02) or by rare earth fluoride micro precipitation onto polypropylene filters (Eichrom Method SPA01).
- 1.2. This method does not address all aspects of safety, quality control, calibration or instrument set-up. However, enough detail is given for a trained radiochemist to achieve accurate and precise results for the analysis of the analyte(s) from the appropriate matrix, when incorporating the appropriate agency or laboratory safety, quality and laboratory control standards.

2. SUMMARY OF METHOD

- 2.1. Thorium is separated by Eichrom TEVA resin prior to measurement by alpha spectrometry. A calcium phosphate precipitation is used to concentrate actinides from water samples. Tracers are used to monitor chemical recoveries and correct results to improve precision and accuracy.

3. SIGNIFICANCE OF USE

- 3.1. This is a rapid, reliable method for measurement of actinides in water samples that is more cost-effective and efficient than traditional ion exchange, solvent extraction and precipitation techniques.

4. INTERFERENCES

- 4.1. Nuclides with unresolvable alpha energies such as ^{241}Am and ^{238}Pu , ^{237}Np and ^{234}U , or ^{232}U and ^{210}Po must be chemically separated to enable measurement. This method separates these isotopes effectively.
- 4.2. Very high levels of phosphate in the sample may lead to reduced recovery of actinides in calcium phosphate precipitation and column separations. Adjusting the amount of phosphate added to co-precipitate the actinides may be necessary in these cases.

4.3. The sample preparation procedure outlined in this method will adequately recover actinides from freshly collected, well preserved, homogenous water samples. Older, poorly preserved samples or samples with significant organic or solid matter may require more aggressive treatment to recover actinides which have precipitated or adsorbed to the walls of the storage container or solid matter. Rinsing the empty storage container with warm HNO_3 , adjusting the HNO_3 concentration of the sample to 1M HNO_3 and boiling, and/or wet-washing the calcium phosphate precipitate may be required for older, poorly preserved samples.

5. APPARATUS

- Analytical balance, 0.0001 g sensitivity
- Beakers, glass
- Centrifuge tubes, 50mL and 250mL
- Centrifuge, with rotor and carriers for 50mL and 250mL tubes
- Column rack, Eichrom Part: AC-103
- Extension Funnels, 25 mL, Eichrom Part: AC-120
- Fume hood
- Hotplate
- Stir rods, glass

6. REAGENTS

Note: Analytical grade or ACS grade reagents are recommended. Evaluation of key reagents, such as aluminum nitrate and ammonium hydrogen phosphate, for contribution to method background levels from naturally occurring radioactive materials is recommended.

<i>Aluminum nitrate nonahydrate, $\text{Al}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$</i>
<i>Ammonium hydrogen phosphate, $(\text{NH}_4)_2\text{HPO}_4$</i>
<i>Ammonium hydroxide(57% NH_4OH or 28% NH_3), concentrated NH_4OH</i>
<i>Ammonium oxalate monohydrate, $(\text{NH}_4)_2\text{C}_2\text{O}_4 \cdot \text{H}_2\text{O}$</i>
<i>Appropriate tracers or standards (Th-229)</i>
<i>Calcium nitrate, CaNO_3</i>
<i>Deionized water, All reagents are prepared with deionized water</i>
<i>Hydrochloric acid (37%), concentrated HCl</i>
<i>Isopropanol, $\text{C}_3\text{H}_7\text{OH}$</i>
<i>Nitric acid (70%), concentrated HNO_3</i>

Oxalic acid dihydrate, H₂C₂O₄·2H₂O

Phenolphthalein pH indicator

TEVA® resin, 2mL prepacked column, 100-150µm, Eichrom Part TE-C50-A

- 6.1. *Ammonium hydrogen phosphate (3.2M)* - Dissolve 104g of (NH₄)₂HPO₄ in 200mL of water. Heat gently to dissolve. Dilute to 250mL with water.
- 6.2. *Calcium nitrate (1.25M)* - Dissolve 51g of Ca(NO₃)₂ in 100mL of water and dilute to 250mL with water.
- 6.3. *Hydrochloric acid (9M)* - Add 750mL of concentrated HCl to 100mL of water. Dilute to 1L with water.
- 6.4. *Nitric acid solution (3M)* - Add 188mL of concentrated HNO₃ to 800mL of water. Dilute to 1L with water.
- 6.5. *Nitric acid solution (3M) - aluminum nitrate (1M)* - Add 188mL of concentrated HNO₃ to 500mL of water. Dissolve 375g Al(NO₃)₃·9H₂O. Dilute to 1L with water.
- 6.6. *Phenolphthalein indicator* - Dissolve 1g of phenolphthalein in 50mL of isopropyl alcohol and add 50mL of water.

7. PROCEDURE

7.1. Water Sample Preparation:

- 7.1.1. If samples larger than 1L are analyzed, evaporate the sample to approximately 1L.
- 7.1.2. Aliquot 500 to 1000mL of the filtered sample (or enough to meet required detection limit) into an appropriate size beaker.
- 7.1.3. Add 5mL concentrated HNO₃.
- 7.1.4. Add appropriate tracers and/or reference standards per lab protocol.
- 7.1.5. Calcium phosphate precipitation:
 - 7.1.5.1. Add 0.5mL of 1.25M Ca(NO₃)₂ to each sample.
 - 7.1.5.2. Place each beaker on a hotplate.
 - 7.1.5.3. Cover each beaker with a watch glass.

- 7.1.5.4. Heat at medium setting for 30-60 minutes.
- 7.1.5.5. Once the samples boil, take the watch glass off the beaker and turn the heat down.
- 7.1.5.6. Add 2-3 drops of phenolphthalein indicator and 1mL of 3.2M $(\text{NH}_4)_2\text{HPO}_4$ solution.
- 7.1.5.7. While stirring, slowly add enough concentrated NH_4OH to reach the phenolphthalein end point and form a calcium phosphate precipitate. Heat samples for another 20-30 minutes.
- 7.1.5.8. Remove samples from the hot plate, cool to room temperature, and allow precipitate to settle until solution can be decanted (30 minutes to 2 hours) or centrifuge.
- 7.1.5.9. Decant supernate and discard to waste.
- 7.1.5.10. Transfer the precipitate to a centrifuge tube and centrifuge the precipitate for approximately 10 minutes at 2000 rpm.
- 7.1.5.11. Decant supernate and discard to waste.
- 7.1.5.12. Wash the precipitate with an amount of water approximately twice the volume of the precipitate. Mix well on a vortex mixer. Centrifuge for 5-10 minutes. Discard the supernate.
- 7.1.5.13. If an ammonia odor persists repeat 7.1.5.12.
- 7.1.5.14. Dissolve precipitate in 5mL concentrated HNO_3 . Transfer solution to a 100mL beaker. Rinse centrifuge tube with 2-3mL of concentrated HNO_3 and transfer to beaker. Evaporate solution to dryness.

7.2. Th Separation using TEVA Resin

- 7.2.1. Dissolve each precipitate with 10mL of 3M HNO_3 -1M $\text{Al}(\text{NO}_3)_3$.

Note: An additional 5-10mL may be necessary if the volume of precipitate is large.

Note: Cool all reagent solutions and samples to room temperature.

- 7.2.2. For each sample dissolved, place a TEVA Resin column in the column rack.

- 7.2.3. Place a beaker below each column, remove the bottom plug from each column, press top frit down to top of resin bed, and allow to drain. Attach column reservoirs to each column.
- 7.2.4. Add 5mL of 3M HNO₃ into each column to condition resin and allow to drain.
- 7.2.5. Place a clean, labeled beaker below each column.
- 7.2.6. Transfer each dissolved sample into the appropriate TEVA resin column reservoir. Allow solution to drain.
- 7.2.7. Add 5mL of 3M HNO₃ to rinse to each sample beaker and transfer each solution into the appropriate TEVA Resin column. Allow solution to drain through each column.

Note: Uranium, americium and Np(V) are removed with the load solution and 3M HNO₃ rinses.

- 7.2.8. Add 30mL of 3M HNO₃ into each column. Allow solution to drain through all columns. Discard eluate as waste.
 - 7.2.9. Place a clean, labeled 50mL beaker below each column.
 - 7.2.10. Add 15mL of 9M HCl into each column to strip Th. Allow solution to drain through each TEVA column.
- 7.3. Prepare samples for the determination of thorium by alpha spectrometry using electrodeposition (Eichrom SPA02) or rare earth fluoride micro precipitation (Eichrom SPA01).

8. CALCULATIONS

Calculate the actinide activity as follows:

Calculate tracer yield:

$$\text{Yield} = \frac{(C_s - B_s)}{E_s \times A_s}$$

where:

- C_s = measured actinide tracer, cpm
B_s = background, cpm
E_s = counting efficiency for tracer
A_s = tracer activity, dpm

Note: If any tracer may be present in the sample, a spiked and unspiked sample must be analyzed to determine chemical yield, where:

$$Y = \frac{(\text{spiked sample tracer cpm} - \text{unspiked sample tracer cpm})}{E \times \text{actinide spike activity, dpm}}$$

$$\text{Percent yield} = \text{Yield} \times 100$$

Calculate actinide isotope activity:

$$\text{Sample dpm/L} = \frac{S - B}{E \times V \times Y}$$

where:

- S = sample activity, cpm
- B = background, cpm
- E = counting efficiency = measured cpm/dpm of isotopic standard
- V = sample volume, L
- Y = yield

Conversion of dpm/g to pCi/g: $\text{pCi/L} = (\text{dpm/L})/2.22$

9. PRECISION AND BIAS

- 9.1. *Precision* - A relative standard deviation of 3.5% for Th in the range of 1 pCi/L to 20 pCi/L has been reported.
- 9.2. *Bias* - A mean recovery of 102% ± 3.5% for Th have been reported. Since results are corrected based on spike recovery, no significant bias exists for the method.

10. REFERENCES

- 1) Horwitz, E.P., et al. "Separation and Preconcentration of Uranium from Acidic Media by Extraction Chromatography." *Analytica Chimica Acta*, 266, 25-37 (1992).
- 2) Horwitz, E.P., et al. "Separation and Preconcentration of Actinides from Acidic Media by Extraction Chromatography." *Analytica Chimica Acta*, 281, 361-372 (1993).
- 3) Maxwell, S.L., "Rapid Column Extraction Methods for Actinides and Sr-89/90 in Water Samples," *Journal of Radioanalytical and Nuclear Chemistry*, 267(3), 537-543 (2006).
- 4) ASTM Method D7282-06, "Standard Practice for Set-up, Calibration, and Quality Control of Instruments Used for Radioactivity Measurements."
- 5) ASTM Method D3648-14, "Standard Practices for the Measurement of Radioactivity."

