

# Rapid Determination of Sr in 1-2 Liter Seawater Samples

**Summary of Method** Strontium is separated and concentrated from 1-2L samples of seawater with a calcium phosphate precipitation, enhanced with 200mg of iron. Strontium is separated from matrix impurities and potentially interfering radionuclides in the sample using two stacked 2mL cartridges of Eichrom Sr Resin. Radiostrontium is measured on a low background gas flow proportional counter or liquid scintillation counter. Chemical yield of strontium is determined by gravimetric recovery of native stable strontium in the seawater or by ICP-AES measurement. Average chemical recovery of strontium is  $89 \pm 5\%$  for 1L samples and  $82 \pm 4\%$  for 2L samples. Measured values of  $^{90}\text{Sr}$  agreed to within 1% and 4% of reference values, for 1L and 2L, respectively, with two hour count times. The minimum detectable activity for  $^{90}\text{Sr}$  for 2L samples with a two hour count time is 9.1Bq/L. A single operator can prepare batches of 12-24 samples for measurement of radiostrontium in less than 8 hours.


## Reagents

Sr Resin, 2mL Cartridges (Eichrom SR-R50-S)  
 Nitric Acid (70%)  
 Ammonium Hydroxide (listed as 28%  $\text{NH}_3$  or 56%  $\text{NH}_4\text{OH}$ )  
 Deionized Water  
 Iron Carrier (50mg/mL Fe, as ferric nitrate)  
 3.2M  $(\text{NH}_4)_2\text{HPO}_4$   
 2M  $\text{Al}(\text{NO}_3)_3$   
 $^{90}\text{Sr}$  standard  
 Oxalic acid

## Equipment


Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)  
 Cartridge Reservoir, 20mL (Eichrom AR-200-RV20)  
 Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)  
 Yellow Outer Tips (Eichrom AR-1000-OT)  
 50mL Centrifuge Tubes  
 250-500mL Centrifuge Tubes  
 Centrifuge  
 Cupped Stainless Steel Planchets (~5mL volume)  
 Gas Flow Proportional Counter  
 Hot Plate  
 Analytical Balance  
 Vacuum Pump

**Figure 1. Sample Preparation**




## Figure 2. Strontium Resin Separation (Optional <sup>90</sup>Y Ingrowth)

- (1) Precondition Sr Resin with 10mL 8M HNO<sub>3</sub>.
- (2) Load sample at 1-2mL/min.
- (3) Rinse sample tube with 5mL 8M HNO<sub>3</sub>.
- (4) Add tube rinse to Sr Resin. Elute at 1-2mL/min.
- (5) Rinse Sr Resin sequentially with:
  - 15 mL 8M HNO<sub>3</sub>
  - 10mL 3M HNO<sub>3</sub> - 0.05 oxalic acid
  - 10mL 8M HNO<sub>3</sub>
- (6) Dispose of (1) to (5) as waste.
- (7) Strip Sr with 20mL 0.05M HNO<sub>3</sub> at 1mL/min.



- Gas Flow Proportional Counting\*:
- (8) Evaporate samples to dryness on tared cupped stainless steel planchets.
  - (9) Rinse Sr sample vials with 2mL 0.05M HNO<sub>3</sub>. Transfer vial rinse to planchets. Evaporate to dryness.



- (10) Weigh planchets on an analytical balance to determine gravimetric yield of stable Sr(NO<sub>3</sub>)<sub>2</sub>.

- (11) Measure radiostrontium in samples on low background gas flow proportional counter.

- \* (Options for <sup>89/90</sup>Sr Discrimination)
- (a) Sr fraction from step (7) can be transferred to a liquid scintillation vial. <sup>89</sup>Sr can be measured by Cerenkov counting (no LSC cocktail). <sup>89/90</sup>Sr can then be measured by adding liquid scintillation cocktail.
  - (b) Sr fraction from step (10) can be dissolved in 10mL 8M HNO<sub>3</sub> after >7 days of <sup>90</sup>Y ingrowth. <sup>89/90</sup>Sr can be removed on Sr Resin. <sup>90</sup>Y will elute in Sr Resin load and can be counted by liquid scintillation or gas flow proportional counting.

### Performance of <sup>90</sup>Sr Method for 1L and 2L Seawater Samples

Sample Replicates	Sample Volume, L	<sup>90</sup> Sr, Reference Value (mBq/L)	<sup>90</sup> Sr, Measured Value (mBq/L)	% Bias	Sr carrier % Recovery
11	1	148	150 ± 11	1.2	89 ± 5
4	2	148	154 ± 5	4.2	82 ± 4

2 hour count times

MDA = 9.1 mBq/L for 2L sample

## References

- 1) Sherrod L. Maxwell, Brian K. Culligan, Robin C. Utsey, "Rapid determination of radiostrontium in seawater samples," *J. Radioanal. Nucl. Chem.*, 298(2), 867-875 (2013).