

# Rapid Determination of Actinides in Animal Tissue Samples

**Summary of Method** Uranium, Plutonium, and Americium-Curium are separated and concentrated from up to 200g tissue samples. Samples are digested with aqua regia, wet ashed with HNO<sub>3</sub>-H<sub>2</sub>O<sub>2</sub> and muffled over night at 550°C to destroy organic content. Actinides are separated from matrix impurities and potentially interfering radionuclides in the sample using stacked 2mL cartridges of Eichrom TEVA, TRU and DGA Resin. Actinides are measured via alpha spectrometry following cerium fluoride microprecipitation onto Eichrom Resolve<sup>®</sup> Filters. Average chemical recoveries of Pu for 100-200g samples are 93-101%. Typical americium recoveries are 93-105%. Typical uranium recoveries are 82-96%. A single operator can complete the sample preparation for 12-24 samples, including 16 hours for muffling, in less than 24 hours.


## Reagents

TEVA Resin, 2mL Cartridges (Eichrom TE-R50-S)  
 TRU Resin, 2mL Cartridges (Eichrom TR-R50-S)  
 DGA Resin, 2mL Cartridges (Eichrom DN-R50-S)  
 Hydrofluoric Acid (49%) or Sodium Fluoride  
 Iron Carrier (50mg/mL Fe, as ferric nitrate)  
<sup>242</sup>Pu (or <sup>236</sup>Pu if meas. Np), <sup>243</sup>Am and <sup>232</sup>U tracers  
 Oxalic acid/Ammonium oxalate  
 Nitric Acid (70%)                      Hydrochloric Acid (37%)  
 Hydrogen Peroxide (30%)            Deionized Water  
 Cerium Carrier (1mg/mL)            2M Al(NO<sub>3</sub>)<sub>3</sub>  
 Sodium nitrite                          Sulfamic acid  
 Ascorbic acid                            10% (w:w) TiCl<sub>3</sub>  
 Denatured Ethanol

## Equipment




Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)  
 Cartridge Reservoir, 20mL (Eichrom AR-200-RV20)  
 Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)  
 Yellow Outer Tips (Eichrom AR-1000-OT)  
 Resolve Filters in Funnel (Eichrom RF-DF25-25PP01)  
 Muffle Furnace  
 Hot Plate  
 Analytical Balance  
 600mL Glass Beakers  
 Stainless Steel planchets with adhesive  
 Vacuum Pump  
 Alpha Spectrometry System  
 Heat Lamp

**Figure 1. Sample Preparation**



Adjust valence states of actinides by adding (mix between each addition):  
 0.5mL 1.5M Sulfamic acid, 10uL 50mg/mL Fe carrier, 1.25mL 1M Ascorbic acid, 1mL 3.5M NaNO<sub>2</sub>, and 1.5mL 70% HNO<sub>3</sub>.

**Figure 2. Actinide Separation on TEVA - TRU - DGA\***

<p>(1) Precondition stacked 2mL TEVA, TRU, DGA with 10mL 3M HNO<sub>3</sub>.                  (2) Load sample solution.                  (3) Rinse sample tube with 5mL 6M HNO<sub>3</sub>.** Add tube rinse to cartridges.                  (4) Rinse cartridges with 10mL 3M HNO<sub>3</sub>.                  (5) Separate TEVA, TRU, and DGA cartridges.</p>		<p>(13) Strip Am and Cm from DGA with 10mL 0.25M HCl.                  (14) Rinse TRU cartridge with 15mL 4M HCl-0.2M HF-0.002M TiCl<sub>3</sub>.                  (15) Strip U from TRU with 15mL of 0.1M ammonium bioxalate.                  (16) Add 0.5mL 10% TiCl<sub>3</sub> to U samples, 0.5mL 30% H<sub>2</sub>O<sub>2</sub> to Pu and 0.2mL 30% H<sub>2</sub>O<sub>2</sub> to Am/Cm samples.</p>	<p>(23) Draw vacuum until filter is dry.                  (24) Remove filter from funnel assembly and mount filter on stainless steel planchet with 2-sided tape.</p>
<p>(6) Rinse TEVA cartridge with:                  -10mL 3M HNO<sub>3</sub>                  -20mL 9M HCl (remove Th)                  -5mL 3M HNO<sub>3</sub></p>	<p>(17) Add 50-100ug Ce carrier to all samples. Mix well. Add 1mL 49% HF. Mix well. Wait 15-20 minutes.                  (18) Set up Resolve® Filter Funnel on vacuum box.                  (19) Wet filter with 3mL 80% ethanol followed by 3mL DI water.</p>	<p>(25) Dry filter under heat lamp for 3-5 minutes.                  (26) Measure actinides by alpha spectrometry.</p>	
<p>(7) Strip Pu(Np) from TEVA with 20mL 0.1M HCl-0.05M HF-0.01M TiCl<sub>3</sub>.                  (8) Rinse DGA with 8mL 0.1M HNO<sub>3</sub>.                  (9) Place TRU cartridge above DGA.                  (10) Strip Am/Cm from TRU onto DGA with 15mL 3M HCl.                  (11) Separate TRU and DGA. Set TRU aside for U recovery.                  (12) Rinse DGA with:                  -5mL 3M HCl                  -3mL 1M HNO<sub>3</sub>                  -15mL 0.05M HNO<sub>3</sub></p>	<p>(20) Filter sample.                  (21) Rinse sample tube with 5mL DI water and add to                  (22) Rinse filter funnel with 3mL DI water and 2mL 100% ethanol.</p>	 	

\*Radiostrontium may also be measured by adding a 2mL + 1mL Sr Resin cartridge below DGA and following separation scheme in Eichrom application note AN-1407, "Rapid Determination of Sr in Animal Tissue Samples."

\*\*Adding 50uL of 30% H<sub>2</sub>O<sub>2</sub> to tube rinse can improve U recoveries and decontamination in Pu/Np samples.

**Method Performance for 100-200g Tissue Samples**

Sample	mass, g	replicates	% Tracer Recovery		
			Pu-236	Am-243	U-232
Beef	100	6	98.7 ± 5.7	97.1 ± 8.4	93.4 ± 4.7
Deer	100	59	99.3 ± 12	93.4 ± 10	90.4 ± 8.0
Fish-Bass	200	72	96.2 ± 14	102 ± 13	95.1 ± 8.1
Fish-Bream	100	57	96.6 ± 12	98.4 ± 7.7	91.1 ± 6.3
Fish-Catfish	200	69	98.3 ± 12	103.7 ± 7.6	89 ± 12
Hog	100	17	93 ± 20	96.4 ± 9.7	86 ± 15
Shellfish	100	5	101.3 ± 2.2	97.4 ± 7.1	81.7 ± 3.2

**Reference** Sherrod L. Maxwell, Donald M. Faison, "Rapid column extraction method for actinides and strontium in fish and other animal tissue samples," *J. Radioanal. Nucl. Chem.*, 275(3), 605-612 (2007).