

# Rapid Determination of Actinides in Emergency Urine Samples

**Summary of Method** Uranium, Plutonium, and Americium-Curium are separated and concentrated from 100mL urine samples using calcium phosphate precipitation. The precipitate is dissolved in HNO<sub>3</sub>-H<sub>2</sub>O<sub>2</sub> and wet ashed to destroy residual organic material. The wet-ashed residue is dissolved in nitric acid and aluminum nitrate. Actinides are separated from matrix impurities and potentially interfering radionuclides in the sample using 2mL cartridges of Eichrom TEVA and TRU Resins. Actinides are measured by alpha spectrometry following source preparation by cerium fluoride microprecipitation onto Eichrom Resolve® Filters. Chemical yields are determined by recovery of <sup>232</sup>U, <sup>243</sup>Am, and <sup>242</sup>Pu (or <sup>236</sup>Pu, if measuring <sup>237</sup>Np) tracers. Typical chemical recoveries are >90%. A single operator can complete the separation method for batches of 12-24 samples in as little as 4-5 hours.


## Reagents

TEVA Resin, 2mL Cartridges (Eichrom TE-R50-S)	
TRU Resin, 2mL Cartridges (Eichrom TR-R50-S)	
Ammonium Hydroxide (listed as 28% NH <sub>3</sub> or 56% NH <sub>4</sub> OH)	
<sup>242</sup> Pu (or <sup>236</sup> Pu if meas. Np), <sup>243</sup> Am and <sup>232</sup> U tracers	
Oxalic acid/Ammonium oxalate	
Hydrofluoric Acid (49%) or Sodium Fluoride	
Nitric Acid (70%)	Hydrochloric Acid (37%)
Hydrogen Peroxide (30%)	Deionized Water
Iron Carrier (50mg/mL)	Cerium Carrier (1mg/mL)
1.25M Ca(NO <sub>3</sub> ) <sub>2</sub>	3.2M (NH <sub>4</sub> ) <sub>2</sub> HPO <sub>4</sub>
2M Al(NO <sub>3</sub> ) <sub>3</sub>	10% (w:w) TiCl <sub>3</sub>
Sulfamic Acid	Ascorbic Acid
Sodium Nitrite	Denatured Ethanol





## Equipment

Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)  
 Cartridge Reservoir, 20mL (Eichrom AR-200-RV20)  
 Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)  
 Yellow Outer Tips (Eichrom AR-1000-OT)  
 Resolve Filters in Funnel (Eichrom RF-DF25-25PP01)  
 50mL and 250mL Centrifuge Tubes  
 Centrifuge  
 Hot Plate  
 Analytical Balance  
 250mL Glass Beakers  
 Alpha Spectrometry System  
 Vacuum Pump  
 Heat Lamp

**Figure 1. Sample Preparation**



**Figure 2. Actinide Separation on TEVA - TRU\***

<p>(1) Precondition stacked 2mL TEVA-TRU with 10mL 3M HNO<sub>3</sub>.</p> <p>(2) Load sample solution.</p> <p>(3) Rinse sample tube with 5mL 6M HNO<sub>3</sub>. Add tube rinse to cartridges.**</p> <p>(4) Rinse cartridges with 5mL 3M HNO<sub>3</sub>.</p> <p>(5) Separate TEVA and TRU cartridges.</p>		<p>(11) Add 0.5mL 30% H<sub>2</sub>O<sub>2</sub> to Pu and 0.2mL 30% H<sub>2</sub>O<sub>2</sub> to Am/Cm samples for additional U decon. during CeF<sub>3</sub> ppt.</p> <p>(12) Add 0.5mL of 10% TiCl<sub>3</sub> to each U sample for CeF<sub>3</sub> ppt.</p> <p>(13) Add 50-100ug Ce carrier to all samples. Mix well. Add 1mL 49% HF. Mix well. Wait 15-20 minutes.</p>	<p>(19) Draw vacuum until filter is dry.</p> <p>(20) Remove filter from funnel assembly and mount filter on stainless steel planchet with 2-sided tape.</p>
<p>(6) Rinse TEVA cartridge with: -15mL 3M HNO<sub>3</sub> -20mL 9M HCl (remove Th) -5mL 3M HNO<sub>3</sub></p> <p>(7) Strip Pu(Np) from TEVA with 20mL 0.1M HCl-0.05MHF-0.01M TiCl<sub>3</sub>.</p>	<p>(14) Set up Resolve® Filter Funnel on vacuum box.</p> <p>(15) Wet filter with 3mL 80% ethanol followed by 3mL DI water.</p>		
<p>(8) Strip Am/Cm from TRU with 15mL 4M HCl. Dilute to 30mL prior to CeF<sub>3</sub> ppt.</p> <p>(9) Rinse TRU with 15mL 4M HCl-0.2M HF. (Th removal)</p> <p>(10) Strip U from TRU with 15mL 0.1M ammonium bioxalate.</p>	<p>(16) Filter sample.</p> <p>(17) Rinse sample tube with 5mL DI water and add to filter.</p> <p>(18) Rinse filter funnel with 3mL DI water and 2mL 100% ethanol.</p>	<p>(21) Dry filter under heat lamp for 3-5 minutes.</p> <p>(22) Measure actinides by alpha spectrometry.</p> 	

\*Strontium may also be measured by adding a 2mL Sr Resin Cartridge below DGA and following the separation scheme in Eichrom application note AN-1410, "Rapid Determination of Sr in Emergency Urine Samples."

\*\*Adding 50uL of 30% H<sub>2</sub>O<sub>2</sub> to the 6M HNO<sub>3</sub> tube rinse can help improve U recoveries and decontamination in the Pu/Np fraction.

## References

- 1) Sherrod L. Maxwell, Brian K. Culligan, "Rapid separation method for emergency water and urine samples," *J. Radioanal. Nucl. Chem.*, 279(3), 901-907(2009).