

Rapid Determination of ^{226}Ra in Water Samples

Summary of Method ^{226}Ra is separated from up to 1 liter water samples and measured by alpha spectrometry. Radium is precipitated from samples with calcium carbonate. The calcium carbonate precipitate is dissolved in hydrochloric acid, and cation exchange chromatography is used to purify radium and barium from matrix ions. Eichrom DGA Resin is used to remove other alpha emitting nuclides from radium. Samples are prepared for radium measurement by alpha spectrometry using barium sulfate micro-precipitation method onto Eichrom® Resolve Filters. Sample preparation, including alpha spectrometry source preparation, for batches of 12-24 samples can be completed by a single operator in as little as 3-4 hours. Yields can be traced with ^{133}Ba by gamma spectrometry or ^{225}Ra (^{229}Th) by alpha spectrometry. If tracing with ^{225}Ra , >8 hours of ingrowth time for the alpha emitting ^{217}At daughter of ^{225}Ra is required prior to measurement by alpha spectrometry.

Reagents

Cation Exchange Resin (Eichrom C8-B500-F-H)	
DGA Resin, Normal 2mL Cartridges (Eichrom DN-R50-S)	
Ammonium Hydroxide (listed as 28% NH_3 or 56% NH_4OH)	
^{133}Ba or ^{225}Ra (^{229}Th) Tracer	
Nitric Acid (70%)	Hydrochloric Acid (37%)
1.25M $\text{Ca}(\text{NO}_3)_2$	2M Na_2CO_3
Barium Carrier (1mg/mL)	Isopropyl Alcohol
Ammonium Sulfate	Ascorbic Acid
Denatured Ethanol	Deionized Water
H_2O_2 (30%)	

Equipment

Plastic Chromatography Column (Eichrom AC-50E-5M)
 Column Extension Funnel (Eichrom AC-20X-20M)
 Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)
 Cartridge Reservoir, 20mL (Eichrom AR-200-RV20)
 Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)
 Yellow Outer Tips (Eichrom AR-1000-OT)
 Resolve Filter in Disposable Funnel (Eichrom RF-DF-25-25PP01)
 50mL and 250mL Centrifuge Tubes
 Centrifuge
 Stainless Steel Planchets with adhesive tape
 Hotplate
 Alpha Spectrometry System
 Gamma Spectrometry System (if ^{133}Ba tracer used)
 150mL Glass beakers
 Vacuum Pump
 Heat Lamp

Figure 1. Sample Preparation









Figure 2. Column Purification and Alpha Source Preparation

<p>(1) Prewash 5.0g 50Wx8 200-400 mesh, cation exchange resin column¹: -10mL deionized water -20mL 6M HCl -10mL 0.5M HCl</p> <p>(2) Load Sample²</p> <p>(3) Rinse 30mL 3M HCl</p> <p>(4) Add 2mL DGA cartridge below cation exchange column.</p>		<p>(8) Add 50ug Ba carrier. Mix well.</p> <p>(9) Add 3g (NH₄)₂SO₄ and 5mL iso-propanol. Mix well.</p> <p>(10) Place in ice bath for 30 minutes.</p> <p>(11) Set up Resolve® Filter Funnel on vacuum box.</p> <p>(12) Wet filter with 3mL 80% ethanol followed by 3mL DI water.</p> <p>(13) Filter sample.</p> <p>(14) Rinse sample tube with 5mL DI water and add to filter.</p> <p>(15) Rinse filter funnel with 3mL DI water.</p> <p>(16) Rinse filter funnel with 1-2mL 100% ethanol.</p> <p>(17) Draw vacuum until filter is dry.</p>	<p>(18) Remove filter from funnel assembly and mount filter on stainless steel planchet with adhesive tape.</p>
<p>(5) Strip Ra/Ba with 25mL 5M HNO₃.</p> <p>(6) Add 2mL 30% H₂O₂. Evaporate to dryness.</p> <p>(7) Dissolve residue in 10mL 1.5M HCl.</p>			<p>(19) Dry filter under heat lamp for 3-5 minutes.</p> <p>(20) Measure ²²⁶Ra and ²²⁵Ra(²¹⁷At) by alpha spectrometry after >8 hours ²¹⁷At ingrowth. (¹³³Ba by gamma, if necessary.)</p>
			

¹If using ¹³³Ba tracer, 3.0g of cation exchange resin and proportionally smaller rinse volumes may be used.

²If tracing with ²²⁹Th, a 20mL 1M HCl-1M H₃PO₄ rinse following the sample load can improve purity of final ²²⁶Ra fraction.

Method Performance ²²⁶ Ra in Water				
Sample	²²⁵ Ra(²¹⁷ At) % Yield*	²²⁶ Ra(mBq/L) Reference	²²⁶ Ra(mBq/L) Measured	% Bias
1	84.8	73.8	69.6	-5.7
2	87.3	73.8	75.7	2.6
3	86.2	73.8	71.3	-3.4
4	98.7	73.8	66.9	-9.3
AVG	89 ± 6	73.8	70.9 ± 3.7	-3.9

*²²⁵Ra tracer is added as ²²⁹Th in equilibrium with its daughters and measured by its alpha emitting ²¹⁷At daughter (7.066MeV) after >8 hr ingrowth.

References

1) Sherrod L. Maxwell, Brian K. Culligan, "Rapid Determination of ²²⁶Ra in Environmental Samples," *J. Radioanal. Nucl. Chem.*, 293(1), 149-155 (2012).