

# Rapid Determination of Actinides in Emergency Air Filter Samples


**Summary of Method** U, Pu, Np, Am and Cm are separated and concentrated from air filters. Samples are digested in Teflon beakers once with HNO<sub>3</sub>-H<sub>2</sub>O<sub>2</sub>-HF and then several times with HNO<sub>3</sub>-H<sub>2</sub>O<sub>2</sub>. After evaporating to dryness from HNO<sub>3</sub>-H<sub>3</sub>BO<sub>3</sub> to complex any residual fluoride, actinides are valence adjusted and separated on stacked 2mL cartridges of Eichrom TEVA and TRU resins. Actinides are measured by alpha spectrometry following CeF<sub>3</sub> microprecipitation onto Eichrom Resolve<sup>®</sup> Filters. Chemical yields of tracers averaged from 94±12% for <sup>242</sup>Pu, 87±6% for <sup>243</sup>Am, and 67±32% for <sup>232</sup>U. Poor <sup>232</sup>U recoveries in some samples were traced to insufficient mass of Ce carrier in the source preparation step. Recovery of <sup>232</sup>U improved upon increasing to 100ug of Ce carrier. Measured values typically agreed to within 10% of reference values. Sample preparation for batches of 12 samples can be completed by a single operator in <8 hours. Alpha spectrometry count times will depend on detection limit and data quality objectives.

## Reagents




- TEVA Resin, 2mL Cartridges (Eichrom TE-R50-S)
- TRU Resin, 2mL Cartridges (Eichrom TR-R50-S)
- Iron carrier (50mg/mL Fe, as ferric iron nitrate)
- <sup>242</sup>Pu (or <sup>236</sup>Pu if meas. Np), <sup>243</sup>Am and <sup>232</sup>U tracers
- Oxalic acid/Ammonium oxalate
- Ce carrier (1mg/mL)
- Deionized water            2M Al(NO<sub>3</sub>)<sub>3</sub>
- 10% (w:w) TiCl<sub>3</sub>            HNO<sub>3</sub> (70%)
- HCl (37%)                    HF (49%) or NaF
- Boric acid                    H<sub>2</sub>O<sub>2</sub> (30%)
- NaNO<sub>2</sub>                        Denatured ethanol
- Sulfamic Acid                Ascorbic Acid

## Equipment

- Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)
- Cartridge Reservoir, 20mL (Eichrom AR-200-RV20)
- Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)
- Yellow Outer Tips (Eichrom AR-1000-OT)
- Resolve Filters in Funnel (Eichrom RF-DF25-25PP01)
- 50mL Centrifuge Tubes
- Centrifuge
- Heat Lamp
- Hot Plate
- Analytical Balance
- 250mL Teflon beakers
- Stainless Steel Planchets with adhesive tape
- Alpha Spectrometry System
- Vacuum Pump



**Figure 2. Actinide Separation on TEVA - TRU\***

<p>(1) Precondition stacked 2mL TEVA-TRU with 10mL 3M HNO<sub>3</sub>.</p> <p>(2) Load sample solution.</p> <p>(3) Rinse sample tube with 5mL 6M HNO<sub>3</sub>. Add tube rinse to cartridges.**</p> <p>(4) Rinse cartridges with 5mL 3M HNO<sub>3</sub>.</p> <p>(5) Separate TEVA and TRU cartridges.</p>		<p>(11) Add 0.5mL of 10% TiCl<sub>3</sub> to each U sample for CeF<sub>3</sub> ppt.</p> <p>(12) Add 50-100ug Ce carrier to all samples. Mix well. Add 1mL 49% HF. Mix well. Wait 15-20 minutes.</p> <p>(13) Set up Resolve® Filter Funnel on vacuum box.</p>	<p>(19) Remove filter from funnel assembly and mount filter on stainless steel planchet with 2-sided tape.</p>
<p>(6) Rinse TEVA cartridge with: -15mL 3M HNO<sub>3</sub> -20mL 9M HCl( remove Th) -5mL 3M HNO<sub>3</sub></p> <p>(7) Strip Pu(Np) from TEVA with 20mL 0.1M HCl-0.05MHF-0.01M TiCl<sub>3</sub>. Add 0.5mL 30% H<sub>2</sub>O<sub>2</sub> for additional U decon. during CeF<sub>3</sub></p>	<p>(14) Wet filter with 3mL 80% ethanol followed by 3mL DI water.</p> <p>(15) Filter sample.</p> <p>(16) Rinse sample tube with 5mL DI water and add to filter.</p>		<p>(20) Dry filter under heat lamp for 3-5 minutes.</p> <p>(21) Measure actinides by alpha spectrometry.</p>
<p>(8) Strip Am/Cm from TRU with 15mL 4M HCl. Dilute to 30mL and add 0.2mL 30% H<sub>2</sub>O<sub>2</sub>.</p> <p>(9) Rinse TRU with 15mL 4M HCl-0.2M HF. (Th removal)</p> <p>(10) Strip U from TRU with 15mL 0.1M ammonium bioxalate.</p>	<p>(17) Rinse filter funnel with 3mL DI water and 2mL 100% ethanol.</p> <p>(18) Draw vacuum until filter is dry.</p>		

<sup>89/90</sup>Sr can also be measured by placing a 2mL Sr Resin cartridge below DGA and following the separation scheme in application note AN-1434

\*\*Adding 50uL of 30% H<sub>2</sub>O<sub>2</sub> to the 6M HNO<sub>3</sub> tube rinse can help improve U recoveries and decontamination in the Pu/Np fraction.

## References

1) Sherrod L. Maxwell, Brian K. Culligan, Gary W. Noyes, "Rapid separation method for actinides in emergency air filter samples," *Applied Radiation and Isotopes*, 68(12), 2125-2131 (2010).