

# Rapid Determination of Np/Pu in 20-50g Soil Samples

**Summary of Method** Plutonium and Neptunium are separated and concentrated from 20-50 gram soil samples. Samples are leached with HNO<sub>3</sub> and HCl. The leachates are evaporated to dryness, and sequential precipitations with Fe/Ti-hydroxide and LaF<sub>3</sub> facilitate matrix removal. Pu-Np are separated on 2mL cartridges of Eichrom TEVA resin. Pu-Np are measured by alpha spectrometry following CeF<sub>3</sub> microprecipitation onto Eichrom Resolve® Filters. Chemical yields of the <sup>236</sup>Pu tracer ranged from 82-96%. Measured values typically agreed to within 10% of reference values. Sample preparation for batches of 12 samples can be completed by a single operator in <8 hours.


## Reagents

TEVA Resin, 2mL Cartridges (Eichrom TE-R50-S)	
Iron carrier (50mg/mL Fe, as ferric iron nitrate)	
<sup>236</sup> Pu tracer	NH <sub>4</sub> OH (28% NH <sub>3</sub> or 56% NH <sub>4</sub> OH)
La carrier (10mg/mL)	Ce carrier (1mg/mL)
Deionized Water	2M Al(NO <sub>3</sub> ) <sub>3</sub>
10% (w:w) TiCl <sub>3</sub>	HNO <sub>3</sub> (70%)
HCl (37%)	HF (49%) or NaF
Boric acid	H <sub>2</sub> O <sub>2</sub> (30%)
NaNO <sub>2</sub>	Denatured ethanol
Sulfamic Acid	Ascorbic Acid

## Equipment



Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)  
 Cartridge Reservoir, 20mL (Eichrom AR-200-RV20)  
 Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)  
 Yellow Outer Tips (Eichrom AR-1000-OT)  
 Resolve Filters in Funnel (Eichrom RF-DF25-25PP01)  
 600mL Glass beakers  
 Stainless Steel Planchets with adhesive tape  
 Alpha Spectrometry System  
 50mL and 250mL Centrifuge Tubes  
 Centrifuge  
 Heat Lamp  
 Hot Plate  
 Analytical Balance  
 Vacuum Pump

**Figure 1. Sample Preparation**



**Figure 2. Actinide Separation on TEVA - TRU - DGA and Source Preparation**

<p>Cool samples to room temp. Fix valence by adding: (mix between steps)</p> <ul style="list-style-type: none"> <li>-0.5mL 1.5M sulfamic acid</li> <li>-40uL 50mg/mL Fe carrier</li> <li>-1.5mL 1M ascorbic acid (Wait 3 min)</li> <li>-1mL 3.5M NaNO<sub>2</sub></li> </ul>	<p>(7) 50ug Ce carrier. Mix. Add 1mL 49% HF. Mix. Wait 15-20 minutes.</p> <p>(8) Set up Resolve® Filter Funnel on vacuum box.</p>	<p>(14) Remove filter from funnel assembly and mount filter on stainless steel planchet with 2-sided tape.</p>
<p>(1) Precondition 2mL TEVA, cartridges with 5mL 8M HNO<sub>3</sub>.</p> <p>(2) Load Sample.</p> <p>(3) Rinse centrifuge tube with 5mL 6M HNO<sub>3</sub>.<sup>*</sup> Add to TEVA.</p> <p>(4) Rinse cartridges with:</p> <ul style="list-style-type: none"> <li>-15mL 3M HNO<sub>3</sub></li> <li>-20mL 9M HCl</li> <li>-5mL 3M HNO<sub>3</sub>.</li> </ul> <p>(5) Strip Pu/Np from TEVA with 20mL 0.1M HCl-0.05M HF-0.03M TiCl<sub>3</sub>.</p> <p>(6) Add 0.5mL 30% H<sub>2</sub>O<sub>2</sub> to each sample for additional uranium decontamination during CeF<sub>3</sub> ppt.</p>	<p>(9) Wet filter with 3mL 80% ethanol followed by 3mL DI water.</p> <p>(10) Filter sample.</p> <p>(11) Rinse sample tube with 5mL DI water and add to filter.</p> <p>(12) Rinse filter funnel with 3mL DI water and 2mL 100% ethanol.</p> <p>(13) Draw vacuum until filter is dry.</p>	<p>(15) Dry filter under heat lamp for 3-5 minutes.</p> <p>(16) Measure actinides by alpha spectrometry.</p>



\*Adding 50uL 30% H<sub>2</sub>O<sub>2</sub> to the tube rinse can improve Uranium decontamination.

Sample, g	replicates	<sup>236</sup> Pu Tracer % Yield	<sup>239</sup> Pu			<sup>238</sup> Pu			<sup>237</sup> Np		
			Reference (mBq)	Measured (mBq)	<sup>239</sup> Pu %Bias	Reference (mBq)	Measured (mBq)	<sup>238</sup> Pu %Bias	Reference (mBq)	Measured (mBq)	<sup>237</sup> Np %Bias
20	6	89 <sub>+6</sub>	116.3	118 <sub>+7</sub>	1.5	63.2	67 <sub>+4</sub>	6.0	37.0	39 <sub>+4</sub>	5.4
20	6	96 <sub>+7</sub>	1.69	2.1 <sub>+0.4</sub>	24	25.3	25 <sub>+2</sub>	-1.2	37.0	35 <sub>+2</sub>	-5.4
30	6	82 <sub>+6</sub>	116.3	121 <sub>+5</sub>	4.0	63.2	68 <sub>+5</sub>	7.6	37.0	39 <sub>+4</sub>	5.4
50	6	88 <sub>+5</sub>	116.3	114 <sub>+3</sub>	-2.0	63.2	64 <sub>+2</sub>	1.3	37.0	21 <sub>+11</sub>	-43

**References**

1) Sherrod L. Maxwell, Brian K. Culligan, Gary W. Noyes, "Rapid separation method for <sup>237</sup>Np and Pu isotopes in large soil samples," *Applied Radiation and Isotopes*, 69(7), 917-925 (2011).