

**Summary of Method** A method for the preparation of  $^{225}\text{Ac}$  ( $t_{1/2} = 10$  days) and  $^{225}\text{Ra}$  ( $t_{1/2} = 14.8$  days) from  $^{229}\text{Th}$  ( $t_{1/2} = 7340$  years) source material is presented. The method employs 2mL cartridges of UTEVA and DGA resins to obtain high purity  $^{225}\text{Ac}$  in small volumes of eluate while preserving valuable  $^{229}\text{Th}$  source material. The method is meant for  $^{225}\text{Ac}$  tracer production from  $^{229}\text{Th}$  containing 5-10mg or less of total Th. For separations from larger masses of Th see the Eichrom website bibliography for other options (Recent Advances in the Recovery and Purification of Actinium Isotopes, Horwitz and McAlister, National Meeting of the American Chemical Society, 2009). The source material, containing  $^{229}\text{Th}$ ,  $^{225}\text{Ac}$ ,  $^{225}\text{Ra}$  and other daughter nuclides in 4M  $\text{HNO}_3$ , is loaded onto stacked 2mL cartridges of UTEVA and DGA resins.  $^{229}\text{Th}$  is retained on UTEVA, while  $^{225}\text{Ac}$  is retained on DGA and  $^{225}\text{Ra}$  passes through both cartridges.  $^{225}\text{Ra}$  can be saved for use as a radiotracer or as an additional source of  $^{225}\text{Ac}$ , following a suitable ingrowth period.  $^{225}\text{Ac}$  is recovered from DGA with a small volume of 0.1-2.0M HCl. The  $^{229}\text{Th}$  source is recovered from UTEVA with a small volume of 0.5M HCl. Following a suitable ingrowth period, the  $^{229}\text{Th}$  can be acidified to 4M  $\text{HNO}_3$  and used to produce additional  $^{225}\text{Ac}$  and  $^{225}\text{Ra}$ . The  $^{229}\text{Th}$  is preserved nearly indefinitely and continuously purified from chemical and radiologic impurities run to run.

## Reagents

UTEVA Resin Cartridges (Eichrom UT-R50-S)

DGA Resin Cartridges (Eichrom DN-R50-S)

$^{229}\text{Th}$  Source

Deionized Water

HCl

$\text{HNO}_3$



## Equipment

Glass vials for storage of  $^{229}\text{Th}$  source.

Glass or plastic vials/bottles for collection of  $^{225}\text{Ac}$ .  $^{225}\text{Ra}$  and waste.

5, 10 or 20mL plastic luer lock syringes

Gamma spectrometry system and/or alpha spectrometry for measurement of  $^{225}\text{Ac}$  ( $^{221}\text{Fr}$ ),  $^{225}\text{Ra}$  and  $^{229}\text{Th}$ .



## $^{225}\text{Ac}/^{225}\text{Ra}/^{229}\text{Th}$ Separation

(1) Precondition stacked 2mL cartridges of UTEVA and DGA with 10mL 4M  $\text{HNO}_3$ .

(2) Acidify  $^{229}\text{Th}$  eluate from previous separation with 5mL  $\text{HNO}_3$ . (If new  $^{229}\text{Th}$  source, dilute to 20mL with 4M  $\text{HNO}_3$ .)

(3) Load  $^{229}\text{Th}$  and daughters in 20mL 4M  $\text{HNO}_3$ . Collect and save eluate containing  $^{225}\text{Ra}$ .\*


(4) Rinse UTEVA/DGA with 10mL 4M  $\text{HNO}_3$ . Collect  $^{225}\text{Ra}$ .\*

(5) Separate UTEVA and DGA cartridges.




(6) Rinse DGA with 10mL 8M  $\text{HCl}$ .



(7) Strip  $^{225}\text{Ac}$  with 10mL 0.1M  $\text{HCl}$ . (For higher  $^{229}\text{Th}$  decontamination, use 10mL of 2M  $\text{HCl}$ .)



(8) Strip  $^{229}\text{Th}$  from UTEVA cartridge with 15mL 0.5M  $\text{HCl}$ . (More complete recovery may be obtained by stripping in opposite direction of load.) Save  $^{229}\text{Th}$  for future use.



\* $^{225}\text{Ra}$  can be used directly as a tracer or as a source of additional  $^{225}\text{Ra}$ .



## References

1) McAlister and Horwitz, "Chromatographic Generator Systems for the actinides and natural decay series elements," *Radiochimica Acta*, 99:1-9 (2011).