

RADIUM-228 IN WATER

(WITH VACUUM BOX SYSTEM)

1. SCOPE

- 1.1. This is a method for the separation and measurement of ^{228}Ra in water via its beta emitting ^{228}Ac daughter.
- 1.2. This method does not address all aspects of safety, quality control, calibration or instrument set-up. However, enough detail is given for a trained radiochemist to achieve accurate and precise results for the analysis of the analyte(s) from the appropriate matrix, when incorporating the appropriate agency or laboratory safety, quality and laboratory control standards.

2. SUMMARY OF METHOD

- 2.1. A barium sulfate precipitation technique is used to concentrate radium from water samples. ^{133}Ba is used to monitor chemical recovery and correct results to improve precision and accuracy. After ingrowth, ^{228}Ac , a daughter of ^{228}Ra , is separated by Ln Resin, prior to measurement by low background gas flow proportional counter.

3. SIGNIFICANCE OF USE

- 3.1. This is a rapid, reliable method for measurement of ^{228}Ra in water samples that is based on chemistry similar to EPA procedure RA-05.

4. INTERFERENCES

- 4.1. Potential beta emitters such as bismuth, yttrium and thorium would be retained on the Ln Resin, while eluting actinium from the resin.
- 4.2. Interferences from other radioactive rare earth elements are eliminated under the stripping conditions of 0.35M HNO_3 .

5. APPARATUS

- Analytical balance- 0.0001 g sensitivity
- Beakers, glass

- Bunsen burner
- Cartridge reservoirs, 10mL (Eichrom Part: AR-200-RV10) or 20mL (Eichrom Part: AR-200-RV20)
- Centrifuge, with rotor and carriers for 15mL tubes
- Centrifuge tubes, 50mL
- Fume hood
- Gamma spectrometry system (for determination of ^{133}Ba recovery)
- Hotplate
- Liquid scintillation vials
- Low background gas flow proportional counter with appropriate sample carriers and planchets.
- Stir rods, glass
- Test tubes, Pyrex glass, 15mL
- Vacuum box system, Eichrom Part: AR-12-BOX or AR-24-BOX
- Vacuum box white inner support tube-PE- Eichrom Part: AR-1000-TUBE-PE
- Vacuum box yellow outer tips- Eichrom Part: AR-1000-OT
- Vacuum pump, 115 V, 60 Hz Fisher Part: 01-092-25 (or equivalent) or house vacuum
- Vortex mixer
- Watch glasses
- Optional: vacuum box inner liner - For collection of load and rinse fractions, Eichrom Part: AR-12-LINER or AR-24-LINER

6. REAGENTS

Note: Analytical grade or ACS grade reagents are recommended.

^{133}Ba tracer (~3000 dpm/mL)
Barium chloride dihydrate, $\text{BaCl}_2 \cdot 2\text{H}_2\text{O}$
Deionized water, all reagents are prepared with deionized water
Hydrogen peroxide (30%), concentrated H_2O_2
$\text{LN}^{\text{®}}$ resin, 2mL prepacked cartridge, 50-100 μm , Eichrom Part LN-R50-S
Nitric acid (70%), concentrated HNO_3
Potassium carbonate, K_2CO_3
Sulfuric acid (96%), concentrated H_2SO_4

- 6.1. *Barium carrier (30 mg/ml)*- Dissolve 13.3g $\text{BaCl}_2 \cdot 2\text{H}_2\text{O}$ in 200mL deionized water. Dilute 250mL with water.

- 6.2. *Nitric acid solution (0.095 M)*- Add 5.9mL of concentrated HNO₃ to 900mL of water. Dilute to 1L with water. Check that pH of solution is ~1.
- 6.3. *Nitric acid solution (0.35 M)*- Add 21.9mL of concentrated HNO₃ to 800mL of water. Dilute to 1L with water.
- 6.4. *Potassium carbonate (50wt%)*- Dissolve 50g K₂CO₃ in 50mL of water.

7. PROCEDURE

- 7.1. Water Sample Preparation:
 - 7.1.1. If required, filter the sample through a 0.45 micron filter.
 - 7.1.2. Aliquot 500 to 1000mL of the sample (or enough to meet required detection limit) into an appropriate size beaker. Add ¹³³Ba tracer and 1mL of barium chloride carrier.
 - 7.1.3. Add 10mL of concentrated H₂SO₄ to sample. Place each beaker on a hot plate.
 - 7.1.4. Cover each beaker with a watch glass. Heat samples at about 90°C for 2 hours.
 - 7.1.5. Allow samples to cool. A fine white precipitate is formed. Allow precipitate to settle until solution can be decanted or centrifuge.
 - 7.1.6. Decant supernate and discard to waste.
 - 7.1.7. Transfer the precipitate to a 15mL Pyrex glass centrifuge tube using water. Centrifuge for 10 minutes at 2000 rpm.
 - 7.1.8. Decant supernate and discard to waste.
 - 7.1.9. Wash the precipitate with 4-5 mL of water. Mix well and centrifuge for 5-10 minutes. Check the pH of the supernate. The pH should be about 6, if necessary decant supernate and repeat wash/centrifuge steps.
- 7.2. Conversion of Ba(Ra)SO₄ to Ba(Ra)CO₃:
 - 7.2.1. Add 1 mL of K₂CO₃ solution and 2-3 mL of water, submerging the precipitate completely.
 - 7.2.2. Heat the bottom of the Pyrex tube over a low flame while stirring the tube periodically. Continue heating, mixing until

volume is reduced to ~1mL. (This step may also be performed using microwave heating or a water bath.)

- 7.2.3. Let the tube cool and repeats step 7.2.1 and 7.2.2 two times.
- 7.2.4. Add 15-20 mL of water. Mix well. Centrifuge and check the pH of the supernate. The pH of at this point should be approximately 12.
- 7.2.5. Discard the supernate and continue washing with water and centrifuging until a pH of 7 is achieved. Discard supernate.

Note: The precipitate now should be primarily BaCO₃. The water washing of the BaCO₃ precipitate removes excess CO₃²⁻ and SO₄²⁻ ions. Any sulfate remaining in the sample can result in incomplete dissolution during the following steps and decreased Ba/Ra yields.

- 7.2.6. Once the pH of 7 is achieved add 5 mL of 0.095M HNO₃ to the tube and gently heat as necessary to dissolve the residue in the tube.

Note: If the entire residue does not dissolve in the load solution, then centrifuge the solution and check the supernatant for ¹³³Ba yield recovery. If acceptable (>80%) recovery is produced then proceed further if not then take the residue (undissolved) and repeat steps 7.2.1 through 7.2.5. Combine the dissolved solutions and check for Ba-133 recovery.

- 7.2.7. Transfer the dissolved solution into a 20mL liquid scintillation vial. Rinse the Pyrex tube with an additional 5 mL of 0.095M HNO₃ and add the tube rinse to the same vial.
- 7.2.8. Count the ¹³³Ba by gamma spectrometry. Record ¹³³Ba recovery.
- 7.2.9. Let the solution sit in the vial for at least 30 hours for ingrowth of ²²⁸Ac from ²²⁸Ra.

7.3. ²²⁸Ac Separation Using Ln Resin:

Note: ²²⁸Ac has a 6.13 hour half-life. The following steps should be performed quickly to achieve the lowest detection limits.

- 7.3.1. Place the inner tube rack (supplied with vacuum box system) into the vacuum box with the centrifuge tubes in the rack. Fit the lid to the vacuum system box. Alternatively, a vacuum box inner liner may be used.

- 7.3.2. Place yellow outer tips into all 12 or 24 openings in the lid of the vacuum box. Fit a white inner support tube into each yellow tip.
- 7.3.3. For each sample solution, fit an LN Resin cartridge on to the inner support tube.
- 7.3.4. Add syringe barrels (funnels/reservoirs) to the top end of each LN Resin cartridge.
- 7.3.5. Connect the vacuum pump to the box. Turn the vacuum pump on and ensure proper fitting of the lid.

Note: The unused openings on the vacuum box should be sealed. Vacuum manifold plugs can be used to plug unused white tips to achieve good seal during the separation. Alternatively, unused vacuum box holes can be sealed with scotch tape.

- 7.3.6. Add 5mL of 0.095M HNO₃ into each LN Resin cartridge reservoir to precondition resin. Adjust vacuum to achieve a flow rate of 1mL/min.
- 7.3.7. Transfer each sample solution from step 7.2.7. into the appropriate LN Resin cartridge reservoir.
- 7.3.8. Allow the load solution to completely pass through the cartridge at 1 mL/min.
- 7.3.9. Add 5mL of 0.095M HNO₃ into the sample vial and transfer this rinse to the appropriate LN Resin cartridge reservoir. Allow the solution to completely pass through each cartridge at 1mL/min. Record the date and time at the start of this rinse for decay correction.
- 7.3.10. Place clean reservoirs above each LN Resin cartridge.
- 7.3.11. Rinse each LN Resin cartridge with 5mL of 0.095M HNO₃. Allow solution to completely pass through each cartridge at 1mL/min.
- 7.3.12. Repeat step 7.3.10.
- 7.3.13. Disengage vacuum. Dispose of eluate to this point as waste (or retain in case issues arise with subsequent analysis of ²²⁸Ra via ²²⁸Ac).
- 7.3.14. Place clean, labeled centrifuge tubes below each LN Resin cartridge. Replacing yellow outer tubes and inner tube supports at this point will help ensure clean ²²⁸Ac fractions in the following steps.

7.3.15. Add 10mL of 0.35M HNO₃ to each LN Resin cartridge reservoir to eluate ²²⁸Ac. Engage vacuum. Allow solution to completely pass through each cartridge at 1mL/min.

7.3.16. Prepare sources for counting by rare earth micro precipitation using Eichrom Method SPA-01.

8. CALCULATIONS

Calculate ²²⁸Ra in pCi/L:

$$^{228}\text{Ra (pCi/L)} = \frac{A}{2.22 \times E \times V \times Y \times e^{-\lambda t_1}} \times \frac{\lambda t_2}{1 - e^{-\lambda t_2}}$$

where:

- A = net count rate, cpm
- E = counting efficiency expressed as fraction
- Y = ¹³³Ba (Ra) yield expressed as fraction
- V = Sample volume (liters)
- t₁ = decay time of ²²⁸Ac, from start of rinse until start of counting (minutes)
- t₂ = counting time (minutes)
- λ = decay constant of ²²⁸Ac (1.88*10⁻³ min⁻¹)

9. PERFORMANCE DATA

Reported by Burnett, et al. (1995)

Sample	Expected Result Ra-228 (pCi/L)	Experimental Result Ra-228 (pCi/L)
EMSL Ra in Water, 7/17/92	16.7 ± 3.3	17.6 ± 0.7
EMSL Standard Ra-228 Solution	59.3 ± 2.1	58.0 ± 1.3
Decontamination Experiment ^a	29.7 ± 1.0	31.7 ± 0.5
EMSL Performance Evaluation "A" 4/19/94 ^b	20.1 ± 5.0	21.6 ± 1.2

- a. The decontamination experiment consisted of using approximately 800 mL of EMSL Performance Evaluation "B" (April 20, 1993) which contained the following nuclides: ⁹⁰Sr(⁹⁰Y) = 397 dpm; ⁶⁰Co = 466 dpm; ¹³⁴Cs = 438 dpm; ¹³⁷Cs = 438 dpm; and ²²⁶Ra 213 dpm. In addition the following were added: approximately 8000 dpm ¹³³Ba, 900 dpm ²⁰⁷Bi, 500 dpm ²¹⁰Pb, and 65.9 dpm ²²⁸Ra.
- b. Ra-226 was analyzed via radon emanation by collection of the sample load and rinse from the Ln Resin column. Their result of 20.1 +/- 1.1 pCi/L compared very well to the expected result of 20.0 +/- 3.0 pCi/L.

10. REFERENCES

- 1) Burnett, W.C., P.H. Cable, and Russ Moser, "Determination of Radium-228 in Natural Waters Using Extraction Chromatographic Resins," *Radioactivity & Radiochemistry*, Vol. 6, No. 3, pp. 36-43 (1995).
- 2) Maxwell, S.L., et al., "Rapid method for determination of 228Ra in water samples," *J. Radioanal. Nucl. Chem.* 295, 2181-2188, (2013).

