


Summary of Method A method for the preparation of ²¹⁰Po ($t_{1/2} = 138.4$ days) from ²¹⁰Pb ($t_{1/2} = 22.26$ years) source material is presented. The method employs 2mL cartridges of DGA and Sr resins to obtain high purity ²¹⁰Po in small volumes of eluate while preserving valuable ²¹⁰Pb source material. The source material, containing ²¹⁰Pb/²¹⁰Bi/²¹⁰Po in 2.67M HCl, is loaded onto stacked 2mL cartridges of DGA and Sr resins. ²¹⁰Po and ²¹⁰Bi are retained on DGA Resin, while ²¹⁰Pb is retained on Sr Resin. The ²¹⁰Pb source is recovered from Sr Resin with a small volume of 8M HCl. Following a suitable ingrowth period, the ²¹⁰Pb can be diluted to 2.67M HCl and used to produce additional ²¹⁰Po. The ²¹⁰Pb is preserved nearly indefinitely and continuously purified from chemical and radiologic impurities run to run. ²¹⁰Po is recovered from DGA resin with 0.05M HNO₃, but should be acidified to 1M HNO₃ to prevent loss of Po to glass vials. The ²¹⁰Bi will remain on the DGA resin during the Po elution, and can be recovered with 10mL of 0.05M ammonium bioxalate. The DGA/Sr Resin chemistry is an improvement over the UTEVA/Sr Resin chemistry previously described (AN-1616a), which required 6M HNO₃ to recover the ²¹⁰Po.

Reagents

Sr Resin Cartridges (Eichrom SR-R50-S)
 DGA, Normal Cartridges (Eichrom DN-R50-S)
 Liquid Scintillation Cocktail
²¹⁰Pb Source
 Deionized Water
 HCl
 HNO₃

Equipment

Glass vials for storage of ²¹⁰Pb source.
 Glass or plastic vials/bottles for collection of ²¹⁰Po, ²¹⁰Bi and waste.
 10, 20 or 30mL plastic luer lock syringes
 Liquid Scintillation System for measurement of ²¹⁰Bi and ²¹⁰Po.
 Gamma Spectrometry System for measurement of ²¹⁰Pb.




$^{210}\text{Pb}/^{210}\text{Bi}/^{210}\text{Po}$ Separation

- (1) Precondition stacked 2mL cartridges of DGA and Sr Resins with 10mL 2M HCl.
- (2) Dilute ^{210}Pb eluate from previous separation with 20mL DI H_2O . (If new ^{210}Pb source, dilute to 20mL with 2M HCl.)*
- (3) Load $^{210}\text{Pb}/^{210}\text{Bi}/^{210}\text{Po}$ in 30mL 2.67M HCl. (20mL 2M HCl).
- (4) Rinse DGA/Sr with 10mL 2M HCl.
- (5) Elute 10mL 8M HCl through DGA/Sr, collecting ^{210}Pb Source material. Save ^{210}Pb for future use.





- (6) Separate DGA/Sr.
- (7) Rinse DGA with 10mL 8M HNO_3 . Discard as waste.
- (8) Strip ^{210}Po from DGA with 15mL 0.05M HNO_3 .*



*Acidification of the ^{210}Pb to 1M HNO_3 is recommended to prevent hydrolysis and loss of the ^{210}Po to the storage vial.

*Adding 1mg of stable Pb to the ^{210}Pb source can help improve ^{210}Pb recovery from Sr Resin (do only once, not each time).



References

- 1) McAlister and Horwitz, "Chromatographic Generator Systems for the actinides and natural decay series elements," *Radiochimica Acta*, 99:1-9 (2011).