

Options for Separation and Measurement of $^{89}\text{Sr}/^{90}\text{Sr}$

Daniel McAlister, Ph.D.

Eichrom Technologies, LLC

Eichrom Workshop, May 22, 2018



*The 63rd Conference on
Radiobioassay & Radiochemical
Measurements*

Portland, Maine USA

eichrom[®]
TECHNOLOGIES

A GCI COMPANY

npo

eichrom

**Yield
tracers
analytes**

Properties of Selected Nuclides								
Nuclide	Half-Life	Decay		Detector Suitable for Measurement				
		Mode	Energy	GFPC	LSC	Cerenkov	MS/AES	Gamma
^{82}Sr	25.35 days	ε						
^{82}Rb	1.25 min	β^+	β^+ mean = 1479 keV γ = 511 keV (190.4%)					
^{85}Sr	64.849 days	ε/γ	γ = 514 keV (96%)	No	Yes	No	No	Yes
^{88}Sr	Stable			No	No	No	Yes	No
^{89}Sr	50.563 days	β^-	β_{max} = 1500 keV β_{mean} = 587 keV	Yes	Yes	Yes	No	No
^{90}Sr	28.79 years	β^-	β_{max} = 546 keV β_{mean} = 196 keV	Yes	Yes	No	Yes	No
^{90}Y	64 hours	β^-	β_{max} = 2280 keV β_{mean} = 934 keV	Yes	Yes	Yes	No	No

Methods and Application Notes

Sample Extraction / Dissolution

Matrix removal / Concentration (ppt, IX)

Sr Resin separation to isolate pure Sr fraction

Counting source preparation

Detailed instructions for $^{89}\text{Sr}/^{90}\text{Sr}$ discrimination



eichrom®

AN-1624-10

Options for ^{89}Sr and ^{90}Sr
Determination

Factors

Quality Objectives

Ratio $^{89}\text{Sr}/^{90}\text{Sr}$

Urgency


Budget

Available Equipment




Disclaimer

Uptake on Sr Resin




Load from
8M HNO_3
to reject
Ba, Na, K.



Uptake on Sr Resin

Rinse with
3M HNO₃ –
0.05M
Oxalic acid



Pb remains
on column
when
stripping Sr
with dilute
HNO₃.


Po
removed
with 8M
HNO₃.



Why is $^{89}\text{Sr}/^{90}\text{Sr}$ Challenging?

Decay of $^{89}\text{Sr}/^{90}\text{Sr}$ + Ingrowth ^{90}Y




Pure Beta emitters






^{90}Y decays into stable ^{90}Zr



^{89}Sr decays into stable ^{89}Y

Beta/Positron Decay



https://en.wikipedia.org/wiki/Beta_decay

Beta emission accompanied by electron antineutrino.


Energy distributed between electron and neutrino.

Neutrinos have very little interaction with matter (only weak force).

Measurement yields spectrum of electron energies with characteristic β_{max}^- .

Measurement Techniques

LSC Spectra



IAEA/AQ/27, "Rapid Simultaneous Determination of ^{89}Sr and ^{90}Sr in Milk: a Procedure Using Cherenkov and Scintillation Counting," IAEA Analytical Quality in Nuclear Applications Series No. 27.

LSC Measurement

Pros:

4 π geometry
High Efficiency

Cons:

Requires Cocktail
Quenching
Luminescence
Little Selectivity
Destructive

Cherenkov Radiation

The speed of light in a vacuum is a universal constant (c).

The speed at which light propagates in water is only $\sim 0.75c$. (function of refractive index)

Matter may be accelerated beyond this speed ($>0.75c$, $<c$) during nuclear reactions.

Cherenkov radiation results when a charged particle, most commonly an **electron**, travels through a dielectric medium with a speed greater than that at which light propagates in the same medium.




Cherenkov Spectra

Pros:

No cocktail = No Quench
No luminescence
Non-destructive
Selective for high energy β^-/β^+ .
(Measure $^{89}\text{Sr}/^{90}\text{Y}$, reject ^{90}Sr)

Cons:

Lower Efficiency than LSC



IAEA/AQ/27, "Rapid Simultaneous Determination of ^{89}Sr and ^{90}Sr in Milk: a Procedure Using Cherenkov and Scintillation Counting,"
IAEA Analytical Quality in Nuclear Applications Series No. 27.

Gas Flow Proportional Counting

Pros:

Low background
Simultaneous counting
Non-destructive



Cons:

Efficiency (2π geometry)
Low selectivity (reject alpha)



ICP-MS, ICP-AES, MP-AES , AA

- Screen samples for native Sr content.
- Yield tracer by stable Sr.
- ^{90}Sr (MS), relatively high levels.



Method Options

Following Matrix Removal/Sr Resin Isolation of Sr:

- 1) MS (^{90}Sr Only)

- 2) Two count methods
 - Different counting techniques
 - Same technique with period of ingrowth (^{90}Y)


- 3) Count, ingrowth, Separate ^{90}Y

Mass Spectrometry

- ^{90}Sr Only, or ^{88}Sr yield monitor/native content.
- Separation: Concentrate, remove matrix + ^{90}Zr isobar.
- Limited to higher activity samples (half-life 28.79 yrs).
- Couple with radiometric detection for ^{89}Sr ???

Feuerstein J, Boulyga S.F., Galler P., Stingeder G., Prohaska T., "Determination of ^{90}Sr in soil samples using inductively coupled plasma mass spectrometry equipped with dynamic reaction cell (ICP-DRC-MS)." *J. Environ. Radioact.* 11, 1764-9 (2008).

Cherenkov/LSC



Count ^{89}Sr immediately by Cherenkov
(Limit impact of ^{90}Y ingrowth)

Add Cocktail.

Count $^{89/90}\text{Sr}$ by LSC


1 day, Bias: $<20\%$ ^{89}Sr and $<40\%$ ^{90}Sr

Ratios of $^{89}\text{Sr}/^{90}\text{Sr}$ activity of up to 35.

Fastest route to ^{89}Sr and ^{90}Sr .

Separate Sr

2-Count Methods



Prepare Gas Flow Planchet or LSC source.

Count immediately for total $^{89/90}\text{Sr}$
(Limit impact of ^{90}Y ingrowth)


Wait 7-14 days for ^{90}Y Ingrowth
(Maximum ^{90}Y Ingrowth/Minimum ^{89}Sr decay)

Second Count. Solve equations to calculate ^{89}Sr
and ^{90}Sr . (Appendix B)

Economical. Verify ^{90}Sr by ^{90}Y . Difficult to
measure low ratios of ^{89}Sr .

Separate Sr

Cherenkov/Cherenkov



Count ^{89}Sr immediately by Cherenkov.
(Limit impact of ^{90}Y ingrowth)

Wait 7-14 days for ^{90}Y Ingrowth.

Acidify sample to 3-8 M HNO_3 .

Separate ^{90}Y . Count by Cherenkov.


Rugged. Best $^{89}\text{Sr}/^{90}\text{Sr}$ discrimination.

Minimizes waste (no cocktail).

Separate Sr

Banavali, A. D. et al. "Strontium-89, 90 Analysis by Eichrom Column Chemistry and Cherenkov Counting". 38th Annual Conference on Bioassay Analytical and Environmental Radiochemistry. Santa Fe, NM. November 1992.

GFPC/Ingrowth-Separation/GFPC



Prepare Gas Flow Planchet.

Count immediately for total $^{89/90}\text{Sr}$.
(Limit impact of ^{90}Y ingrowth)

Wait 7-14 days for ^{90}Y Ingrowth

Dissolve sample with 3-8 M HNO_3 .




Separate ^{90}Y . Prepare planchet and count ^{90}Y .

Rugged. Good $^{89}\text{Sr}/^{90}\text{Sr}$ discrimination.
Often lowest background method.

ASTM D5811-08, "Standard Test Method for Strontium-90 in Water."


Separate Sr

Uptake on Sr Resin




Uptake on DGA Resin, Normal

k' on DGA Resin vs HCl
50-100 μm , 2 h, 21(1) $^\circ\text{C}$




DGA Resin, Normal

50-100 μm , 1 h Contact Time, 25(2) $^\circ\text{C}$



3-8M HNO_3
5 mg Stable Sr
 $^{89/90}\text{Sr} + ^{90}\text{Y}$



Rinse, Rinse, Rinse!

Conclusions

- There are many options for the discrimination of ^{89}Sr and ^{90}Sr
- $^{89}\text{Sr}/^{90}\text{Sr}$ discrimination may be achieved by selective measurement techniques, ingrowth of ^{90}Y , and/or separation of ingrown ^{90}Y .
- The method choice must take into account factors such as:
 - Quality Objectives
 - Ratio $^{89}\text{Sr}/^{90}\text{Sr}$
 - Urgency
 - Budget
 - Available Equipment

Questions???