

Rapid Determination of Actinides in Vegetation Samples

Summary of Method U, Pu, Am and Cm are separated and concentrated from 5-10 gram vegetation samples. Samples are muffled in zirconium crucibles 2-4 hours to destroy organic content. The residue is wet ashed with HNO_3 - H_2O_2 and then fused with 15g NaOH at 600°C for ten minutes. The fusion cakes are dissolved in water, transferred to 250mL centrifuge tubes and precipitated twice to facilitate matrix removal. Actinides are separated on stacked 2mL cartridges of Eichrom TEVA, TRU and DGA resins. Actinides are measured by alpha spectrometry following CeF_3 microprecipitation onto Eichrom Resolve® Filters. Chemical yields of tracers ranged from 90-101% for ^{242}Pu , 84-93% for ^{243}Am , and 81-87% for ^{232}U . Measured values agreed to within 1-3% of reference values for Pu isotopes, 3-9% for Am and Cm isotopes, and 2-15% for U isotopes for 16 hour count times. A single operator can prepare batches of 12 samples for the measurement of actinides in less than 8 hours.

Reagents

TEVA Resin, 2mL Cartridges (Eichrom TE-R50-S)
 TRU Resin, 2mL Cartridges (Eichrom TR-R50-S)
 DGA Resin, 2mL Cartridges (Eichrom DN-R50-S)
 Iron Carrier (50mg/mL Fe, as ferric nitrate)
 Lanthanum and Cerium Carriers (10mg/mL)
 ^{242}Pu (or ^{236}Pu if meas. Np), ^{243}Am and ^{232}U tracers
 Oxalic acid/Ammonium oxalate
 Hydrofluoric Acid (49%) or Sodium Fluoride
 3.2M $(\text{NH}_4)_2\text{HPO}_4$ 2M $\text{Al}(\text{NO}_3)_3$
 10% (w:w) TiCl_3 Boric acid
 Sodium Hydroxide Sodium Nitrite
 Denature Ethanol Sulfamic Acid
 Ascorbic Acid Hydrogen Peroxide (30%)
 Nitric Acid (70%) Hydrochloric Acid (37%)
 Deionized Water 1.25M $\text{Ca}(\text{NO}_3)_2$

Equipment

Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)
 Cartridge Reservoir, 20mL (Eichrom AR-200-RV20)
 Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)
 Yellow Outer Tips (Eichrom AR-1000-OT)
 Resolve Filters in Funnel (Eichrom RF-DF25-25PP01)
 50mL and 250mL Centrifuge Tubes
 250mL Zirconium crucibles with zirconium lids
 Alpha Spectrometry System
 Centrifuge Muffle Furnace
 Hot Plate Heat Lamp
 Analytical Balance Vacuum Pump

Figure 1. Sample Preparation







Figure 2. Actinide Separation on TEVA - TRU - DGA*

<p>(1) Precondition stacked 2mL TEVA, TRU, DGA with 10mL 3M HNO₃.</p> <p>(2) Load sample solution.</p> <p>(3) Rinse sample tube with 5mL 6M HNO₃. ** Add tube rinse to cartridges.</p> <p>(4) Rinse cartridges with 10mL 3M HNO₃.</p> <p>(5) Separate TEVA, TRU, and DGA cartridges.</p>		<p>(13) Strip Am and Cm from DGA with 10mL 0.25M HCl.</p> <p>(14) Rinse TRU cartridge with 15mL 4M HCl-0.2M HF-0.002M TiCl₃.</p> <p>(15) Strip U from TRU with 15mL of 0.1M ammonium bioxalate.</p> <p>(16) Add 0.5mL 10% TiCl₃ to U samples, 0.5mL 30% H₂O₂ to Pu, and 0.2mL 30% H₂O₂ Am/Cm samples.</p>	<p>(23) Draw vacuum until filter is dry.</p> <p>(24) Remove filter from funnel assembly and mount filter on stainless steel planchet with 2-sided tape.</p>
<p>(6) Rinse TEVA cartridge with: -10mL 3M HNO₃ -20mL 9M HCl (remove Th) -5mL 3M HNO₃</p> <p>(7) Strip Pu(Np) from TEVA with 20mL 0.1M HCl-0.05M HF-0.01M TiCl₃.</p>	<p>(17) Add 50-100ug Ce carrier to all samples. Mix well. Add 1mL 49% HF. Mix well. Wait 15-20 minutes.</p> <p>(18) Set up Resolve® Filter Funnel on vacuum box.</p> <p>(19) Wet filter with 3mL 80% ethanol followed by 3mL DI water.</p>	<p>(25) Dry filter under heat lamp for 3-5 minutes.</p> <p>(26) Measure actinides by alpha spectrometry.</p>	
<p>(8) Rinse DGA with 8mL 0.1M HNO₃.</p> <p>(9) Place TRU cartridge above DGA.</p> <p>(10) Strip Am/Cm from TRU onto DGA with 15mL 3M HCl.</p> <p>(11) Separate TRU and DGA. Set TRU aside for U recovery.</p> <p>(12) Rinse DGA with: -5mL 3M HCl -3mL 1M HNO₃ -15mL 0.05M HNO₃</p>	<p>(20) Filter sample.</p> <p>(21) Rinse sample tube with 5mL DI water and add to filter.</p> <p>(22) Rinse filter funnel with 3mL DI water and 2mL 100% ethanol.</p>	 	

*Radiostrontium may also be measured by adding a 2mL + 1mL Sr Resin cartridge below DGA and following separation scheme in Eichrom application note AN-1405, "Rapid Determination of Sr in Vegetation Samples."

**Adding 50uL of 30% H₂O₂ to the tube rinse can help improve U recovery and decontamination in Pu/Np fractions.

Performance of Actinides in Vegetation Method

5 gram Samples						10 gram Samples					
Nuclide	Replicates	Reference (mBq/g)	Measured (mBq/g)	% Bias	% Tracer Recovery	Nuclide	Replicates	Reference (mBq/g)	Measured (mBq/g)	% Bias	% Tracer Recovery
²³⁸ Pu	6	29.4	30.1 ± 3.7	2.4	101 ± 6	²³⁸ Pu	2	27.4	28.1 ± 0.4	2.6	90 ± 15
²³⁹ Pu	6	56.8	57.0 ± 4.8	0.3	101 ± 6	²³⁹ Pu	2	32.8	32.4 ± 0.9	-1.2	90 ± 15
²⁴¹ Am	6	48.0	48.5 ± 4.6	1.0	93 ± 7	²⁴¹ Am	2	31.2	30.8 ± 0.0	-1.3	84 ± 12
²⁴⁴ Cm	6	6.28	5.9 ± 0.6	-6.1	93 ± 7	²³⁴ U	2	41.6	41.3 ± 1.3	-0.7	81 ± 12
²³⁴ U	6	69.2	81 ± 7	17	87 ± 7	²³⁸ U	2	43.2	42.0 ± 0.3	-2.8	81 ± 12
²³⁸ U	6	71.8	83 ± 10	16	87 ± 7						

References

1) Sherrod L. Maxwell, Brian K. Culligan, Gary W. Noyes, "Rapid separation of actinides and radiostrontium in vegetation samples," *J. Radioanal. Nucl. Chem.*, 286(1), 273-282 (2010).