

# Rapid Determination of Sr in Emergency Urine Samples

**Summary of Method** Strontium is separated and concentrated from 100mL urine samples using calcium phosphate precipitation. An optional wet-ashing step with  $\text{HNO}_3\text{-H}_2\text{O}_2$  destroys residual organic material. The precipitate or wet-ashed residue is dissolved in nitric acid and aluminum nitrate. Strontium is then separated from matrix impurities and potentially interfering radionuclides in the sample using a 2mL cartridge of Eichrom Sr Resin. Radiostrontium is measured on a low background gas flow proportional counter or liquid scintillation counter. Chemical yield of strontium is determined by gravimetric recovery of stable strontium or ICP-AES measurement. Typical chemical recovery of strontium is >80%. Measured values of  $^{90}\text{Sr}$  agreed to within 1.7% of reference values for 10 minute count times, although longer count times can be used to improve detection limits and uncertainty. A single operator can complete the separation method for batches of 12-24 samples in as little as 3-4 hours.


## Reagents

Sr Resin, 2mL Cartridges (Eichrom SR-R50-S)  
 Nitric Acid (70%)  
 Ammonium Hydroxide (listed as 28%  $\text{NH}_3$  or 56%  $\text{NH}_4\text{OH}$ )  
 Hydrogen Peroxide (30%)  
 Deionized Water  
 1.25M  $\text{Ca}(\text{NO}_3)_2$   
 3.2M  $(\text{NH}_4)_2\text{HPO}_4$   
 Sr Carrier (10mg/mL)  
 2M  $\text{Al}(\text{NO}_3)_3$   
 $^{90}\text{Sr}$  standard  
 Oxalic acid

## Equipment


Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)  
 Cartridge Reservoir, 20mL (Eichrom AR-200-RV20)  
 Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)  
 Yellow Outer Tips (Eichrom AR-1000-OT)  
 50mL and 250mL Centrifuge Tubes  
 Centrifuge  
 Cupped Stainless Steel Planchets (~5mL volume)  
 Gas Flow Proportional Counter  
 Hot Plate  
 Analytical Balance  
 250mL Glass Beakers  
 Vacuum Pump

**Figure 1. Sample Preparation**




## Figure 2. Load Solution Preparation and Strontium Separation

- (1) Precondition Sr Resin with 10mL 8M HNO<sub>3</sub>.
- (2) Load sample at 1-2mL/min.
- (3) Rinse sample tube with 3mL 3M HNO<sub>3</sub>.
- (4) Add tube rinse to Sr Resin. Elute at 1-2mL/min.
- (5) Rinse Sr Resin sequentially with:
  - 10 mL 8M HNO<sub>3</sub>
  - 5mL 3M HNO<sub>3</sub> - 0.05 oxalic acid
  - 5mL 8M HNO<sub>3</sub>
- (6) Dispose of (1) to (5) as waste.
- (7) Strip Sr with 15mL 0.05M HNO<sub>3</sub> at 1mL/min.



### Gas Flow Proportional Counting:\*

- (8) Evaporate samples to dryness on tared cupped stainless steel planchets.
- (9) Rinse Sr sample vials with 2mL 0.05M HNO<sub>3</sub>. Transfer vial rinse to planchets. Evaporate to dryness.



- (10) Weigh planchets on an analytical balance to determine gravimetric yield of stable Sr(NO<sub>3</sub>)<sub>2</sub>.

- (11) Measure radiostrontium in samples on low background gas flow proportional counter.

### \*(Options for <sup>89/90</sup>Sr Discrimination)

(a) Sr fraction from step (7) can be transferred to a liquid scintillation vial. <sup>89</sup>Sr can be measured by Cerenkov counting (no LSC cocktail). <sup>89/90</sup>Sr may then be measured after adding liquid scintillation cocktail.

(b) Sr fraction from step (10) can be dissolved in 10mL 8M HNO<sub>3</sub> after >7 days of <sup>90</sup>Y ingrowth. <sup>89/90</sup>Sr can be removed on Sr Resin. <sup>90</sup>Y will elute in Sr Resin load and can be counted by liquid scintillation or gas flow proportional counting.

Actinides may also be measured by adding a 2mL TEVA, TRU and DGA cartridges above Sr Resin and following the separation scheme in Eichrom application note AN-1412, "Rapid Determination of Actinides in Emergency Urine Samples."

## References

- 1) Sherrod L. Maxwell, Brian K. Culligan, "Rapid separation method for emergency water and urine samples," *J. Radioanal. Nucl. Chem.*, 279(3), 901-907 (2009).