

# Rapid Determination of $^{90}\text{Sr}$ in 10g Concrete Samples

**Summary of Method**  $^{90}\text{Sr}$  is determined by the direct separation of its daughter  $^{90}\text{Y}$  from 10 gram concrete samples. Samples are finely ground and fused in a zirconium crucible for 30 minutes at  $600^\circ\text{C}$  with 40 grams of sodium hydroxide. The fusion cake is dissolved in water, and strontium is concentrated and separated from the matrix using a ferric hydroxide precipitate. A secondary precipitation with Y/Ca-fluoride removes additional matrix (including silicates) and decreases the volume of precipitate. The Y/Ca-fluoride precipitate is dissolved with nitric acid-boric acid-aluminum nitrate to form the load solution.  $^{90}\text{Y}$  is separated from remaining matrix and potentially interfering radionuclides using stacked 2mL TRU and DGA Resin cartridges.  $^{90}\text{Y}$  is measured by gas flow proportional counting following microprecipitation onto Resolve® Filters. Chemical yields are determined by ICP-AES analysis. Batches of 12-24 samples can be prepared for analysis in less than 8 hours. This method is only suitable for aged samples, where the shorter lived  $^{89}\text{Sr}$  ( $t_{1/2} = 50.55$  days) and fission products such as  $^{91}\text{Y}$  are unlikely to be present. For samples not meeting this criterion,  $^{89/90}\text{Sr}$  can be determined from up to 5g concrete samples, using Eichrom Application Note AN-1605.

## Reagents

DGA Resin, 2mL Cartridges (Eichrom DN-R50-S)  
 TRU Resin, 2mL Cartridges (Eichrom TR-R50-S)  
 Yttrium Carrier (10mg/mL)  
 Iron Carrier (50mg/mL Fe, as ferric nitrate)  
 $^{90}\text{Sr}$  standard                      HF(49%)  
 Nitric Acid (70%)                      Sodium Hydroxide  
 Hydrochloric Acid (37%)              Deionized Water  
 1.25M  $\text{Ca}(\text{NO}_3)_2$                       3.2M  $(\text{NH}_4)_2\text{HPO}_4$   
 2M  $\text{Al}(\text{NO}_3)_3$                           Oxalic acid  
 Boric acid

## Equipment

Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)  
 Cartridge Reservoir, 20mL (Eichrom AR-200-RV20)  
 Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)  
 Yellow Outer Tips (Eichrom AR-1000-OT)  
 Resolve Filters with funnel (Eichrom RF-DF25-25PP01)  
 50mL and 500mL Centrifuge Tubes  
 Stainless Steel Planchets with two sided tape  
 250mL Zirconium crucibles with zirconium lids  
 Centrifuge                              Gas Flow Proportional Counter  
 Muffle Furnace                      Hot Plate/Heat Lamp  
 Analytical Balance                  Vacuum Pump

**Figure 1. Sample Preparation**

10g finely ground sample in zirconium crucible

↓  
 Add 2mg Y carrier.

Heat samples to dryness on hot plate.

↓  
 Add 40g of NaOH.

Cover crucibles with zirconium lid and place in muffle furnace at  $600^\circ\text{C}$  for 30 minutes.

↓  
 Carefully remove samples from furnace and cool in fume hood.

↓  
 Add 50-100mL of water and heat on hot plate to dissolve fusion cake.

↓  
 Transfer to a 500mL centrifuge tube. Rinse crucible with water. Dilute to 450mL with water.

Cool to room temp. Add 125mg Fe Mix. Centrifuge 10min. Decant supernate.



↓  
 Rinse precipitate with 150mL pH ~9 NaOH. Centrifuge. Decant Supernate. Repeat.

↓  
 Dissolve precipitate in 200mL 1.5M HCl. Add 50mL 0.01M HCl and 15mL 49%HF. Mix. Cool in ice bath 10min. Centrifuge 10min. Decant supernate.

↓  
 Dissolve precipitate in 7mL 3M  $\text{HNO}_3$ -0.25M Boric acid, 7mL conc.  $\text{HNO}_3$ , 7mL 8M  $\text{HNO}_3$  and 7mL 2M  $\text{Al}(\text{NO}_3)_3$ . Warm as needed.

**Figure 2. TRU-DGA Separation and Gas Flow Proportional Counting**

<p>(1) Precondition 2mL TRU* + DGA Resin with 10mL 8M HNO<sub>3</sub>.</p> <p>(2) Load sample at 1-2mL/min.</p> <p>(3) Rinse sample tube with 5mL 8M HNO<sub>3</sub>.</p> <p>(4) Add tube rinse TRU + DGA. Elute 1-2 mL/min.</p> <p>(5) Rinse sequentially with:</p> <ul style="list-style-type: none"> <li>- 5mL 6M HNO<sub>3</sub></li> <li>- 10mL 3M HNO<sub>3</sub></li> <li>- 12mL 4M HCl</li> </ul> <p>(6) Dispose of (1) to (5) as waste.</p> <p>(7) Discard TRU Resin.</p> <p>(8) Rinse DGA sequentially with:</p> <ul style="list-style-type: none"> <li>- 10mL 8M HNO<sub>3</sub> (Ca)</li> <li>- 15mL 0.1M HNO<sub>3</sub> (U)</li> <li>- 25mL 3M HNO<sub>3</sub>-0.25M HF (Th)</li> </ul> <p>(9) Strip Y with 18mL 0.25M HCl.</p> <p>(10) Dilute to 20mL with water.</p> <p>(11) Take 1mL aliquot and dilute to 20mL for ICP-AES yield measurement.</p>	<p>(8) Add 100ug Ce carrier to samples.</p> <p>(9) Mix well.</p> <p>(10) Add 2mL 49% HF.</p> <p>(11) Mix well.</p> <p>(12) Wait 15-20 minutes.</p> <p>(13) Set up Resolve® Filter Funnel on vacuum box.</p> <p>(14) Wet filter with 3mL 80% ethanol followed by 3mL DI water.</p> <p>(15) Filter sample.</p> <p>(16) Rinse sample tube with 5mL DI water and add to filter.</p> <p>(17) Rinse filter funnel with 3mL DI water and 2mL 100% ethanol.</p>	<p>(18) Draw vacuum until filter is dry.</p> <p>(19) Remove filter from funnel assembly and mount filter on stainless steel planchet with 2-sided tape.</p> <p>(20) Dry filter under heat lamp for 3-5 minutes.</p> <p>(21) Measure Y-90 by gas flow proportional counting.</p>
---	---	---



\*TRU Resin improves decontamination factors for U, Th and Bi isotopes, which could interfere with the measurement of <sup>90</sup>Y by gas flow proportional counting.

**Method Performance (10g gram Concrete, TRU/DGA Resin Method)**

Sample	% Y tracer recovery	<sup>90</sup> Sr Bq/g reference	<sup>90</sup> Sr Bq/g measured	% bias
1	81.7	0.0327	0.031	-5.4
2	83.3	0.0327	0.033	1.2
3	83.7	0.0327	0.031	-5.0
4	86.3	0.0327	0.033	-0.6
AVG	84 ± 2		0.032 ± 0.001	

**References**

1) Maxwell, Culligan, Hutchinson, Utsey, Sudowe, McAlister, "Rapid Method to Determine 89/90Sr in Large Concrete Samples," *J. Radioanal. Nucl. Chem.* accepted (2016). DOI 10.1007/s10967-016-4787-4