

# Rapid Determination of Pu in Steel Samples

**Summary of Method** Plutonium is separated and measured from 1-2 gram steel samples. Samples are digested with concentrated nitric, hydrochloric, and hydrofluoric acids. The digestate is evaporated to dryness, the residue dissolved in HNO<sub>3</sub>/H<sub>3</sub>BO<sub>3</sub>, and a CaF<sub>2</sub>/LaF<sub>3</sub> precipitate is used to concentrate the Pu and remove matrix. An optional NaOH fusion may also be performed, post sample digestion, to dissolve concrete or stone included in the sample and deal more rigorously with refractory Pu. Plutonium is separated from matrix impurities and potentially interfering radionuclides in the sample using 2 mL cartridges of Eichrom TEVA Resin. Plutonium is measured by alpha spectrometry following rare earth fluoride microprecipitation onto Eichrom Resolve filters. The chemical recovery of Pu, determined by <sup>242</sup>Pu tracer, was 90–99%. Measured values of Pu typically agreed to within 7-8% of reference values for 16 hour count times. The minimum detectable activity for Pu in 2 g samples with 16 hour count times was 0.25 mBq/g. A single operator can prepare batches of 12 samples for the measurement of Pu in less than 8 hours.

## Reagents


TEVA Resin, 2 mL Cartridges (Eichrom TE-R50-S)  
 Nitric Acid (70%)  
 Hydrochloric Acid (37%)  
 Hydrofluoric Acid (49%) or Ammonium Bifluoride  
 Lanthanum Carrier (10 mg/mL)  
 Cerium Carrier (10 mg/mL)  
 Deionized Water 1.25M Ca(NO<sub>3</sub>)<sub>2</sub>  
 2M Al(NO<sub>3</sub>)<sub>3</sub> <sup>242</sup>Pu Tracer  
 Boric acid NaNO<sub>2</sub>  
 Ascorbic Acid 30% H<sub>2</sub>O<sub>2</sub>  
 10-20% (w:w) TiCl<sub>3</sub> in HCl  
 3.2M (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub>\*  
 Sodium Hydroxide\*

## Equipment


Vacuum Pump  
 Centrifuge  
 Muffle Furnace\*  
 Hot Plate  
 Analytical Balance  
 Teflon Beakers (Zr Crucibles\*)  
 50 mL and 250 mL Centrifuge Tubes  
 Alpha Spectrometry System  
 Vacuum Box (Eichrom AR-24-BOX or AR-12-BOX)  
 Cartridge Reservoir, 20 mL (Eichrom AR-200-RV20)  
 Inner Support Tubes-PE (Eichrom AR-1000-TUBE-PE)  
 Yellow Outer Tips (Eichrom AR-1000-OT)  
 Resolve Filters in Funnel (Eichrom RF-DF25-25PP01)

**Figure 1. Sample Preparation**


Add 1-2 g steel sample to Teflon beaker\*.  
 \*If using optional fusion, omit HF/H<sub>3</sub>BO<sub>3</sub> and use Zr crucible.




**Figure 2. Load Solution Preparation and Plutonium Separation**

<p>(1) Ensure samples are cooled to room temperature.</p> <p>(2) Add 1.25 mL of 1.5M ascorbic acid. Mix well. Wait 5-10 minutes.</p> <p>(3) Add 1 mL of 3M NaNO<sub>2</sub>. Mix well. Wait 5-10 minutes.</p> <p>(4) Precondition 2 mL TEVA Resin with 5 mL 3M HNO<sub>3</sub>.</p> <p>(6) Load sample.</p> <p>(7) Rinse sample tube with 5 mL 8M HNO<sub>3</sub> + 50 µL 30% H<sub>2</sub>O<sub>2</sub>.</p> <p>(8) Add tube rinse to TEVA Resin.</p> <p>(9) Rinse TEVA Resin sequentially with:</p> <ul style="list-style-type: none"> <li>- 15 mL 3M HNO<sub>3</sub> (U decon.)</li> <li>- 20 mL 9M HCl (Th)</li> <li>- 5 mL 3M HNO<sub>3</sub></li> </ul>		<p>(10) Dispose of (4) to (9) as waste.</p> <p>(11) Strip Pu with 15 mL 0.1M HCl-0.05M HF-0.01M TiCl<sub>3</sub>.*</p> <p>*If preparing alpha sources by electrodeposition, replace TiCl<sub>3</sub> with Rongalite or hydroxylamine.</p> <p>(12) Add 50 µg Ce carrier and 0.5 mL 30% H<sub>2</sub>O<sub>2</sub> to all samples. Mix well.</p> <p>(13) Add 1 mL 49% HF. Mix well. Wait 15-20 minutes.</p> <p>(14) Set up Resolve® Filter Funnel on vacuum box.</p> <p>(15) Wet filter with 3 mL 80% ethanol followed by 3 mL DI water.</p> <p>(16) Filter sample.</p> <p>(17) Rinse sample tube with 5 mL DI water and add to filter.</p>	<p>(18) Rinse filter funnel with 3 mL DI water and 2 mL 100% ethanol.</p> <p>(19) Draw vacuum until filter is dry.</p> <p>(20) Remove filter from funnel assembly and mount filter on stainless steel planchet with 2-sided tape.</p> <p>(21) Dry filter under heat lamp for 3-5 minutes.</p> <p>(22) Measure actinides by alpha spectrometry.</p>
---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Filter assembly with 25mm, 0.1µm Resolve™ polypropylene filter





**Method Performance for Pu in Steel Samples**

Details	Sample replicates	Reference (mBq/sample)	Measured (mBq/sample)	Average % Diff.	242Pu tracer % Yield
238Pu in 2 g Steel	5	37.0	37.7 ± 1.6	4.2	89.3 ± 2.3
239Pu in 2 g Steel	5	24.5	24.4 ± 1.6	6.6	96.5 ± 3.4
239Pu (refractory) in 2 g Steel	5	24.5	23.4 ± 0.9	3.8	98.9 ± 6.6
239 Pu in 5 g Steel	4	37.0	38.3 ± 1.0	2.6	92 ± 14

**References**

1) Sherrod L. Maxwell, Brian K. Culligan, Jay B. Hutchison, Robin. C. Utsey, Ralf Sudowe, Daniel R. McAlister, "Rapid method to determine plutonium isotopes in steel samples," *J. Radioanal. Nucl. Chem.*, 314(2), 1103-1111 (2017).