

**eichrom**<sup>®</sup>  
TECHNOLOGIES

A GCI COMPANY

npo

eichrom



# Resolve Filters: New Polyethylene Material

Madeleine Eddy

1 November 2023

RRMC



PRODUCTS

**TECHNICAL INFO**

RESOURCES

+1 (630) 963-0320

**CONTACT**

- APPLICATIONS NOTES**
- AVAILABLE METHODS
- BIBLIOGRAPHY SEARCH
- USER WORKSOPS
- NEWSLETTER ARCHIVE

## Application Notes

Reference No.	Title	Download
	All Application Notes 2014-2018	<a href="#">Download</a> >
AN-1401	Rapid Determination of 226Ra in Emergency Urine and Water	<a href="#">Download</a> >
AN-1402	Rapid Determination of Sr in Emergency Milk Samples	<a href="#">Download</a> >
AN-1403	Rapid Determination of Sr in 50g Soil Samples	<a href="#">Download</a> >

https://www.eichrom.com/eichrom/applications-notes/



A GCI COMPANY




www.eichrom.com

# Common Recipes for CeF<sub>3</sub> Microprecipitates


Analyte	Matrix	15 mPre-treatment*	Ce(ug)	HF (mL)
Ac	10 mL 0.35 M HNO <sub>3</sub>	None	100	1.0
Ac	15 mL 2 M HCl	None	100	1.0
<b>Am/Cm</b>	<b>15 mL 4 M HCl</b>	<b>Dilute 2x</b>	<b>50</b>	<b>3.0</b>
Am/Cm	15 mL 0.25 M HCl	None	50	1.0
Am/Cm	20 mL 1 M HCl	None	50	1.0
Np/Pu	10 mL 0.1 M NH <sub>4</sub> bioxalate	None	50	1.0
Np/Pu	20 mL 0.1 M HCl – 0.01 M HF – 0.01 M TiCl <sub>3</sub>	None	50	1.0
Pu	25 mL 0.05 M HNO <sub>3</sub> – 0.05 M HF – 0.02 M TiCl <sub>3</sub>	None	50	1.0
<b>Th</b>	<b>15 mL 9 M HCl</b>	<b>Dilute to 40 mL</b>	<b>40</b>	<b>3.0</b>
U	15 mL 1 M HCl	None	100	1.0
U	10 mL 0.1 M NH <sub>4</sub> bioxalate	None	100	1.0
Y	15 mL 8 M HNO <sub>3</sub>	Dilute 2x	100	3.0

# CeF<sub>3</sub> Microprecipitate Method Overview


- 1) Dilute samples as necessary and add Ce Carrier (See Table I).
- 2) **U Samples**, Add 0.5mL 10% TiCl<sub>3</sub>
- 3) **Th, Np, Pu, Am/Cm samples requiring additional U decontamination:** Add 50 μL 30% H<sub>2</sub>O<sub>2</sub>.
- 4) Swirl to mix.
- 5) Add 1-3.5mL conc. HF (see Table I).
- 6) Swirl to mix.
- 7) Wait 20-30 minutes before filtration.




- 8) Ensure tight fit of filter assembly.
- 9) Engage vacuum.
- 10) Wet filter with 3-5mL 80% ethanol.
- 11) Wet filter with 3-5mL DI water.
- 12) Add sample.
- 13) Rinse tube with 5mL DI water and add to filter assembly.
- 14) Allow all fluid to pass through filter.
- 15) Rinse filter funnel with 3-5mL DI Water.
- 16) Rinse filter funnel with 2-3mL ethanol.
- 17) Filter until dry.



- 18) Remove filter (lifting on edges).



- 19) Mount filter to stainless steel planchet.
- 20) Place mounted filter into petri dish.
- 21) Dry under heat lamp.
- 22) Place lid on petri dish. Store until alpha measurement.



# 2022 Conclusions

- **New Filter Material**
  - Polyethylene is a viable replacement for Resolve® Filters
  - PE performance is equivalent to PP for standard actinide precipitate conditions
  - PE does experience slight curling
- **Review of RE Precipitate Method**
  - Precipitate quality is unaffected by different amounts of rare earth (25-200 ug) or HF (1-3 mL) added
  - Peroxide is detrimental to high HCl samples
  - Finishing with 100% ethanol reduces filter curling
  - Finishing/mounting/drying has little impact on spectra quality but does impact handling

# Overview of Customer Comments/Questions

- **Physical**

- Filters not centered on funnel
- Oval-shaped filters
- Wrinkling/ballooning of filters in funnel

- **Chemical**

- High FWHM with tailing for Am/Cm samples
  - Specifically for TRU resin eluents
- High FWHM for Th samples

- **Kinetic**

- How long do precipitation reactions need?


# 2023 Experiments

- **Physical Handling**
  - Wrinkled/off-center filter performance
  - Filter durability during long-term storage
- **Aqueous Phase Chemistry**
  - Clean acid vs. column eluent
  - Varied [acid], and acid neutralization
- **Precipitation Kinetics**
  - Filtration of Am/Cm samples at short time-points


# Examples of Distorted Filters



Oval-shaped and wrinkled loose filters



Filters placed off-center on funnels




Ballooned filter removed from funnel




# Resolve Funnel vs. Pall Gelman Funnel



Comparison of filter edge width



Pall Gelman filter with (L) and without (R) mesh insert



Oval-cut filter on Pall Gelman


# Physical Handling “Worst Case Scenario” Tests

## Experimental Conditions:

- 200 dpm  $^{243}\text{Am}$ , 100 dpm  $^{241}\text{Am}$ , and  $^{244}\text{Cm}$
- 10 mL 1 M HCl
- 50 ug Ce
- 1 mL conc. HF
- 20 min precipitation time

Relative Percent (%) recoveries of physically handled sample sets vs. control group (center-placed PP fliters)

Sample Set	$^{243}\text{Am}$	$^{241}\text{Am}$	$^{244}\text{Cm}$
PE – off-center	120 ± 20	105 ± 4	102 ± 4
PP – dropped	105 ± 3	105 ± 2	104 ± 4
PE – dropped	104 ± 3	103 ± 2	102 ± 2
PP – Pall Gelman	100 ± 3	98 ± 2	100 ± 3
PE – Pall Gelman	102 ± 5	93 ± 4	95 ± 4
PE – Pall Gelman w/mesh	102 ± 2	98 ± 3	100 ± 2





A GCI COMPANY

[www.eichrom.com](http://www.eichrom.com)

# Durability/Storage Testing

## Experimental Conditions:

- 200 dpm  $^{243}\text{Am}$ , 100 dpm  $^{241}\text{Am}$ , and  $^{244}\text{Cm}$
- 10 mL 1 M HCl
- 50 ug Ce
- 1 mL conc. HF
- 20 min precipitation time



# Durability/Storage Testing

Percent Change in Activity for PP Samples

Sample	<sup>243</sup> Am	<sup>241</sup> Am	<sup>244</sup> Cm
Initial Activity (dpm)	348	144	138
Two Weeks	+0.9%	-0.2%	-0.4%
Two Weeks – shaken	-3.9%	-2.7%	-4.0%
Three Months	-5.3%	-5.8%	-6.9%
Six Months	-5.2%	-5.9%	-6.1%
Eight Months	-1.0%	-1.9%	-2.2%

Percent Change in Activity for PE Samples

Sample	<sup>243</sup> Am	<sup>241</sup> Am	<sup>244</sup> Cm
Initial Activity (dpm)	242	141	137
Two Weeks	+0.7%	-1.2%	-0.7%
Two Weeks – shaken	+0.0%	-0.7%	-1.4%
Three Months	+0.5%	-2.7%	-2.9%
Six Months	-1.1%	-2.7%	-2.3%
Eight Months	-1.0%	-2.7%	-2.3%

# Durability/Storage Testing

## Percent Change in FWHM for PP Samples

Sample	<sup>243</sup> Am	<sup>241</sup> Am	<sup>244</sup> Cm
Initial FWHM (keV)	37.9	39.7	29.9
Two Weeks	+5.4%	+3.4%	+16%
Two Weeks – shaken	+4.8%	+3.0%	+13%
Three Months	+38%	+34%	+82%
Six Months	+16%	+11%	+55%
Eight Months	+16%	+12%	+50%

## Percent Change in FWHM for PE Samples

Sample	<sup>243</sup> Am	<sup>241</sup> Am	<sup>244</sup> Cm
Initial FWHM (keV)	39.4	41.1	34.9
Two Weeks	+6.9%	+8.2%	+17%
Two Weeks – shaken	+6.0%	+6.1%	+9.5%
Three Months	+32%	+33%	+67%
Six Months	+14%	+16%	+48%
Eight Months	+15%	+13%	+50%


# Physical Handling – Conclusions and Process Changes

## Conclusions

- Minimal change in % recovery unless gaps in filter
  - Can identify gaps easily during wetting by “whooshing” of EtOH
- Pall Gelman funnels larger active area – more important to fully center filters
- Minimal activity loss but significant increase in FWHM with long-term storage

## Process Changes

- Warning labels
- Less compression to funnel stacks
- Adjust dye cutting pressure/method to reduce jagged edges




# Clean Acid vs. Column Eluent


TEVA	UTEVA	TRU	DGA
PC/L/R 30 mL 3 M HNO <sub>3</sub>	PC/L/R 30 mL 3 M HNO <sub>3</sub>	PC/L/R 30 mL 3 M HNO <sub>3</sub>	PC/L/R 30 mL 3 M HNO <sub>3</sub>
<b>Th – 15 mL 9 M HCl</b>	Rinse – 20 mL 5 M HCl/0.05 M oxalic acid	<b>Am – 15 mL 4 M HCl</b>	Rinse – 20 mL 0.25 M HNO <sub>3</sub>
<b>Pu/Np – 15 mL 0.1 M HCl/0.05 M HF/0.01 M TiCl<sub>3</sub></b>	<b>U – 15 mL 1 M HCl</b>	Rinse – 4 M HCl/0.25 M HF	<b>Am – 0.1 M HCl</b>
		<b>U or Pu/Np – 15 mL 0.1 M ammonium bioxalate</b>	

# Clean Acid vs. Column Eluent


### 50 dpm Th: TEVA




### 20 dpm U: UTEVA and TRU



### 100 dpm Am/Cm: TRU and DGA



### 100 dpm Np/Pu: TEVA and TRU





# Clean Acid vs. Column Eluent



TEVA		UTEVA		TRU		DGA	
<b>Th:</b> 15 mL 9 M HCl diluted to 40 mL		<b>U:</b> 15 mL 1 M HCl		<b>Am:</b> 15 mL 4 M HCl diluted to 30 mL		<b>Am:</b> 15 mL 0.1 M HCl	
<i>dpm</i> -6.5%	<i>FWHM</i> +4.1%	<i>dpm</i> -0.4%	<i>FWHM</i> -8.5%	<i>dpm</i> -25%	<i>FWHM</i> -2.2%	<i>dpm</i> -33%	<i>FWHM</i> -5.1%
<b>Pu/Np:</b> 15 mL 0.1 M HCl/0.05 M HF/ 0.01 M TiCl <sub>3</sub>				<b>U:</b> 15 mL 0.1 M ammonium bioxalate			
<i>dpm</i> -8.4%	<i>FWHM</i> -2.8%			<i>dpm</i> -4.7%	<i>FWHM</i> -11%		
				<b>Pu/Np:</b> 15 mL 0.1 M ammonium bioxalate			
				<i>dpm</i> +3.1%	<i>FWHM</i> +12%		

# Am/Cm Acid Dependence

## Experimental Conditions:


- 20 dpm  $^{243}\text{Am}$ , 10 dpm  $^{241}\text{Am}$ , and  $^{244}\text{Cm}$
- 50 ug Ce
- 1 mL conc. HF
- 20 min precipitation time

Volume (mL)	[HCl] (M)
15	0.1
50	1.2
30	2.0
15	4.0
15 + 7.6 mL conc $\text{NH}_4\text{OH}$	4.0 + 3.9 M $\text{NH}_4\text{OH}$




# Am/Cm Acid Dependence


## Am-243 FWHM for PP v PE with varied [HCl]



## Am-241 FWHM for PP v PE with varied [HCl]



## Cm-244 FWHM for PP v PE with varied [HCl]





# Th Acid Dependence

## Experimental Conditions:

- 200 dpm  $^{230}\text{Th}$ ,  $^{229}\text{Th}$ , and  $^{228}\text{Th}$
- 40 ug Ce
- 3 mL conc. HF
- 20 min precipitation time

Volume (mL)	[HCl] (M)
30	2
45	3
15	4
15 + 16.9 mL conc $\text{NH}_4\text{OH}$	9.0 + 8.9 M $\text{NH}_4\text{OH}$

## Th FWHM for PP v PE with varied [HCl]






# Aqueous Phase Chemistry Conclusions

- **Acid vs eluent**
  - No significant difference in yield for acid vs eluent
  - Some increase in FWHM for Np/Pu-TRU sample
- **Am/Cm vs [HCl]**
  - FWHM PE > FWHM PP
  - Significant increase in FWHM with increasing [HCl]
  - Neutralizing HCl restores original spectral quality
- **Th vs [HCl]**
  - FWHM PE > FWHM PP
  - No significant dependence on [HCl] but PE > PE for all conditions
  - Neutralizing HCl greatly improves spectral quality

# Short-term Delayed Filtration Test for Am/Cm

## Experimental Conditions:

- 200 dpm  $^{243}\text{Am}$ , 100 dpm  $^{241}\text{Am}$ , and  $^{244}\text{Cm}$
- 10 mL 1 M HCl
- 50  $\mu\text{g}$  Ce
- 1 mL conc. HF
- **VARIED** precipitation time



# Current Work: Investigation of Thicker PE Membrane

- **Motivation**
  - Determine if thicker filters are viable alternative to new PE materials
  - Hope that the thicker membrane may have fewer physical issues related to curling and shifting in funnels
- **Experiments**
  - Perform standard QC to determine product quality
  - Test filter curling
  - Test Am/Cm sample acid dependence
- **Long-term Objectives**
  - If we proceed with this new material, we will monitor customer comments related to filters shifting and ballooning in funnels to determine if we see a decrease in frequency. We hope the thicker material will be sturdier and create a tighter fit in the funnels which will prevent it from slipping during shipping and handling.

# Questions?

For more information on alpha spectrometry please join Eichrom at ORTEC's Alpha Spectrometry Training Course from October 14-18, 2024, at GEL Laboratories in Charleston, SC

<https://www.ortec-online.com/service-and-support/training/alpha-spectrometry>



Laboratories





# Acid vs Elution Multi-day Test

- Conduct mock columns to collect “eluent” samples
- Prepare 10x samples for each acid and eluent solution with tracers
- Add Ce and HF to 5x samples
- Delay precipitation for other 5x samples
- Each day prepare one delayed filtration and delayed precipitation sample for each solution
- Monitor how DPM and FWHM are affected over time


TEVA	UTEVA	TRU	DGA
PC/L/R 30 mL 3 M HNO <sub>3</sub>	PC/L/R 30 mL 3 M HNO <sub>3</sub>	PC/L/R 30 mL 3 M HNO <sub>3</sub>	PC/L/R 30 mL 3 M HNO <sub>3</sub>
<b>Th – 15 mL 9 M HCl</b>	Rinse – 20 mL 5 M HCl/0.05 M oxalic acid	<b>Am – 15 mL 4 M HCl</b>	Rinse – 20 mL 0.25 M HNO <sub>3</sub>
<b>Pu/Np – 15 mL 0.1 M HCl/0.05 M HF/0.01 M TiCl<sub>3</sub></b>	<b>U – 15 mL 1 M HCl</b>	Rinse – 4 M HCl/0.25 M HF	<b>Am – 0.1 M HCl</b>
		<b>U or Pu/Np – 15 mL 0.1 M ammonium bioxalate</b>	

## TEVA






# Uranium

## TRU





## UTEVA





# Americium

## TRU





## DGA





# Neptunium/Plutonium

## TRU




## TEVA



# Kinetics Conclusions

- **Short-term**
  - Activity unaffected, but slight increase in FWHM over time
- **Multi-day**
  - No relationship between FWHM and day of filtration/precipitation
  - Lower FWHM for Th-TEVA acid samples
  - Decrease in activity/yield for U-UTEVA eluent sample over time
  - Am-DGA samples lower FWHM than Am-TRU
    - Related to [HCl]

# Current Work: Investigation of Thicker PE Membrane



# Investigations of Thicker PE Membrane

## FWHM for thicker PE with varied [HCl]

