


**eichrom®**

A BRAND OF  
EICHROM TECHNOLOGIES



## Eichrom 101: Separations Basics

# https://www.eichrom.com/eichrom/presentations



The screenshot shows a web browser displaying the Eichrom website. The URL in the address bar is <https://www.eichrom.com/eichrom/presentations/>. The website has a blue header with the Eichrom logo and navigation links: PRODUCTS, TECHNICAL INFO (highlighted with a red box), and RESOURCES. On the right side of the header, there is a phone number +1 (630) 963-0320 and a CONTACT button. Below the header, there are several menu items: APPLICATIONS NOTES, AVAILABLE METHODS, BIBLIOGRAPHY SEARCH, and USER WORKSOPS (highlighted with a red box). The main content area features a table with two columns: 'Year' and 'Presentation'. The table lists presentations from 2014 to 2023.

Year	Presentation
2023	<a href="#">12th International Symposium on Targeted Alpha Therapy</a>
2022	<a href="#">2022 RRM</a>
2021	<a href="#">2021 Virtual UGM</a>
2019	<a href="#">64th RRM, Santa Fe, NM, October 27 - November 1, 2019</a>
2018	<a href="#">63rd RRM, Portland, ME May 20-25, 2018</a>
2017	<a href="#">62nd RRM, Honolulu, HI February 6-10, 2017</a>
2015	<a href="#">61st RRM, Iowa City, IA October 26-30, 2015</a>
2014	<a href="#">60th RRM, Knoxville, TN October 27-31, 2014</a>

# Separation Types


Gas – Liquid	Gas – Solid	Liquid – Liquid	Liquid – Solid
Disk	Adsorption	Solvent Extraction	Precipitation
Gas Chromatograph	Sublimation	Exclusion	Fractional Crystallization
	Molecular Sieves		Ion Exchange
	Gas Chromatograph		Extraction Chromatography
			Adsorption
			Ion Exclusion

Phase 1: Aqueous

Phase 2:  
Solid or Liquid


Elements or groups of elements or compounds

# Separation Types




**Aqueous Phase**

**Phase  
Boundary**




$PuF_4 (s)$

$Pu(NO_3)_4 \cdot (DAAP)_2 \cdot S_x$



**Non-Aqueous Phase**

# Aqueous Phase




M = metal ion

n = solvation number


z = charge

Solubility related to strength  
of solvation vs interaction  
between cations and anions



# Hydrolysis

Hydrolysis - Splitting of  $\text{H}_2\text{O}$  into  $\text{OH}^-$  and  $\text{H}^+$



# ADSORPTION

7

Point of zero charge – pH at which surface charge is 0.

At higher pH


- Surface charge is negative/Cations adsorb

At lower pH

- Surface charge is positive/Anions adsorb

Adsorption normally increases with charge

$M^+ < M^{2+} < M^{3+} < M^{4+}$



Borosilicate Glass

1L Bottle  $10^{-7}$ - $10^{-8}$  mole H<sup>+</sup>

# Storage Conditions

- Aqueous samples often stored in  $\text{HNO}_3$  at  $\text{pH} < 2$ .
- Radioactive tracers may be stored in 0.1 – 3.0 M  $\text{HNO}_3$  or  $\text{HCl}$  to limit hydrolysis and adsorption.
- Glass vs plastic.
- **Carrier** may be added to limit adsorption.



# Carrier

Stable element added to a system to provide additional mass.

- Same element or different element with similar chemistry
- prevent adsorption to container walls.
- provide mass to exceed solubility to promote precipitation (Stable Sr/Ca for Sr-89/90)
- gravimetric yield monitor (Stable Sr)
- prevent adsorption to precipitate (hold-back carrier)
  - Adding Co, Cs when precipitating Ni

# Common Yield Tracers

Analyte	Tracer	Tracer Measurement
Fe, Ni, Sr, Pb	stable element	Gravimetric, ICP-MS, ICP-AES
$^{99}\text{Tc}$	Stable Re, $^{99\text{m}}\text{Tc}$	Re: ICP-MS, ICP-AES, $^{99\text{m}}\text{Tc}$ : Gamma
$^{210}\text{Po}$	$^{209}\text{Po}$ , $^{208}\text{Po}$	Alpha
$^{226}\text{Ra}$	$^{133}\text{Ba}$ , $^{225}\text{Ra}$ ( $^{229}\text{Th}$ )	$^{133}\text{Ba}$ : Gamma, $^{225}\text{Ra}$ : Alpha
Th	$^{229}\text{Th}$	Alpha
U	$^{232}\text{U}$ (self-cleaning)	Alpha
Np	$^{236}\text{Pu}$ , $^{239}\text{Np}$ ( $^{243}\text{Am}$ )	$^{236}\text{Pu}$ : Alpha, $^{239}\text{Np}$ : Gamma
Pu	$^{242}\text{Pu}$ , $^{236}\text{Pu}$	Alpha
Am/Cm	$^{243}\text{Am}$	Alpha

# Oxidation State Adjustments

Element	Common Oxidation States in Acidic Aqueous Solution	Adjustments
Ac	Ac(III)	
Th	Th(IV)	High acidity to reduce hydrolysis
Pa	Pa(V)	Complex with F <sup>-</sup> to reduce hydrolysis
U	U(IV), <b>U(VI)O<sub>2</sub><sup>2+</sup></b>	Reduce to U(IV) with TiCl <sub>3</sub>
Np	Np(IV), <b>Np(V)O<sub>2</sub><sup>+</sup></b> , Np(VI)O <sub>2</sub> <sup>2+</sup>	Reduce to Np(IV) with Fe(II)
Pu	Pu(III), <b>Pu(IV)</b> , Pu(VI)O <sub>2</sub> <sup>2+</sup>	Reduce to Pu(III) with Fe(II) Oxidize to Pu(IV) with NO <sub>2</sub> <sup>-</sup>
Am/Cm	Am(III)/Cm(III)	

# Precipitation

## Beginning of methods

- Concentrate large aqueous samples
- Reject unwanted matrix ions
- Easy to dissolve in acidic matrix for further separations

## End of methods

- Provide low mass, concentrated sample for efficient measurement
- Reject unwanted interferences
- Difficult to dissolve to make stable sources

# Hydroxide precipitates

Actinides form acidic cations in aqueous solution and readily hydrolyze to precipitate as hydroxides.

La(OH)<sub>3</sub> and Fe(OH)<sub>3</sub> are often used as carriers for actinides in all oxidation states (rejects Sr):

III (k<sub>sp</sub> ~10E-20)

IV (k<sub>sp</sub> ~10E-54)

V (k<sub>sp</sub> ~10E-10)

VI (k<sub>sp</sub> ~10E-25)

Hydroxides are normally easy to dissolve in acid for further processing.

# Phosphate precipitates

Ca phosphate and Fe(III) phosphate are often used as carriers for actinides in all oxidation states + Sr.

Phosphates are normally easy to dissolve in acid for further processing.

Phosphate may complex actinides (Th/Pu/Np), interfering with further separation steps. Addition of  $\text{Al}(\text{NO}_3)_3$  during precipitate dissolution can help reduce the impact of phosphate.

## Fluoride and Oxalate precipitates

Rare earth fluorides and oxalates are often used as carriers for actinides in the III and IV oxidation states.

Oxalates are normally easy to dissolve in acid for further processing.

Fluorides can be dissolved by the addition of Boric acid (complexes  $F^-$  to form fluoboric acid).

Fluoride precipitates are particularly useful for removal of silicates from soil/rock/concrete samples.

## Rare Earth Fluoride micro. ppt


CeF<sub>3</sub> or NdF<sub>3</sub> are alternatives to electrodeposition for alpha source preparation.

CeF<sub>3</sub> or NdF<sub>3</sub> microprecipitates (50-150ug) can be prepared directly from HNO<sub>3</sub> or HCl solution.

CeF<sub>3</sub> or NdF<sub>3</sub> microprecipitates (50-150ug) can be used to provide additional decontamination from U.


- Add H<sub>2</sub>O<sub>2</sub>, U(VI) not carried on CeF<sub>3</sub>
- Add TiCl<sub>3</sub>, U(IV) carried on CeF<sub>3</sub>





### CeF<sub>3</sub> filter

- Faster
- Simpler
- Adequate resolution
- Less durable (contamination)
- Additional U purification




### Electrodeposition

- Slower
- More complex
- Superior resolution
- More durable

# Extraction Chromatography


# Extraction Chromatographic Resin




- Solvent extraction reagents on inert chromatographic support.
- High selectivity.
- Less waste.
- Less Labor.
- Rapid.
- Can be applied in glovebox or hot cell.

# Multi-dimensional separations

- Different column materials can be stacked
- Analysis of multiple analytes from a single sample preparation




# Acid Dependency Curve




- Data often presented as  $k'$  vs acid concentration
- $k'$  = capacity factor
  - Free column volumes to peak maximum
  - Relative retention strength
  - Simple way to view selectivity
  - $\Delta k'$  is one factor in resolution

# Column Chromatography Principles



# Elution Chromatogram




- $\Delta k'$  is one factor in resolution
  - Resin Type
  - Mobile Phase
- Peak width also affects resolution
  - Particle Size
  - Temperature

# Particle Size

Comparison of Elution Curves for  $\text{Sr}^{2+}$  for Two Particle sizes of Sr Resin

Eluent: 3.2 M  $\text{HNO}_3$ , 23-24°C




- Smaller particle size
  - Narrower elution bands
  - Higher back pressure
- Particle size affects
  - Breakthrough
  - Tailing




# Resin Particle Sizes (Grades)


A-grade 100-150 $\mu$



S-grade 50-100 $\mu$



F-grade 25-53 $\mu$



- Part number: XX-XXX-A
  - Gravity Flow ‘columns’
- Part number: XX-XXX-S
  - Vacuum assisted ‘cartridges’
- Part number: XX-XXX-F
  - Specialty applications
  - Small cartridges
  - Rare Earth separations

Questions????



[www.eichrom.com](http://www.eichrom.com)