

Rapid Determination of Ra-226 in Environmental Samples

S. L. Maxwell, B.K. Culligan, and P. J. Shaw

Savannah River National Laboratory
Aiken, SC
November 3, 2011

57th Radiobioassay and Radiochemical Measurements Conference

Destin, FL

Background

- Need for rapid radiochemical methods
 - Emergency response
 - IND, RDD, nuclear accident
 - Rapid turnaround times
- Ra-226
 - EPA and DOE FRMAC radionuclide concern list
 - alpha emitter, 4.78 MeV (94.5%), 4.61 MeV (5.55%)
 - 1600 year half-life (alpha spectrometry and ICP-MS)
 - radiotoxic, follows calcium in food chain into bones
- Newer, more rapid analytical methods are needed
 - Solid samples
 - Radon emanation/ gamma methods take time (ingrowth)
 - Yield/uncertainty issues?
 - Urine
 - Can be time-consuming



Literature: Water Methods


- Moon DS, Burnett WC, Nour S, Horwitz P, Bond A. 2003. Preconcentration of radium isotopes from natural waters using MnO₂ resin. *Applied Radiation and Isotopes* 59:255–62.
- W.C. Burnett, et al, Radium-228 determination of natural waters via concentration on manganese dioxide and separation using Diphonix ion exchange resin, *Appl. Radiation and Isotopes*, 2004, Dec 61, 1173
- S.L. Maxwell, Rapid method for ²²⁶Ra and ²²⁸Ra analysis in water samples *Journal of Radioanalytical and Nuclear Chemistry*, Vol. 270, No.3 (2006) 651–655

Ra-226 using MnO₂ resin



Ba-133 tracer

Ln Resin and DGA Resin for Ra-226 and Ac-228 separation



Journal of Radioanalytical and Nuclear Chemistry, Vol. 270, No.3 (2006) 651–655
Rapid method for ²²⁶Ra and ²²⁸Ra analysis in water samples
S. L. Maxwell

Literature: Water Methods

- Galit Sharabia, Boaz Lazarb, Yehoshua Kolodny, Nataliya Teplyakov, Ludwik Halicz, High precision determination of ^{228}Ra and $^{228}\text{Ra}/^{226}\text{Ra}$ isotope ratio in natural waters by MC-ICPMS
 - Radium is pre-concentrated from 60 L of Lake Kinneret water (or 2 L of the saline springs water) by co-precipitation with MnO_2 .
 - Separation of radium from the matrix elements was by large BioRad® AG50W-X8 and Sr Resin
- Rapid Radiochemical Methods for Selected Radionuclides in Water for Environmental Restoration Following Homeland Security Events, EPA 402-R-10-001 February 2010 (RRMC 2010 presentation. R. Litman, S. Workman
 - MnO_2 resin and Ra-225 tracer

Literature: Solid Samples

- Chabaux, F., Othman, D.B., Birck, J.L., 1994. A new Ra–Ba chromatographic separation and its application to Ra mass-spectrometric measurement in volcanic rocks., *Chemical Geology* 114, 191–197.
 - Use of multiple, large cation resin columns
 - Sr Resin to separate Ba from Ra
- Larivière, D. K. Brownell, V. N. Epov, R. J. Cornett, and R. D. Evans, Determination of ^{226}Ra in sediments by ICP-MS: A comparative study of three sample preparation approaches
 - Tested various combinations of large cation resin, sulfate precipitation and manganese dioxide
 - Multiple columns
 - 25 ml cation resin (**345 ml 6M HCl total rinse volume on 1st column** in Chabaux option)
- M.T. Crespo. On the determination of ^{226}Ra in environmental and geological samples by α -spectrometry using ^{225}Ra as yield tracer, *Applied Radiation and Isotopes* 53 (2000) 109-114
 - Th-229(Ra-225) tracer/At-217

Ra-226 in urine


- K. Kehagia, C. Potiriadis, S. Bratakos, V. Koukouliou and G. Drikos, Determination of ^{226}Ra in Urine Samples by Alpha Spectrometry, Radiation Protection Dosimetry (2007), pp. 1–4
 - Evaporation, wet-ash, overnight in furnace, redissolution, lead sulfate ppt.
 - Anion resin, cation resin, EDTA, electrodeposition
 - 2 days to prep 4 samples, wait 4 days to begin count
 - 7200 minute count time, 50-70% yields, Ra-225 (At-217)

Ra-226 in Emergency Response Samples

- Ra-226 on emergency response radionuclide list for U.S. EPA and DOE FRMAC
- Challenges
 - Difficulty separating Ca from Ra/Ba
 - Adequate tracer? (Ba-133)
 - Barium interference on alpha spectrometry source preparation
 - Poor alpha resolution
 - Isobaric interferences using ICP-MS (Sr+Ba, Sr + some lanthanides)
 - MnO₂ resin can be used for waters but Ra precipitates with Fe(OH)₃ at high pH for solid samples containing Fe
 - MnO₂ separates Ca in water but Ca may ppt when used as a MnO₂ ppt from solids
- SRS has developed new methods for alpha or ICP-MS
 - Concrete, brick, vegetation, air filters, urine
 - Ra-225 tracer (At-217 alpha)

Loss of Ra when Fe Is present (High pH)


Fe(OH)₃ Precipitation Recoveries vs pH (NaOH)



Courtesy Dr. Dan McCalister, P&G Research Foundation


Ra-225 Tracer Decay

^{225}Ra (Including Parent) Decay Scheme



Rapid Radiochemical Methods for Selected Radionuclides in Water for Environmental Restoration Following Homeland Security Events, EPA 402-R-10-001 February 2010

Alpha Spectra



Rapid Radiochemical Methods for Selected Radionuclides in Water for Environmental Restoration Following Homeland Security Events, EPA 402-R-10-001 February 2010

SRS Approach

- Use Ra-225 (At-217)
- Use sodium hydroxide fusion for solid samples
 - refractory particles, fast
- Use calcium carbonate precipitation
 - overcome $\text{Fe}(\text{OH})_3$ loss of Ra at high pH and column flow problems using MnO_2 resin at high pH from solid matrices
 - no sulfate to carbonate conversion
 - simplify column procedures
- Remove Ba using Sr Resin
 - to eliminate alpha peak resolution interference
 - to allow large sample aliquots and reduced count times

Rapid Ra-226 Solid Sample Preparation



Air Filter
Add Th-229 (Ra-225) tracer

5g Vegetation
Add Th-229 (Ra-225) tracer

1g Concrete/Brick
Add Th-229 (Ra-225) tracer

Heat at 700°C for 2 hours in Furnace

Wet Ash with HNO₃/H₂O₂ on Hot Plate


Heat for ~10 min. at 600°C in Furnace

Fuse in Zr Crucible for 15 min. with 10g NaOH at 600°C

Carbonate Precipitation (10mL conc. HCl, Ca*, 25mL 2M Na₂CO₃)

Redissolve in 20mL 1.5M HCl Add 3mL 1.5M Ascorbic Acid

Column Load Solution




*Ca addition:
Air Filter: 125mg
Vegetation: 50mg
Concrete/Brick: 25mg




Rapid Ra-226 Aqueous Sample Preparation

Wet-ash urine ppt



*Ca addition:
Water: 100-150mg
Urine: 100mg





Ba Removal

Figures 2 and 3

Acid dependency of k' for various ions at 23-25°C.

Sr Resin




Ba/Sr retained

Ra eluted

3M HNO_3 <10 free column volumes


Horwitz, et al., (HP292)

Ln Resin



Courtesy Dr. Dan McCalister, P&G Research Foundation


Rapid Ra-226 Column Separation




Option: ICP-MS

Ra-226 Column Separation

Ca Removal




Sr Resin




Only if Ba removal needed

Ln Resin



Ready for barium sulfate microprecipitation or ICP-MS

Ra-226 Spectra using At -217 tracer



Ra-226 in Vegetation Samples

Sample ID	²¹⁷ At Yield (%)	²²⁶ Ra Reference Value (mBq smp ⁻¹)	²²⁶ Ra Measured Value* (mBq smp ⁻¹)	Difference (%)
1	91.5	73.8	70.8	-4.0
2	88.3	73.8	73.8	0.0
3	93.1	73.8	69.8	-5.4
4	82.2	73.8	68.5	-7.2
5	80.2	73.8	81.4	10.3
Avg	87.1		72.8	-1.2
SD	5.7		5.1	
% RSD	6.5		7.1	

8hr count time

At -217 ingrowth to mid-point-9.03 hrs.

*corrected for 9.17 mBq Ra-226 native content

Ra-226 in Concrete Samples

Sample ID	²¹⁷ At Yield (%)	²²⁶ Ra Reference Value (mBq smp ⁻¹)	²²⁶ Ra Measured Value * (mBq smp ⁻¹)	Difference (%)
1	88.2	184.5	173.5	-6.0
2	90.8	184.5	188.3	2.1
3	81.2	184.5	172.8	-6.3
4	84.2	184.5	192.4	4.3
5	72.9	184.5	177.6	-3.7
6	90.4	184.5	178.8	-3.1
Avg	84.6		180.6	-2.1
SD	6.8		8.0	
% RSD	8.1		4.4	

4hr count time

At -217 ingrowth to mid-point-11.86 hrs.

*corrected for 26.8 mBq Ra-226 native content

Ra-226 in Brick Samples

Sample ID	²¹⁷ At Yield (%)	²²⁶ Ra Reference Value (mBq smp ⁻¹)	²²⁶ Ra Measured Value* (mBq smp ⁻¹)	Difference (%)
1	77.9	73.8	79.5	7.8
2	88.4	73.8	84.5	14.5
3	86.3	73.8	80.3	8.9
4	90.6	73.8	88.5	19.9
5	79.7	73.8	74.7	1.3
6	91.5	73.8	75.9	3.0
7	94.9	73.8	71.7	-2.8
Avg	87.0		79.3	7.5
SD	6.2		5.8	
% RSD	6.5		6.5	

4hr count time

At -217 ingrowth to mid-point-12.15 hrs.

*corrected for 29.5 mBq Ra-226 native content

Ra-226 in Air Filter Samples

Sample ID	²¹⁷ At Yield (%)	²²⁶ Ra Reference Value (mBq smp ⁻¹)	²²⁶ Ra Measured Value * (mBq smp ⁻¹)	Difference (%)
1	80.7	73.8	70.5	-4.4
2	79.9	73.8	80.8	9.6
3	78.6	73.8	77.0	4.4
4	73.0	73.8	79.5	7.8
5	71.5	73.8	77.7	5.3
Avg	76.7			
SD	4.2		77.1	4.5
% RSD	1.4			

8hr count time

At -217 ingrowth to mid-point-12.47 hrs.

*corrected for 3.7 mBq Ra-226 in blank air filter

Ra-226 in Urine

Sample ID	²¹⁷ At Yield (%)	²²⁶ Ra Reference Value (mBq smp ⁻¹)	²²⁶ Ra Measured Value * (mBq smp ⁻¹)	Difference (%)
1	91.2	73.8	75.5	2.4
2	88.3	73.8	88.4	19.9
3	84.7	73.8	82.7	12.1
4	85.0	73.8	76.8	4.1
5	94.1	73.8	78.3	6.1
6	87.7	73.8	78.7	6.6
7	99.2	73.8	70.8	-4.0
8	92.3	73.8	64.8	-12.1
Avg	90.3			
SD	4.9		77.0	4.4
% RSD	3.4			

8hr count time

At -217 ingrowth to mid-point-17.75 hrs.

*corrected for 5.0 mBq Ra-226 in blank urine

Cation Resin + Ln Resin

Summary

- **New method for Ra226 has been developed at SRS**
 - for solid samples and urine
- **Rapid fusion**
 - Ca /Ra carbonate precipitation
 - Ba Removal (if needed)
 - Ra-225 (At-217)
- **Cation + (Sr resin) + Ln resin**
 - Good yields, column flow
 - Effective removal of Ac-225
 - Alpha spectrometry
 - ICP-MS option