





Laboratoire de  
**Radioécologie**

## SEPARATION AND ANALYSIS OF Sr-90 AND Zr-90 FOR NUCLEAR FORENSIC APPLICATIONS



**Dominic Larivière, Université Laval**  
**Derek McLain, Argonne National Laboratory**



UNIVERSITÉ  
**LAVAL**  
Département de chimie

**Argonne**   
NATIONAL LABORATORY

# Acknowledgments

---



**Ana Paula Zattoni**  
MSc in May 2015  
Researcher at Cambrian Solutions Inc



**Agilent Technologies**

A portion of this material is based upon work supported by the U.S. Department of Homeland Security under contract number HSHQDN-16-X-00011. The views and conclusions contained in this document are those of the authors and should not be interpreted as representing the official policies, either expressed or implied, of the U.S. Department of Homeland Security.



# Availability of $^{90}\text{Sr}$ sources

## RTG

- $^{90}\text{SrTiO}_3$
- Fairly common in Russia
- Activity  $\sim 10^{15}$  Bq/RTG

## Nuclear waste

Fission product of  $^{235}\text{U}$  (reactors)

## Medical isotopes

- Production of  $^{90}\text{Y}$
- Cancer treatment



# Radiological properties

Parameters	Information
Appearance	Type of material (e.g. powder)
Dimensions and isotopic composition	Type of reactor
Impurities	Process, Geolocation
Surface	Production plant
$^{18}\text{O}/^{16}\text{O}$	Geolocation
Microstructure	Process
Age	Production date



# Radiological properties


Parameters	Information
Appearance	Type of material (e.g. powder)
Dimensions and isotopic composition	Type of reactor
Impurities	Process, Geolocation
Surface	Production plant
$^{18}\text{O}/^{16}\text{O}$	Geolocation
Microstructure	Process
Age	Production date



**Origin**



# Decay scheme of $^{90}\text{Sr}$




**Radiometric approaches**

**Mass spectrometric approaches**



# Radiochronometry



## Age

Time since the last purification

**Production or purification  $\rightarrow t = 0$**

100 %  $^{90}\text{Sr}$



0%  $^{90}\text{Zr}$

$$t = \frac{28,9}{\ln 2} \ln \left( 1 + \frac{[^{90}\text{Zr}]}{[^{90}\text{Sr}]} \right)$$




# Choices


Ion-exchange resin → Charged active site



Extraction resins → Chelating agents



**Sr Resin**



**DGA**

F. W. E. Strelow, *Analytical Chemistry* **1959**, 31, 1974-1977

L. R. Bunney, N. E. Ballou, J. Pascual, S. Foti, *Analytical Chemistry* **1959**, 31, 324-326


E. P. Horwitz, R. Chiarizia, M. L. Dietz, *Solvent Extraction and Ion Exchange* **1992**, 10, 313-336

D. W. Pawlak, J. L. Parus, T. D. A. Muklanowicz, R. Mikolajczak, *Talanta* **2013**, 114, 1-4






# Separation scheme




- E1 → *Partitioning coefficient very different*
- Preliminary tests → Stable isotopes
- Zr separation → Interference removal
- Zr recovery → Radioactive source determination



# Ion-exchange resins



# Extraction resins




## Zr in 1st fraction

- Faster
- Limited exposure to radioactivity

Sr Resin


— Sr  
—■ Zr





DGA Resin



# Limitations



# Limitations



# Zr recovery

## Sample preparation







Addition of  $F^-$  to complex Zr

## Anionic exchange resin

- Neutral or positive complexes
- Zr less retain  $\rightarrow$  effect on separation

## DGA Resin

- Neutral complex  $\rightarrow$  Zr more retained
- Increased affinity on the resin  $\rightarrow$  effect on recovery yield




**Which species will be favored?**



# F<sup>-</sup> - Dowex1-X8


Sample preparation → E1 + HF 0,01M  
Zr recovery yield =  $124 \pm 11$  %



Neutral complex forming!!!



# F<sup>-</sup> DGA-resin





# Comparison

Resin	Sr removal	Y removal	HF
AG50W-X8	✓	X	X
Dowex1-X8	✓	✓	X
Sr	✓	X	X
DGA	✓	✓	✓


## Adding HF in the sample

DGA resin → Excellent alternative for proper separation-recovery




# New Procedure

- DGA Resin
  - Wash
    - DIW
  - Pre-Condition
    - 1 M HNO<sub>3</sub>/0.2 M HF
  - Load
    - 1 M HNO<sub>3</sub>/0.2 M HF
  - Elute Zr
    - 1 M HNO<sub>3</sub>
  - Elute Sr
    - 0.05 M HNO<sub>3</sub>




# Separation Results – DGA Resin

DGA Resin Elution Profile



# Separation Results – DGA Resin

DGA Resin Elution Profile



# Recovery Comparison

**Sr Resin Separation % Recovery**

Element	Load Fraction	Rinse Fraction	Sr Elution Fraction
Zr	0.00 ± 0.00	100.19 ± 2.72	0.01 ± 0.02
Y	0.00 ± 0.00	99.91 ± 2.40	0.00 ± 0.00
Sr	0.01 ± 0.02	0.03 ± 0.01	97.33 ± 1.09

- Minor Sr cross-contamination
- Minor Zr cross-contamination
- Quantitative recovery of Zr/Y
- Nearly quantitative recovery of Sr
- RSD for calculated age: 2.99%

**DGA Resin Separation % Recovery**

Element	Load Fraction	Zr Elution Fraction	Sr Elution Fraction
Zr	0.00 ± 0.01	99.88 ± 0.43	0.12 ± 0.19
Y	0.93 ± 3.21	0.00 ± 0.00	0.00 ± 0.00
Sr	0.00 ± 0.00	0.00 ± 0.01	100.94 ± 0.57

- Minor Zr cross-contamination
- Quantitative recovery of Zr
- Quantitative recovery of Sr
- RSD for calculated age: 1.51%



# DGA Resin – High Fidelity Measurements

Sample ID	Sr in Sample (ng)	Y in Sample (ng)	Zr in Sample (ng)	Zr/Sr Ratio
Process Blank - Zr	<1.0	<1.0	<1.0	-
Process Blank - Sr	<1.0	<1.0	<1.0	
DGA Resin 1 - Zr	<1.0	<1.0	102.2	0.979
DGA Resin 1 - Sr	104.4	<1.0	<1.0	
DGA Resin 2 - Zr	<1.0	<1.0	100.7	0.959
DGA Resin 2 - Sr	105.1	<1.0	1.1	
DGA Resin 3 - Zr	<1.0	<1.0	101.9	0.986
DGA Resin 3 - Sr	103.3	<1.0	<1.0	
DGA Resin 4 - Zr	<1.0	<1.0	102.7	0.980
DGA Resin 4 - Sr	104.8	<1.0	<1.0	
DGA Resin 5 - Zr	<1.0	<1.0	101.1	0.984
DGA Resin 5 - Sr	102.7	<1.0	<1.0	
			Average	0.978
			SD	0.011
DGA Resin Spike 1	99.9	43.8	98.2	0.983
DGA Resin Spike 2	100.2	43.9	98.4	0.983
DGA Resin Spike 3	100.3	43.9	98.3	0.980
			Average	0.982
			SD	0.002



# Conclusions/Future Plans

---

- Conclusions

- DGA Resin procedure works at least as well as Sr Resin separation

- Load and Elute Zr in 1 M  $\text{HNO}_3$ /0.2 M HF

- Resin requires cleaning

- ~1.5 ng Zr/g DGA resin

- Future Work

- Test re-usability

- Likely possible with 0.05 M HCl



