

Certificate of Analysis

No. CA 18106

Product: Eichrom's Anion Exchange Resin
Lot #: 070318A4M
Date of Manufacture: 3-JUL-2018

Eichrom's Anion Exchange resin Lot # 070318A4M has been manufactured according to Eichrom standard operating procedures. This resin has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

| <u>Test Value</u> | <u>Specification</u> | <u>Result</u> |
|---------------------|------------------------------|---------------|
| Appearance | White to light Yellow | Pass |
| Extractables | < 1 mg (residue) / g (resin) | n.d. |
| pH | ≥ 4 | Pass |
| Wet Volume Capacity | ≥ 1.0 meq/mL | 1.2 |



Jill Bryant
Quality System Coordinator

July 6, 2018