

Certificate of Analysis
No. CA 18129

Product: Eichrom's Cation Exchange Resin; 50Wx8, 50-100 mesh, Hydrogen Form


Date of Manufacture: 31-JUL-2018

Lot #: 073118C8C

Eichrom's Cation Exchange resin, Lot # 073118C8C has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

Test Value	Specification	Result
Extractables	< 1 mg(residue) / g(resin)	0.08
pH	≥ 4.0	pass
Wet Volume Capacity	min. 1.7 meq/mL	2.1



Jill Bryant
Quality System Coordinator

July 31, 2018