




**Certificate of Analysis**  
No. CA 18105

Product: Eichrom's LN-3 Resin Cartridges  
Date of Manufacture: 19-APR-2016  
Lot #: L3SR16A

Eichrom's Ln3 Resin cartridge Lot L3SR16A has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, L3S110409, has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

**ANALYSIS RESULTS**

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
%La in load/rinse	95-105	> 99 %
%Be in load/rinse	95-105	> 99 %
%U in load/rinse	< 3	< 0.1
%La in (NH <sub>4</sub> ) <sub>2</sub> CO <sub>3</sub> strip	< 1	< 0.1
%Be in (NH <sub>4</sub> ) <sub>2</sub> CO <sub>3</sub> strip	< 1	< 0.1
%U in (NH <sub>4</sub> ) <sub>3</sub> CO <sub>3</sub> strip	95-105	98

  
\_\_\_\_\_  
Jill Bryant  
Quality System Coordinator

July 5, 2018