



# Certificate of Analysis

No. CA 18131

Product: Eichrom's Pb Resin  
Date of Manufacture: 31-JUL-2018  
Lot #: PBS073118

Eichrom's Pb Resin Lot # PBS073118, has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
% Lead in Load + Rinse	< 0.1 %	0.0 %
% Lead in Strip	> 99 %	100 %
Gravimetric Yield	80-120 %	96.7 %
Chemical Yield	80-120 %	92.5 %

*Jill Bryant*  
\_\_\_\_\_  
Jill Bryant  
Quality System Coordinator

August 6, 2018

