

Certificate of Analysis

No. CA 18188

Product: Eichrom's Anion Exchange Resin
 Lot #: 110218A4M
 Date of Manufacture: 2-NOV-2018

Eichrom's Anion Exchange resin Lot # 110218A4M has been manufactured according to Eichrom standard operating procedures. This resin has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White to light Yellow	Pass
Extractables	< 1 mg (residue) / g (resin)	n.d.
pH	< 4	Pass
Wet Volume Capacity	≥ 1.0 meq/mL	1.2


 Jill Bryant
 Quality System Coordinator

November 5, 2018