

Certificate of Analysis


No. CA 18227

Product: Eichrom's Anion Exchange Resin
Date of Manufacture: 3-DEC-2018
Lot #: 120318A8I

Eichrom's Anion Exchange resin Lot 120318A8I has been manufactured according to Eichrom standard operating procedures. This resin has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White to light Orange	Pass
Wet Volume Capacity	≥ 1.5 meq/mL	1.5



Jill Bryant
Quality System Coordinator

December 18, 2018