

# Certificate of Analysis

No. CA 24151

Product: Eichrom's Cation Exchange Resin Cartridges; 50Wx8, 100-200 mesh, H+ Form

Lot #: C8MR24D

Date of Manufacture: 16-JUL-2024

Eichrom's Cation Exchange resin cartridge Lot C8MR24D has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot, 060324C8M, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

Test Value	Specification	Result
Extractables	< 1 mg(residue) / g(resin)	0.16 mg/g
pH	> 4.0	pass
Wet Volume Capacity	min. 1.7 meq/mL	2.0 meq/mL (Per vendor cert)

  
Jose Munoz  
Quality Analyst

July 16, 2024