

Certificate of Analysis


No. CA 24155

Product: Eichrom's TEVA[®] Discs
Lot #: TEFD24A
Date of Manufacture: 19-JUL-2024

Eichrom's TEVA[®] Discs, Lot # TEFD24A has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	47mm, round	Pass
Tc-99 yield in pH 2 (100-200 mL/min)	90 – 110 %	99.6 %



Jose Munoz
Quality Analyst

July 19, 2024