

# Certificate of Analysis

No.

**Product:**

**Part Number:**

**Lot Number:**

**Manufactured Date:**

This product has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, \_\_\_\_\_, has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

All products come with a 1-year warranty from purchase, unless labeled differently. Many products may be used beyond the 1-year warranty. It is suggested that older material be reevaluated according to customer's policies and procedures to determine suitability for continued use.

<u>Test</u>	<u>Analysis</u>	<u>Specification</u>	<u>Result</u>
Appearance		White Beads	Pass
% <sup>233</sup> U stripped in 0.01M HNO <sub>3</sub>			
Max. acceptable mean		2.0	
Max. standard deviation		1.0	
% <sup>241</sup> Am stripped in 0.01M HNO <sub>3</sub>			
Min. acceptable mean		90.0	
Max. standard deviation		2.5	
% <sup>233</sup> U stripped in 0.5M HNO <sub>3</sub>			
Min. acceptable mean		95.0	
Max. standard deviation		2.5	
% <sup>241</sup> Am stripped in 0.5M HNO <sub>3</sub>			
Max. acceptable mean		2.0	
Max. standard deviation		1.0	