

Certificate of Analysis

No.

Product:

Part Number:

Lot Number:

Manufactured Date:

This product has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, _____, has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

All products come with a 1-year warranty from purchase, unless labeled differently. Many products may be used beyond the 1-year warranty. It is suggested that older material be reevaluated according to customer's policies and procedures to determine suitability for continued use.

Analysis

| <u>Test</u> | <u>Specification</u> | <u>Result</u> |
|--|----------------------|---------------|
| Appearance | White Beads | Pass |
| Gravimetric Yield (%) | 80-120 | |
| Chemical Yield (%) | 80-120 | |
| % Pb in DI H ₂ O Strip | ≥ 99 | |
| % Pb in 1M HNO ₃ Load + 0.1M HNO ₃ Rinse | ≤ 0.1 | |